PERFORMING NO-GLARE BRAND ACCORDING TO INDEPENDENT EYECARE PROFESSIONALS*

NEW Crizal® Easy Pro

EASIER CLEANING

is amongst the top-5 criteria for consumers when choosing prescription eyeglasses⁽¹⁾ On average, eyeglass wearers wipe their lenses **20,000 TIMES** before renewing them⁽²⁾

That's almost **60 HOURS** every year⁽³⁾

THE TECHNOLOGY IN NEW CRIZAL® EASY PROVIDES OUTSTANDING EASINESS-TO-CLEAN TO EASILY REPEL SMUDGES, DUST, AND WATER.



Reinforces the resistance to smudges and water on both sides of the lens for easier and optimal cleaning.

Anti-Static Layer

The **anti-static layer** reduces static on the lens, allowing it to repel dust and other particulates and keeps lenses cleaner longer.



Crizal[®] Easy^{Pro} lenses went through a rigorous protocol of tests to assess the performance technically and in real-life simulations.

CONTACT ANGLE TEST

To validate hydrophobic & oleophobic properties the contact angle measurement is used. The higher it is, the better the smudge efficiency.

STAIN RELEASE TEST

20,000 wiping motion with ink splatter on the lens every 5,000 times.

To demonstrate the superior resistance of Crizal[®] Easy^{Pro} against water, oil, and dust vs. standard anti-reflective coatings, a series of situations were reproduced in real-life simulations.

Standard coating



WATER



OIL



DUST

Crizal® EasyPro



Crizal® Easy^{Pro} coating has high water droplet repellence performance based on contact angle criteria.



Crizal® Easy^{Pro} coating has superior ability to repel oil and smudges compared to main competitors⁽⁴⁾.



The anti-static properties of Crizal® Easy^{Pro} coating help to repel dust.

Significantly higher contact angle vs main competitors⁽⁴⁾ 3x easier to clean after extensive smudging vs main competitors⁽⁴⁾ Optimal protection against reflections, dust and UV rays

Crizal[®] Easy^{Pro} is easy to clean and thanks to its oleophobic, hydrophobic, and anti-static properties, the lenses stay cleaner longer.

CRIZAL® EASYPRO EASIER CLEANING, SIMPLER LIVING

Contact your Essilor Account Executive for more details.

