





Crizal®

Crizal® FOR THE CLEAREST VISION POSSIBLE™

CRIZAL® NO-GLARE LENSES PRODUCT GUIDE

With *Crizal* No-Glare lenses, patients also get a two-year warranty with up to two replacements in two years.*

	<i>Crizal</i> SAPPHIRE ^{360°UV}	<i>Crizal</i> PREVENCIA	<i>Crizal</i> AVANCÉ ^{UV}	<i>Crizal</i> ALIZÉ ^{UV}	<i>Crizal</i> easy ^{UV}
CLARITY	GLARE REDUCTION Virtually reflection-free for better vision during the day and safer driving at night.	BEST-IN-CLASS CLARITY	EXCEPTIONAL CLARITY		
		VERY FAINT BLUE HUE	BLUE-VIOLET REFLECTION	FAINT GREEN HUE	
DURABILITY	DUST REPELLENT Prevents particles from adhering to the lens for less frequent cleaning.	ANTI-PARTICULATE LAYER FOR PROTECTION AGAINST DUST AND DIRT			
	WATER REPELLENT Repels moisture for convenient outdoor wear.	SUPERIOR CLEANABILITY: HIGH SURFACE DENSITY™ TECHNOLOGY MAKES LENS SUPER HYDROPHOBIC + THE PATENTED PAD-CONTROL SYSTEM™			EASY TO CLEAN WITH PATENTED PAD-CONTROL SYSTEM™
	SCRATCH AND SMUDGE RESISTANCE Protects lenses from the wear and tear of everyday use and cleaning. Reduces smears and fingerprints for clarity without distractions.	CONTACT ANGLE OF 116°			CONTACT ANGLE OF 112°
		2X MORE SCRATCH RESISTANCE WITH SR BOOSTER™ LAYER†		BETTER SCRATCH RESISTANCE THAN ORDINARY SCRATCH-COATED LENSES	
PROTECTION	UV PROTECTION Essential daily protection against UV light.				
	HARMFUL BLUE LIGHT** PROTECTION Crizal® Prevenia® selectively filters out Harmful Blue Light for constant protection.	100% UV PROTECTION ON FRONT SIDE OF LENS*** WITH ADDITIONAL PROTECTION FROM UV REFLECTION OFF BACK SIDE OF LENS			
			FILTRATION OF UP TO 20% OF HARMFUL BLUE-VIOLET LIGHT		

With *Crizal* lenses

Without *Crizal* lenses



SEE MORE. DO MORE.



*Lenses are guaranteed for up to two years from original order date from product failure resulting in atypical scratching or crazing from normal use.


Harmful Blue Light is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells. *With the exception of clear plastic (1.50) substrates. †As compared to Crizal Alize.
©2019 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries. These products may be protected by one or more patents listed at www.essilorusa.com/patents. 54930_PRO_ZAL SHK/SSX 10/19

Simulated Images

Crizal lenses combine multiple layers and technologies to offer the best combination of clarity, durability, and protection.

FOR THE CLEAREST VISION POSSIBLE™

CRIZAL FOR KIDS



Crizal *SAPPHIRE* 360° UV

- REDUCES GLARE
- RESISTS SCRATCHES
- REPELS DUST
- REPELS WATER
- RESISTS SMUDGES
- OPTIMAL UV PROTECTION (E-SPF 35)
- Best in Clarity, Scratch Resistance, Cleanability, and UV Protection

Crizal *PREVENCIA*

- REDUCES GLARE
- RESISTS SCRATCHES
- REPELS DUST
- REPELS WATER
- RESISTS SMUDGES
- HARMFUL BLUE LIGHT* FILTRATION
- OPTIMAL UV PROTECTION (E-SPF 35)
- Best in UV Protection, Harmful Blue Light Filtration, and Cleanability

Crizal *AVANCE* UV™

- REDUCES GLARE
- RESISTS SCRATCHES
- REPELS DUST
- REPELS WATER
- RESISTS SMUDGES
- OPTIMAL UV PROTECTION (E-SPF 35)
- Best in Cleanability and Scratch Resistance

Crizal *ALIZE* UV™

- REDUCES GLARE
- RESISTS SCRATCHES
- REPELS WATER
- RESISTS SMUDGES
- COMPREHENSIVE UV PROTECTION (E-SPF 25)
- Best in Cleanability

Crizal *easy* UV™

- REDUCES GLARE
- RESISTS SCRATCHES
- REPELS WATER
- RESISTS SMUDGES
- COMPREHENSIVE UV PROTECTION (E-SPF 25)
- Everyday Essential No-Glare Lens

Crizal *PREVENCIA* kids

- REDUCES GLARE
- RESISTS SCRATCHES
- REPELS WATER
- RESISTS SMUDGES
- HARMFUL BLUE LIGHT* FILTRATION
- OPTIMAL UV PROTECTION
- Best in UV Protection and Harmful Blue Light Filtration for Kids

Crizal *kids* UV™

- RESISTS SCRATCHES
- REDUCES GLARE
- REPELS WATER
- RESISTS SMUDGES
- COMPREHENSIVE UV PROTECTION
- Everyday Essential No-Glare Lens for Kids

Crizal for Kids No-Glare lenses are a no-glare treatment, polycarbonate lens material package. They are **impact resistant** and **virtually unbreakable**—perfect for kids' active lives!

Image represents the Crizal Sapphire® 360° UV No-Glare lens treatment stack. Thickness and order of layers is for demonstration purposes only.

In a clear lens, E-SPF® is an index rating of the overall UV protection of a lens. E-SPF® was developed by Essilor International and endorsed by third-party experts. The E-SPF® index relates to lens performance only and excludes direct eye exposure that depends on external factors (wearer's morphology, frame shape, position of wear). *Harmful Blue Light is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells.

Recommended for ages 5-12. Crizal Kids UV™ is available in polycarbonate lens material only.