

A 100% DUALOPTIX™ DIGITAL FAMILY OF PRODUCTS



The right product for your customers

| | DEFINITY® STANDARD | DEFINITY ³ CUSTOMIZED | DEFINITY ³ ^{PLUS} PERSONALIZED |
|---|-----------------------|-------------------------------------|---|
| Computer automatically selects between Regular and Short design | | ✓ | ✓ |
| Fixed split of the ADD: +0.75 on the back, remaining ADD power on the front | ✓ | | |
| Dynamic split of the ADD individually tailored to the patient's Rx | | ✓ | ✓ |
| Optimized to frame choice with A, B, PD, & Seg Height fitting measurements | | ✓ | ✓ |
| Incorporates 3 additional fitting parameters for a personalized fit: Panto, Wrap and Vertex | | | ✓ |

Keep your patients satisfied as presbyopia progresses, with quite possibly the only PAL they'll ever need

- 35% less unwanted astigmatism vs competing PALs
- 4 out of 5 patients adapt immediately
- Outperforms the competition, especially at high ADD powers

Discover a distinctively different experience with DEFINITY lenses—unique designs, superior support and an unparalleled community of dispensers. Visit definitycommunity.com.

DEFINITY® DIGITAL LENSES—AVAILABILITY CHART

NEW DEFINITY® 3 PLUS

Minimum fitting height: 15mm. Max Ground View Advantage™ benefits at fitting heights ≥ 19mm

| Material | 1.50 Plastic | TREXA™ Trivex® | Airwear® Polycarbonate | Thin&Lite® | | |
|-------------|-----------------|-------------------|---------------------------|---------------|-----------------|-------------------|
| | | | | 1.60 | 1.67 | 1.74 ¹ |
| Power Range | +6.00 to -10.00 | +6.00 to -10.00 | +6.00 to -10.00 | Not Available | +9.00 to -12.00 | +9.00 to -14.00 |
| Cylinder | up to -6.00 | up to -6.00 | up to -6.00 | | up to -6.00 | up to -6.00 |
| ADD Power | +1.00 to +3.50 | +1.00 to +3.50 | +1.00 to +3.50 | | +1.00 to +3.50 | +1.00 to +3.50 |

LENS OFFERINGS

| | | | | | | |
|------------------------------------|--------------|--------------|--------------|---------------|--------------|-----|
| Clear | ● | ● | ● | Not Available | ● | ● |
| Xperio UV™ Polarized ² | N/A | N/A | Gray | | N/A | N/A |
| Transitions® Signature™ VII Lenses | Gray & Brown | Gray & Brown | Gray & Brown | | Gray & Brown | N/A |

NEW DEFINITY® 3

Minimum fitting height: 15mm. Max Ground View Advantage benefits at fitting heights ≥ 19mm

| Material | 1.50 Plastic | TREXA™ Trivex® | Airwear® Polycarbonate | Thin&Lite® | | |
|-------------|-----------------|-------------------|---------------------------|---------------|-----------------|-------------------|
| | | | | 1.60 | 1.67 | 1.74 ¹ |
| Power Range | +6.00 to -10.00 | +6.00 to -10.00 | +6.00 to -10.00 | Not Available | +9.00 to -12.00 | +9.00 to -14.00 |
| Cylinder | up to -6.00 | up to -6.00 | up to -6.00 | | up to -6.00 | up to -6.00 |
| ADD Power | +1.00 to +3.50 | +1.00 to +3.50 | +1.00 to +3.50 | | +1.00 to +3.50 | +1.00 to +3.50 |

LENS OFFERINGS

| | | | | | | |
|------------------------------------|--------------|--------------|--------------|---------------|--------------|-----|
| Clear | ● | ● | ● | Not Available | ● | ● |
| Xperio UV™ Polarized ² | N/A | N/A | Gray | | N/A | N/A |
| Transitions® Signature™ VII Lenses | Gray & Brown | Gray & Brown | Gray & Brown | | Gray & Brown | N/A |

DEFINITY® / DEFINITY SHORT™

Regular: min fitting height: 18mm. Max Ground View Advantage benefits at fitting heights ≥ 22mm
Short corridor: min fitting height: 15mm. Max Ground View Advantage benefits at fitting heights ≥ 19mm

| Material | 1.50 Plastic | TREXA™ Trivex® | Airwear® Polycarbonate | Thin&Lite® | | |
|-------------|-----------------|-----------------------|---------------------------|-----------------|-----------------------|------------------------|
| | | | | 1.60 | 1.67 | 1.74 ¹ |
| Power Range | +6.00 to -10.00 | +6.00 to -10.00 | +4.00 to -7.00 | +8.00 to -12.00 | +10.00 to -12.00 | +9.00 to -14.00 |
| Cylinder | up to -4.00 | up to -4.00 | up to -4.00 | up to -4.00 | up to -4.00 | up to -4.00 |
| ADD Power | +1.00 to +3.00 | +1.00 to +3.50 & 4.00 | +1.00 to +3.00 | +1.00 to +3.00 | +1.00 to +3.50 & 4.00 | +1.00 to +3.00 & +3.50 |

LENS OFFERINGS

| | | | | | | |
|------------------------------------|--------------|------|--------------|---------------|------|-----|
| Clear | ● | ● | ● | ● | ● | ● |
| Xperio UV™ Polarized ² | Gray & Brown | N/A | Gray & Brown | Not Available | N/A | N/A |
| Transitions® Signature™ VII Lenses | Gray | Gray | Gray | | Gray | N/A |

¹ DEFINITY Thin&Lite 1.74 lenses will systematically include your choice of Crizal Sapphire UV™, Crizal Avancé UV™, or Crizal Alizé UV™

² May also be ordered in Essilor Polarized. Crizal SunShield UV™ is recommended with every DEFINITY polarized job.

To learn more visit definitycommunity.com

DEFINITY®
DISTINCTIVELY DIFFERENT®