

### IDEAL FAMILY: ESSILOR IDEAL, IDEAL FIT, IDEAL SHORT, IDEAL ADVANCED, IDEAL ADVANCED FIT, IDEAL SPORT, IDEAL COMPUTER, IDEAL DRIVING

MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER*	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25
	1.50 Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® Vantage® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	2.00, 3.25, 4.25, 6.25, 8.25
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions Signature VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions XTRActive (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50
	Airwear Transitions Vantage (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50
THIN & LITE®	Thin & Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin & Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin & Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin & Lite 1.67	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin & Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin & Lite 1.67 Transitions VI (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin & Lite 1.67 Transitions XTRActive (Gray)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin & Lite 1.67 Xperio UV (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00
	Thin & Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00
	Thin & Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00

\*Max add for Essilor Ideal Short is 3.00

### IDEAL FAMILY: IDEAL ADVANCED WRAP

MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE
PLASTIC	1.50 Index Plastic	-7.00 to +6.00	up to -6.00	+0.75 to +3.50	5.25, 6.25, 7.25, 8.75
	1.50 Transitions® Signature VII™ (Gray/Brown/GG)	-7.00 to +6.00	up to -6.00	+0.75 to +3.50	5.25, 6.25, 7.25, 8.75
	1.50 Transitions® VI (Gray/Brown)	-7.00 to +6.00	up to -6.00	+0.75 to +3.50	5.25, 6.25, 8.25
	1.50 Xperio UV™ (Gray/Brown)	-7.00 to +6.00	up to -4.00	+0.75 to +3.50	5.25, 6.25, 7.25, 8.75
AIRWEAR® (POLYCARBONATE)	Airwear®	-8.00 to +6.00	up to -6.00	+0.75 to +3.50	4.75, 5.50, 7.00, 8.50
	Airwear® Transitions® Signature VII™ (Gray/Brown/GG)	-8.00 to +6.00	up to -6.00	+0.75 to +3.50	4.75, 5.50, 7.00, 8.50
	Airwear® Transitions® VI (Gray/Brown)	-8.00 to +6.00	up to -4.00	+0.75 to +3.50	5.50, 7.00, 8.50
	Airwear® Xperio UV™ (Gray/Brown)	-8.00 to +6.00	up to -6.00	+0.75 to +3.50	4.25, 5.50, 7.00, 8.50
THIN & LITE®	Thin & Lite® 1.67	-9.00 to +6.00	up to -6.00	+0.75 to +3.50	4.75, 6.00, 6.75, 8.00
	Thin & Lite 1.67 Transitions® Signature VII™ (Gray/Brown/GG)	-9.00 to +6.00	up to -6.00	+0.75 to +3.50	4.75, 6.00, 6.75, 8.00
	Thin & Lite 1.67 Transitions® XTRActive (Gray)	-9.00 to +6.00	up to -6.00	+0.75 to +3.50	4.75, 6.00, 6.75, 8.00

### IDEAL FAMILY-DESIGN ENGRAVINGS

LENS DESIGN	1.50 PLASTIC	AIRWEAR®	TREXA™	THIN & LITE® 1.60	THIN & LITE® 1.67	THIN & LITE® 1.74	MFHmm
ESSILOR IDEAL®	E	E P	E *	E 6	E 67	E 74	14
ESSILOR IDEAL® FIT	EF	EF P	EF *	EF 6	EF 67	EF 74	14
ESSILOR IDEAL® SHORT	E	E P	E *	E 6	E 67	E 74	14
ESSILOR IDEAL® ADVANCED	E3	E3 P	E3 *	E3 6	E3 67	E3 74	14
ESSILOR IDEAL® ADVANCED WRAP	E3 WR	E3 P WR	E3 * WR	E3 6 WR	E3 67 WR	E3 74 WR	14
ESSILOR IDEAL® ADVANCED FIT	E3F	E3FP	E3F *	E3F 6	E3F 67	E3F 74	14
ESSILOR IDEAL® DRIVING	ED	ED P	ED *	ED 6	ED 67	ED 74	14
ESSILOR IDEAL® SPORT	EO	EO P	EO *	EO 6	EO 67	EO 74	18
ESSILOR IDEAL® COMPUTER	EC	EC P	EC *	EC 6	EC 67	EC 74	13