



Eyezen™ +

Fitting Eyezen™+ Lenses

Eyezen+ lenses require the following standard measurements:

- Monocular Distance PD
- Fitting Height: Simply measure from center pupil to the lowest point of the lens to acquire fitting height. For **Eyezen+**, minimum fitting height is **15mm**.

Eyezen+ lenses with embedded Harmful Blue Light* protection are also recommended with Crizal® No-Glare treatment for our most comprehensive eye protection.

The following laser engravings located on the nasal side of each lens will identify the **Eyezen** design received.

Eyezen+ 0



Eyezen+ 1



Eyezen+ 2



Eyezen+ 3



Eyezen+ 4



Embedded **Harmful Blue Light*** Protection.

* *Eyezen* lenses block at least 20% of Harmful Blue Light, which is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells.

EYEZEN™+ LENSES AVAILABILITY

for Eyecare Professionals

Effective **DECEMBER 2019**

| MATERIAL | LIGHT MANAGEMENT | INDEX | COLOR | POWER RANGE | BASE CURVE | MAX CYLINDER |
|--------------|---------------------------------------|-------|--|------------------|--|--------------|
| PLASTIC | Clear | 1.5 | | -10.00 to +7.50 | 2.00, 3.50, 4.25, 6.25, 7.25, 8.75 | -6.00 |
| | Essilor Polarized | 1.5 | Gray, Brown, Gray-Green, Gray-1, Brown-1 | -10.00 to +7.50 | 2.00, 3.25, 4.25, 6.25, 7.25, 8.25 | -6.00 |
| | Xperio UV™ Polarized | 1.5 | Gray, Brown, Gray-Green, Gray-1, Brown-1 | -10.00 to +7.50 | 2.00, 3.25, 4.25, 6.25, 7.25, 8.25 | -6.00 |
| | Transitions® Signature® GEN 8™ | 1.5 | Gray, Brown, Graphite Green | -10.00 to +7.50 | 2.00, 3.50, 4.25, 6.25, 7.25, 8.75 | -6.00 |
| | Transitions® Signature® Style Colors | 1.5 | Amber, Amethyst, Emerald, Sapphire | -10.00 to +7.50 | 2.00, 3.50, 4.25, 6.25, 7.25, 8.75 | -6.00 |
| | Transitions® Vantage® | 1.5 | Gray | -10.00 to +7.50 | 2.25, 4.25, 5.25, 8.25 | -6.00 |
| | Transitions® XTRActive® Gray* | 1.5 | Gray | -10.00 to +7.50 | 2.25, 4.25, 5.25, 6.25, 8.25 | -6.00 |
| | Transitions® XTRActive® Brown | 1.5 | Brown | -10.00 to +7.50 | 2.00, 3.50, 4.25, 6.25, 7.25, 8.75 | -6.00 |
| TREXA™ | Transitions® XTRActive® Style Mirrors | 1.5 | Blue, Green, Gold, Pink, Red, Silver Shadow | -10.00 to +7.50 | 2.00, 3.50, 4.25, 6.25, 7.25, 8.75 | -6.00 |
| | Clear § | 1.53 | | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25, 8.75 | -6.00 |
| AIRWEAR® | Transitions® Signature® GEN 8™ § | 1.53 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25, 8.75 | -6.00 |
| | Clear | 1.59 | | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | -6.00 |
| | Essilor Polarized | 1.59 | Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50 | -6.00 |
| | Xperio UV™ Polarized | 1.59 | Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50 | -6.00 |
| | Transitions® Signature® GEN 8™ | 1.59 | Gray, Brown, Graphite Green | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | -6.00 |
| | Transitions® Signature® Style Colors | 1.59 | Amber, Amethyst, Emerald, Sapphire | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | -6.00 |
| | Transitions® Vantage®* | 1.59 | Gray | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | -6.00 |
| | Transitions® XTRActive® | 1.59 | Gray, Brown | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | -6.00 |
| THIN & LITE® | Transitions® XTRActive® Style Mirrors | 1.59 | Blue, Green, Gold, Pink, Red, Silver Shadow | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | -6.00 |
| | Colors™* | 1.59 | | -10.00 to +8.00 | 2.25, 4.25, 6.50, 7.50, 8.50 | -6.00 |
| | Clear | 1.6 | | -12.00 to +8.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.50 | -6.00 |
| | Transitions® Signature® GEN 8™ | 1.6 | Gray, Brown | -12.00 to +8.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.50 | -6.00 |
| | Clear | 1.67 | | -14.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | -6.00 |
| | Essilor Polarized | 1.67 | Gray, Brown, Gray-Green | -14.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.75, 8.00 | -4.00 |
| | Xperio UV™ Polarized | 1.67 | Gray, Brown, Gray-Green | -14.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.75, 8.00 | -4.00 |
| | Transitions® Signature® GEN 8™ | 1.67 | Gray, Brown, Graphite Green | -14.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | -6.00 |
| | Transitions® Signature® Style Colors | 1.67 | Amber, Amethyst, Emerald, Sapphire | -14.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | -6.00 |
| | Transitions® XTRActive® | 1.67 | Gray, Brow | -14.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | -6.00 |
| | Transitions® XTRActive® Style Mirrors | 1.67 | Blue, Green, Gold, Pink, Red, Silver Shadow | -14.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | -6.00 |
| | Clear | 1.74 | | -20.00 to +13.00 | 0.75,1.75,2.75,3.75,5.25,6.50,8.00,10.00 | -6.00 |
| | Transitions® Signature® GEN 8™ | 1.74 | Gray | -20.00 to +13.00 | 0.75,1.75,2.75,3.75,5.25,6.50,8.00,10.00 | -6.00 |

* Transitions Vantage, AIRWEAR® Colors™, 1.50 plastic Transitions® XTRActive® Gray available on Eyezen+ 0-3 only

EYEZEN LENS DESIGNS

| | |
|------------------|------------------------------------|
| Eyezen+ 0 | patients age 17 and under |
| Eyezen+ 1 | patients age 18 to 34 |
| Eyezen+ 2 | patients age 35 to 39 |
| Eyezen+ 3 | patients age 40 to 44 |
| Eyezen+ 4 | patients age 45 to 50 ¹ |

There are five unique **Eyezen+** designs. Although the recommended age ranges for each design provided below are based on research, as always, the Eyecare Professional is free to use discretion in prescribing **Eyezen+** enhanced single vision lenses.

¹ Eyezen+ lenses are not adapted to patients that have been prescribed progressive lenses.



EyezenTM START

Fitting EyezenTM Lenses

Eyezen START lenses require the following standard measurements:

- Monocular Distance PD
- Fitting Height: Simply measure from center pupil to the lowest point of the lens to acquire fitting height.
There is no minimum fitting height

Eyezen START lenses with embedded Harmful Blue Light* protection are always recommended with Crizal[®] No-Glare treatment for our most comprehensive eye protection.

The following laser engraving, located on the nasal side of each lens, will identify the **Eyezen START** design received.

Eyezen Start



Embedded **Harmful Blue Light* Protection.**

* Eyezen START lenses block at least 20% of Harmful Blue Light, which is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells.

EYEZEN™ START LENSES AVAILABILITY

for Eyecare Professionals

Effective **DEC 2019**

| MATERIAL | LIGHT MANAGEMENT | INDEX | COLOR | POWER RANGE | BASE CURVE | MAX CYLINDER |
|----------------------------------|--|-------|--|---|---|--------------|
| PLASTIC | Clear | 1.5 | | -10.00 to +6.00 | 2.00, 3.50, 4.25, 5.25, 6.25, 7.25 | -6.00 |
| | Essilor Polarized § | 1.5 | Gray, Brown, Gray-Green, Gray-1, Brown-1 | -10.00 to +6.00 | 2.00, 3.25, 4.25, 6.25, 7.25 | -4.00 |
| | Xperio UV™ Polarized § | 1.5 | Gray, Brown, Gray-Green, Gray-1, Brown-1 | -10.00 to +6.00 | 2.00, 3.25, 4.25, 6.25, 7.25 | -4.00 |
| | Transitions® Signature® GEN 8™ | 1.5 | Gray, Brown, Graphite Green | -10.00 to +6.00 | 2.00, 3.50, 4.25, 5.25, 6.25, 7.25 | -6.00 |
| | Transitions® Signature® GEN 8™ lenses style colors § | 1.5 | Amber, Amethyst, Emerald, Sapphire | -10.00 to +6.00 | 2.00, 3.50, 4.25, 5.25, 6.25, 7.25 | -6.00 |
| | Transitions® XTRActive® § | 1.5 | Brown | -10.00 to +6.00 | 2.00, 3.50, 4.25, 5.25, 6.25, 7.25 | -6.00 |
| | Transitions® XTRActive® Style Mirrors § | 1.5 | Blue, Green, Gold, Pink, Red, Silver Shadow | -10.00 to +6.00 | 2.00, 3.50, 4.25, 5.25, 6.25, 7.25 | -6.00 |
| TREXA™ | Clear § | 1.53 | | -10.00 to +6.00 | 2.00, 3.50, 4.25, 5.25, 6.25, 7.25 | -6.00 |
| | Transitions® Signature® GEN 8™ § | 1.53 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.50, 4.25, 5.25, 6.25, 7.25 | -6.00 |
| AIRWEAR® | Clear | 1.59 | | -10.00 to +6.00 | 1.50, 3.25, 4.25, 4.75, 5.50, 7.00 | -6.00 |
| | Essilor Polarized § | 1.59 | Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby | -10.00 to +6.00 | 1.50, 3.25, 4.25, 5.50, 7.00 | -6.00 |
| | Xperio UV™ Polarized § | 1.59 | Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby | -10.00 to +6.00 | 1.50, 3.25, 4.25, 5.50, 7.00 | -6.00 |
| | Transitions® Signature® GEN 8™ | 1.59 | Gray, Brown, Graphite Green | -10.00 to +6.00 | 1.50, 3.25, 4.25, 4.75, 5.50, 7.00 | -6.00 |
| | Transitions® Signature® GEN 8™ lenses style colors § | 1.59 | Amber, Amethyst, Emerald, Sapphire | -10.00 to +6.00 | 1.50, 3.25, 4.25, 4.75, 5.50, 7.00 | -6.00 |
| | Transitions® XTRActive® | 1.59 | Gray, Brown | -10.00 to +6.00 | 1.50, 3.25, 4.25, 4.75, 5.50, 7.00 | -6.00 |
| | Transitions® XTRActive® Style Mirrors | 1.59 | Blue, Green, Gold, Pink, Red, Silver Shadow | -10.00 to +6.00 | 1.50, 3.25, 4.25, 4.75, 5.50, 7.00 | -6.00 |
| THIN & LITE® | Clear | 1.6 | | -12.00 to +8.00 | 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50 | -6.00 |
| | Transitions® Signature® GEN 8™ | 1.6 | Gray, Brown | -12.00 to +8.00 | 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50 | -6.00 |
| | Clear | 1.67 | | -14.00 to +9.00 | 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00 | -6.00 |
| | Essilor Polarized § | 1.67 | Gray, Brown, Gray-Green | -12.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.75, 8.00 | -4.00 |
| | Xperio UV™ Polarized § | 1.67 | Gray, Brown, Gray-Green | -12.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.75, 8.00 | -4.00 |
| | Transitions® Signature® GEN 8™ | 1.67 | Gray, Brown | -14.00 to +9.00 | 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00 | -6.00 |
| | Transitions® Signature® GEN 8™ lenses style colors § | 1.67 | Amber, Amethyst, Emerald, Sapphire | -14.00 to +9.00 | 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00 | -6.00 |
| | Transitions® XTRActive® § | 1.67 | Gray, Brown | -14.00 to +9.00 | 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00 | -6.00 |
| | Transitions® XTRActive® Style Mirrors § | 1.67 | Blue, Green, Gold, Pink, Red, Silver Shadow | -14.00 to +9.00 | 1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00 | -6.00 |
| | Clear § | 1.74 | | -20.00 to +13.00 | 0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00 | -6.00 |
| Transitions® Signature® GEN 8™ § | 1.74 | Gray | -20.00 to +13.00 | 0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00 | -6.00 | |

EYEZEN START LENS DESIGNS

| LENS ENGRAVING | |
|----------------|---|
| PRODUCT | SYMBOL |
| Eyezen Start |  |

§ Available on March 10, 2020