

ESSILOR  
#1

LENS BRAND RECOMMENDED  
BY EYECARE PROFESSIONALS  
WORLDWIDE\*

# Life has evolved. Why not your lenses?

Effective April 2022

## EYEZEN® ENHANCED SINGLE VISION LENSES AVAILABILITY CHART



ESSILOR

Eyezen®

\* Euromonitor, Eyewear 2022 edition; Essilor International company; worldwide retail value sales at RSP. Prescription lenses defined as Spectacle Lenses in Passport.

SEE MORE. DO MORE.

## FITTING EYEZEN® LENSES

Eyezen lenses require the following standard measurements:

- Monocular Distance PD
- Fitting Height: Simply measure from center pupil to the lowest point of the lens to acquire fitting height. *Eyezen+* 1, 2, 3, and 4 lenses have a minimum fitting height of 15mm. There is no minimum fitting height for *Eyezen Start* and *Eyezen Kids* lenses.

Eyezen lenses come with embedded Harmful Blue Light<sup>1</sup> protection and should also be paired with a Crizal® No-Glare coating for additional clarity, durability, and protection.

## EYEZEN LENS ENGRAVINGS

The following laser engravings located on the nasal side of each lens will identify the *Eyezen* lens design received.

### EYEZEN START LENS DESIGN

Eyezen Start



### EYEZEN KIDS LENS DESIGN

Eyezen Kids



### EYEZEN+ LENS DESIGNS

Eyezen+ 0



Eyezen+ 1



Eyezen+ 2



Eyezen+ 3



Eyezen+ 4



There are five unique *Eyezen+* lens designs. While the recommended age ranges for each design are based on research and provided below, the Eyecare Professional is free to use discretion when prescribing *Eyezen+* enhanced single vision lenses.

<i>Eyezen+</i> 0	patients age 17 and under
<i>Eyezen+</i> 1	patients age 18 to 34
<i>Eyezen+</i> 2	patients age 35 to 39
<i>Eyezen+</i> 3	patients age 40 to 44
<i>Eyezen+</i> 4	patients age 45 to 50 <sup>2</sup>

<sup>1</sup> Eyezen lenses filter at least 20% of Harmful Blue Light, which is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells.

<sup>2</sup> Eyezen+ lenses are not adapted to patients that have been prescribed progressive lenses.

©2022 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries. Transitions, Transitions Signature and XTRActive are registered trademarks of Transitions Optical, Inc. used under license by Transitions Optical Ltd. Photochromic performance is influenced by temperature, UV exposure and lens material. These products may be protected by one or more patents listed at [www.essilorusa.com/patents](http://www.essilorusa.com/patents). 218500\_PRO\_ZEN\_CVD\_03/22



ESSILOR

Eyezen®

# EYEZEN® START™ LENSES AVAILABILITY for Eyecare Professionals

Effective April 2022

MATERIAL	LIGHT MANAGEMENT	INDEX	COLOR	POWER RANGE	BASE CURVE	MAX CYLINDER
PLASTIC	Clear	1.5		-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Essilor Polarized	1.5	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 7.25	-4.00
	Xperio UV™ Polarized	1.5	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 7.25	-4.00
	Transitions® Signature® GEN 8™	1.5	Gray, Brown, Graphite Green	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® Signature® GEN 8™ lenses style colors	1.5	Amber, Amethyst, Emerald, Sapphire	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® XTRActive® new generation	1.5	Gray, Brown, Graphite Green	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® XTRActive® Polarized™		Gray	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
TREXA™	Clear	1.53		-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® Signature® GEN 8™	1.53	Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® XTRActive® new generation	1.53	Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
AIRWEAR®	Clear	1.59		-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Essilor Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50	-6.00
	Xperio UV™ Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50	-6.00
	Transitions® Signature® GEN 8™	1.59	Gray, Brown, Graphite Green	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Transitions® Signature® GEN 8™ lenses style colors	1.59	Amber, Amethyst, Emerald, Sapphire	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Transitions® XTRActive® new generation	1.59	Gray, Brown, Graphite Green	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Transitions® XTRActive® Style Mirrors	1.59	Blue, Green, Gold, Pink, Red, Silver Shadow	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Transitions® XTRActive® Polarized™	1.59	Gray	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
THIN & LITE®	Clear	1.6		-12.00 to +8.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50	-6.00
	Transitions® Signature® GEN 8™	1.6	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50	-6.00
	Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Essilor Polarized	1.67	Gray, Brown, Gray-Green	-12.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-4.00
	Xperio UV™ Polarized	1.67	Gray, Brown, Gray-Green	-12.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-4.00
	Transitions® Signature® GEN 8™	1.67	Gray, Brown, Graphite Green	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Transitions® Signature® GEN 8™ lenses style colors	1.67	Amber, Amethyst, Emerald, Sapphire	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive® new generation	1.67	Gray, Brown, Graphite Green	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive® Style Mirrors	1.67	Blue, Green, Gold, Pink, Red, Silver Shadow	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive® Polarized™	1.67	Gray	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Clear**	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00
	Transitions® Signature® GEN 8™**	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00

# EYEZEN® KIDS LENSES AVAILABILITY for Eyecare Professionals

Effective April 2022

MATERIAL	LIGHT MANAGEMENT	INDEX	COLOR	POWER RANGE	BASE CURVE	MAX CYLINDER
AIRWEAR®	Clear	1.59		-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Transitions® Signature® GEN 8™	1.59	Gray, Brown, Graphite Green	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00

# EYEZEN®+ LENSES AVAILABILITY for Eyecare Professionals

Effective April 2022

MATERIAL	LIGHT MANAGEMENT	INDEX	COLOR	POWER RANGE	BASE CURVE	MAX CYLINDER
PLASTIC	Clear	1.5		-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Essilor Polarized	1.5	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	-6.00
	Xperio UV™ Polarized	1.5	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	-6.00
	Transitions® Signature® GEN 8™	1.5	Gray, Brown, Graphite Green	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® Signature® lenses style colors	1.5	Amber, Amethyst, Emerald, Sapphire	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® XTRActive® new generation	1.5	Gray, Brown, Graphite Green	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® XTRActive® Polarized™	1.5	Gray	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
TREXA™	Clear	1.53		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® Signature® GEN 8™	1.53	Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® XTRActive® new generation		Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
AIRWEAR®	Clear	1.59		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Essilor Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	-6.00
	Xperio UV™ Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	-6.00
	Transitions® Signature® GEN 8™	1.59	Gray, Brown, Graphite Green	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® Signature® lenses style colors	1.59	Amber, Amethyst, Emerald, Sapphire	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® XTRActive® new generation	1.59	Gray, Brown, Graphite Green	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® XTRActive® Style Mirrors	1.59	Blue, Green, Gold, Pink, Red, Silver Shadow	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® XTRActive® Polarized™	1.59	Gray	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
THIN & LITE®	Clear	1.6		-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	Transitions® Signature® GEN 8™	1.6	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Essilor Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-4.00
	Xperio UV™ Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-4.00
	Transitions® Signature® GEN 8™	1.67	Gray, Brown, Graphite Green	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Transitions® Signature® lenses style colors	1.67	Amber, Amethyst, Emerald, Sapphire	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive® new generation	1.67	Gray, Brown, Graphite Green	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive® Style Mirrors	1.67	Blue, Green, Gold, Pink, Red, Silver Shadow	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive® Polarized™	1.67	Gray	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Clear**	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00
	Transitions® Signature® GEN 8™**	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00

\*\* Thin & Lite 1.74 lenses come systematic with your choice of Crizal Sapphire® 360° UV™, Crizal® Previncia®, Crizal® Rock™ or Crizal Alizé UV™.

\*\*\* AIRWEAR® Colors™ available on Eyezen+0-3 only