

Fitting Eyezen™ + Lenses

Eyezen+ lenses require the following standard measurements:

- Monocular Distance PD
- Fitting Height: Simply measure from center pupil to the lowest point of the lens to acquire fitting height. Recommended minimum fitting height is **15mm**.

Eyezen+ lenses with embedded Harmful Blue Light* protection are also recommended with Crizal® No-Glare treatment for the most comprehensive eye protection.

The following laser engravings located on the nasal side of each lens will identify the **Eyezen+** design received.

Eyezen+ 0



Eyezen+ 1



Eyezen+ 2



Eyezen+ 3



Embedded **Harmful Blue Light*** Protection.

*Eyezen+ lenses block at least 20% of Harmful Blue Light, which are the high energy waves found between 415 and 455nm (Blue-Violet light) and are suspected to be the most toxic to retinal cells.

Visit [EyezenPro.com/resources](https://www.eyezenpro.com/resources) for additional information and materials.

EYEZEN™+ LENSES AVAILABILITY

for Eyecare Professionals

Effective JULY 2017

MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	MAX CYLINDER
PLASTIC	1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	-6.00
	1.50 Index Essilor Polarized	1.50	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	-4.00
	1.50 Index Xperio UV™ Polarized	1.50	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	-4.00
	1.50 Index Transitions® Signature® VII	1.50	Gray, Brown, Graphite Green	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	-6.00
	1.50 Index Transitions® Vantage®	1.50	Gray	-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	-6.00
	1.50 Index Transitions® XTRActive®	1.50	Gray	-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	-6.00
AIRWEAR® (POLYCARBONATE)	1.50 Index Transitions® XTRActive®	1.50	Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	-6.00
	Clear	1.59		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Essilor Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Xperio UV™ Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® Signature® VII	1.59	Gray, Brown, Graphite Green	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® XTRActive®	1.59	Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
TREXA™ (TRIVEX®)	Airwear® Transitions® Vantage®	1.59	Gray	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Airwear Colors™	1.59	Gray, Brown, Gray-Green	-10.00 to +8.00	2.25, 4.25, 6.50, 7.50, 8.50	-6.00
	Clear	1.53		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	-6.00
THIN & LITE®	Transitions® Signature® VII	1.53	Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	-6.00
	1.60 Clear	1.60		-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	1.60 Transitions® Signature® VII	1.60	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	1.67 Essilor Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-6.00
	1.67 Xperio UV™ Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-6.00
	1.67 Transitions® Signature® VII	1.67	Gray, Brown, Graphite Green	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	1.67 Transitions® XTRActive®	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	1.74 Clear*	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00
1.74 Transitions® Signature® VII*	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00	

NOTES: *Thin & Lite 1.74 lenses come systematic with choice of Crizal® Previncia, Crizal Sapphire UV™ Crizal Avancé UV™ or Crizal Alizé UV™
Please include fitting height when ordering Eyezen+ lenses. Recommended minimum fitting height is 15mm.

EYEZEN+ LENS DESIGNS

Eyezen+ 0	patients age 17 and under
Eyezen+ 1	patients age 18 to 34
Eyezen+ 2	patients age 35 to 44
Eyezen+ 3	patients age 45 to 50**

There are four unique **Eyezen+** designs. Although the recommended age ranges for each design provided below are based on research, as always, the Eyecare Professional is free to use discretion in prescribing **Eyezen+** enhanced single vision lenses.

**It is recommended that you consider progressive lenses for your patients age 45-50.



Better Sight.
Better Life.

