

Eyezen® Start™

Fitting Eyezen® Lenses

Eyezen Start lenses require the following standard measurements:

- Monocular Distance PD
- Fitting Height: Simply measure from center pupil to the lowest point of the lens to acquire fitting height. There is no minimum fitting height

Eyezen Start lenses with embedded Harmful Blue Light* protection are always recommended with Crizal® No-Glare treatment for our most comprehensive eye protection.

The following laser engraving, located on the nasal side of each lens, will identify the Eyezen Start design received.

Eyezen Start



Embedded Harmful Blue Light* Protection.

* Eyezen Start lenses filter at least 20% of Harmful Blue Light, which is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells.

EYEZEN® START™ LENSES AVAILABILITY for Eyecare Professionals

Effective January 2021

MATERIAL	LIGHT MANAGEMENT	INDEX	COLOR	POWER RANGE	BASE CURVE	MAX CYLINDER
PLASTIC	Clear	1.5		-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Essilor Polarized	1.5	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 7.25	-4.00
	Xperio UV™ Polarized	1.5	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 7.25	-4.00
	Transitions® Signature® GEN 8™	1.5	Gray, Brown, Graphite Green	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® Signature® GEN 8™ lenses style colors	1.5	Amber, Amethyst, Emerald, Sapphire	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® XTRActive®	1.5	Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® XTRActive® Style Mirrors	1.5	Blue, Green, Gold, Pink, Red, Silver Shadow	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
TREXA™	Clear	1.53		-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
	Transitions® Signature® GEN 8™	1.53	Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-6.00
AIRWEAR®	Clear	1.59		-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Essilor Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	-6.00
	Xperio UV™ Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	-6.00
	Transitions® Signature® GEN 8™	1.59	Gray, Brown, Graphite Green	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Transitions® Signature® GEN 8™ lenses style colors	1.59	Amber, Amethyst, Emerald, Sapphire	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Transitions® XTRActive®	1.59	Gray, Brown	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
	Transitions® XTRActive® Style Mirrors	1.59	Blue, Green, Gold, Pink, Red, Silver Shadow	-10.00 to +6.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00	-6.00
THIN & LITE®	Clear	1.6		-12.00 to +8.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50	-6.00
	Transitions® Signature® GEN 8™	1.6	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.50	-6.00
	Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Essilor Polarized	1.67	Gray, Brown, Gray-Green	-12.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-4.00
	Xperio UV™ Polarized	1.67	Gray, Brown, Gray-Green	-12.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-4.00
	Transitions® Signature® GEN 8™	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Transitions® Signature® GEN 8™ lenses style colors	1.67	Amber, Amethyst, Emerald, Sapphire	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive®	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive® Style Mirrors	1.67	Blue, Green, Gold, Pink, Red, Silver Shadow	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	-6.00
	Clear**	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00
	Transitions® Signature® GEN 8™**	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00

** Thin & Lite 1.74 lenses come systematic with your choice of Crizal Sapphire® 360° UV™, Crizal® Prevencia®, Crizal® Rock™, Crizal Avancé UV™ or Crizal Alizé UV™.

EYEZEN START LENS DESIGNS

LENS ENGRAVING	
PRODUCT	SYMBOL
Eyezen Start	



Eyezen® +

Fitting Eyezen® + Lenses

Eyezen+ lenses require the following standard measurements:

- Monocular Distance PD
- Fitting Height: Simply measure from center pupil to the lowest point of the lens to acquire fitting height. For **Eyezen+**, minimum fitting height is **15mm**.

Eyezen+ lenses with embedded Harmful Blue Light* protection are also recommended with Crizal® No-Glare treatment for our most comprehensive eye protection.

The following laser engravings located on the nasal side of each lens will identify the **Eyezen** design received.

Eyezen+ 0



Eyezen+ 1



Eyezen+ 2



Eyezen+ 3



Eyezen+ 4



Embedded Harmful Blue Light* Protection.

* Eyezen lenses filter at least 20% of Harmful Blue Light, which is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells.

EYEZEN®+ LENSES AVAILABILITY for Eyecare Professionals

Effective January 2021

MATERIAL	LIGHT MANAGEMENT	INDEX	COLOR	POWER RANGE	BASE CURVE	MAX CYLINDER
PLASTIC	Clear	1.5		-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Essilor Polarized	1.5	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	-6.00
	Xperio UV™ Polarized	1.5	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	-6.00
	Transitions® Signature® GEN 8™	1.5	Gray, Brown, Graphite Green	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® Signature® Style Colors	1.5	Amber, Amethyst, Emerald, Sapphire	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® Vantage®	1.5	Gray	-10.00 to +7.50	2.25, 4.25, 5.25, 8.25	-6.00
	Transitions® XTRActive® Gray*	1.5	Gray	-10.00 to +7.50	2.25, 4.25, 5.25, 6.25, 8.25	-6.00
	Transitions® XTRActive® Brown	1.5	Brown	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® XTRActive® Style Mirrors	1.5	Blue, Green, Gold, Pink, Red, Silver Shadow	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
TREXA™	Clear	1.53		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	Transitions® Signature® GEN 8™	1.53	Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
AIRWEAR®	Clear	1.59		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Essilor Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	-6.00
	Xperio UV™ Polarized	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	-6.00
	Transitions® Signature® GEN 8™	1.59	Gray, Brown, Graphite Green	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® Signature® Style Colors	1.59	Amber, Amethyst, Emerald, Sapphire	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® Vantage®*	1.59	Gray	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® XTRActive®	1.59	Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® XTRActive® Style Mirrors	1.59	Blue, Green, Gold, Pink, Red, Silver Shadow	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Colors™**	1.59		-10.00 to +8.00	2.25, 4.25, 6.50, 7.50, 8.50	-6.00
THIN & LITE®	Clear	1.6		-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	Transitions® Signature® GEN 8™	1.6	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Essilor Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-4.00
	Xperio UV™ Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-4.00
	Transitions® Signature® GEN 8™	1.67	Gray, Brown, Graphite Green	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Transitions® Signature® Style Colors	1.67	Amber, Amethyst, Emerald, Sapphire	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive®	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Transitions® XTRActive® Style Mirrors	1.67	Blue, Green, Gold, Pink, Red, Silver Shadow	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	Clear**	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00
	Transitions® Signature® GEN 8™**	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00

* Transitions Vantage, AIRWEAR® Colors™, 1.50 plastic Transitions® XTRActive® Gray available on Eyezen+ 0-3 only

** Thin & Lite 1.74 lenses come systematic with your choice of Crizal Sapphire® 360° UV™, Crizal® Previncia®, Crizal® Rock™, Crizal Avancé UV™ or Crizal Alizé UV™.

EYEZEN+ LENS DESIGNS

Eyezen+ 0	patients age 17 and under
Eyezen+ 1	patients age 18 to 34
Eyezen+ 2	patients age 35 to 39
Eyezen+ 3	patients age 40 to 44
Eyezen+ 4	patients age 45 to 50 ¹

There are five unique **Eyezen+** designs. Although the recommended age ranges for each design provided below are based on research, as always, the Eyecare Professional is free to use discretion in prescribing **Eyezen+** enhanced single vision lenses.

¹ Eyezen+ lenses are not adapted to patients that have been prescribed progressive lenses.