

HERITAGE relief



Benefits

An innovative single vision lens design for patients who are not quite presbyopic. Generally from ages 25 to 40 but younger patients with heavy digital device use will experience benefits too!

It prevents the patient's eyes from becoming tired and strained if they spend much of the day focusing on the same distance such as using a computer.

Supports near vision focusing by having a small amount of addition that allows the eyes to use less of their accommodation and therefore allows the wearer to continue doing near vision tasks for longer without tiring their eyes.

Built In Technology

FULL BACK
SIDE

EYE
VIEW

Corridor
Additions

14mm
+0.66D and +1.00D

Performance

Distance
Intermediate
Reading
Comfort

