

2023

AVAILABILITY CHART



EFFECTIVE MARCH 2023

MIRROR OPTIONS:



Blue Mirror

High Efficiency Polarized Filter
14% Light Transmission
91% Blue Light Filtered



Silver Mirror

High Efficiency Polarized Filter
14% Light Transmission
90% Blue Light Filtered



Gold Mirror

High Efficiency Polarized Filter
17% Light Transmission
94% Blue Light Filtered



Orange Mirror

High Efficiency Polarized Filter
17% Light Transmission
94% Blue Light Filtered



Green Mirror

High Efficiency Polarized Filter
14% Light Transmission
89% Blue Light Filtered



Ice Blue Mirror

High Efficiency Polarized Filter
12% Light Transmission
95% Blue Light Filtered



Red Mirror

High Efficiency Polarized Filter
20% Light Transmission
96% Blue Light Filtered



Pink Mirror

High Efficiency Polarized Filter
17% Light Transmission
90% Blue Light Filtered



Purple Mirror

High Efficiency Polarized Filter
18% Light Transmission
96% Blue Light Filtered



Black Mirror

High Efficiency Polarized Filter
10% Light Transmission
92% Blue Light Filtered



Chrome Mirror

High Efficiency Polarized Filter
12% Light Transmission
89% Blue Light Filtered



Rose Gold

High Efficiency Polarized Filter
16% Light Transmission
97% Blue Light Filtered



Champagne

High Efficiency Polarized Filter
14% Light Transmission
87% Blue Light Filtered

MATTE MIRROR OPTIONS



Sunny Mist

High Efficiency Polarized Filter
14% Light Transmission
97% Blue Light Filtered



Sea Mist

High Efficiency Polarized Filter
15% Light Transmission
88% Blue Light Filtered



Charcoal Mist

High Efficiency Polarized Filter
15% Light Transmission
88% Blue Light Filtered

SOLID COLOR OPTIONS:



Gray-3

High Efficiency Polarized Filter
15% Light Transmission
86% Blue Light Filtered*
Also available non-polarized



Brown-3

High Efficiency Polarized Filter
14% Light Transmission
96% Blue Light Filtered
Also available non-polarized



Gray-Green

High Efficiency Polarized Filter
13% Light Transmission
89% Blue Light Filtered
Also available non-polarized



Gray-1

High Efficiency Polarized Filter
35% Light Transmission
67% Blue Light Filtered



Brown-1

High Efficiency Polarized Filter
27% Light Transmission
85% Blue Light Filtered



Copper

High Efficiency Polarized Filter
18% Light Transmission
94% Blue Light Filtered



Ruby

High Efficiency Polarized Filter
15% Light Transmission
91% Blue Light Filtered



Grass

High Efficiency Polarized Filter
35% Light Transmission
79% Blue Light Filtered



Amber/Melanin

High Efficiency Polarized Filter
24% Light Transmission
92% Blue Light Filtered



Sunflower

Low Efficiency Polarized Filter
87% Light Transmission
62% Blue Light Filtered



Violet

Mid Efficiency Polarized Filter
41% Light Transmission
26% Blue Light Filtered



Rose

Low Efficiency Polarized Filter
44% Light Transmission
50% Blue Light Filtered



Ocean

Mid Efficiency Polarized Filter
50% Light Transmission
34% Blue Light Filtered



Quartz

Non-Polarized
41% Light Transmission
65% Blue Light Filtered

GRADIENT OPTIONS:



Gray Gradient

High Efficiency Polarized Filter
5% at darkest, 35% lightest
82% Blue Light Filtered



Brown Gradient

High Efficiency Polarized Filter
5% at darkest, 27% lightest
95% Blue Light Filtered



Ruby Gradient

High Efficiency Polarized Filter
10% at darkest, 31% lightest
82% Blue Light Filtered



Blue Gradient

High Efficiency Polarized Filter
17% at darkest, 31% lightest
73% Blue Light Filtered



Green Gradient

High Efficiency Polarized Filter
7.5% at darkest, 31% lightest
83% Blue Light Filtered



www.kbco.net

AVAILABILITY CHART 2023



Plastic 1.55 HC/UC	Blank	Base	Add													
Single Vision	75	.50,2,4,6,8		S	S	M	M	M	M	M	M	M	M	M	M	M
Single Vision	85	4,6,8		M	M	M	M	M	M			M				
FT-28 Bifocal	75	.50,2,4,6,8	1.25D-3.25D	S	S	M	M	M	M	M	M	M	M	M	M	M
FT-35 Bifocal	75	4,6,8	1.25D-3.25D	M	M	M	M	M	M	M	M	M	M	M	M	M
7X28 Trifocal	75	4,6,8	1.75D-3.25D	S	S	M	M	M	M	M	M	M	M	M	M	M
8X35 Trifocal	75	4,6,8	1.75D-3.25D	M	M	M	M	M	M	M	M	M	M	M	M	M
iRx Pro™ PAL MFH 18mm	76/80	2,4,6,7	1.00D-3.50D	M	M	M	M	M	M	M	M	M	M	M	M	M
Polycarbonate HC	Blank	Base	Add													
Plano Rectangle 1.5mm thick	60x70	6,8		S	S	S										
Plano Rectangle 2.2 mm thick	75	6,8		S	S	S										
Plano Round 1.5mm thick	71	6		S	S	S										
Single Vision (DS)	81	1,3,4,5,7,8		S	S	S										
Single Vision (DS)	76	1,2,3,4,5,6,7,8		S	S	S	S	S	S	S						
Single Vision	82	6,7,8		S	S	S										
FT-28 Bifocal	76	2,4,6,8	1.00D-3.00D	S	S	S										
FT-35 Bifocal	76	2,4,6,8	1.50D-3.00D	S	S											
Xplore™ PAL MFH 17mm	81	1,3,4,5,6	1.00D-3.00D	S		S										
High-Index 1.60 HC/UC	Blank	Base	Add													
Single Vision HC (DS)	75	1,2,3,4,5,6,7,8,9		S	S											
Single Vision UC (DS)	81	1,2,3,4,5,6,8		S	S	S										
High-Index 1.67 HC/UC	Blank	Base	Add													
Single Vision	75	.50,2,4,6,8					S	S	S	S	S		S	S	S	S
Single Vision (DS)	81	1.25,1.75,2,3,4,5,6,6.75,8		S	S	S										
Single Vision Aspheric	75	.50,2,4,6,8		S	S											
FT-28 Bifocal	75	.50,2,4,6,8	1.25D-4.00D	S	S	M	M	M	M	M	M		M	M	M	M
FT-35 Bifocal	75	4,6,8	1.25D-4.00D	M	M	M	M	M	M	M	M		M	M	M	M
7X28 Trifocal	75	4,6,8	2.00D-4.00D	S	S	M	M	M	M	M	M		M	M	M	M
8X35 Trifocal	75	4,6,8	2.00D-4.00D	M	M	M	M	M	M	M	M		M	M	M	M
iRx Pro™ PAL MFH 18mm	75	2,4,5,7	1.00D-3.00D			M	M	M	M	M	M		M	M	M	M

KBCO NON-POLARIZED LENSES

Polycarbonate Airwear Colors™ HC Non-Polar	Blank	Base	Add	Gray-3	Brown-3	Gray-Green	Quartz	Blue
Single Vision	76	2,4,6,7,8		S	S	S	S	S

XPERIO UV™ LENSES

Go to XperioUV.com for complete listing of colors & lens styles.

Plastic 1.50	Gray-3	Brown-3	Gray-Green	Gray-1	Brown-1	Copper	Ruby
Varilux® Stylistic®	●	●	●	●	●		
Varilux Physio DRx™	●	●	●	●	●		
Varilux Physio Short DRx™	●	●	●	●	●		
Varilux Comfort DRx™	●	●	●	●	●		
Varilux Comfort Short DRx™	●	●	●	●	●		
Single Vision	●	●	●	●	●		
FT-28	●	●	●				
FT-35	●	●					
7x28	●	●					
8x35	●	●					
Irx Pro™	●	●					
Xplore™	●	●	●				
Plastic 1.55	Gray-3	Brown-3	Gray-Green	Gray-1	Brown-1	Copper	Ruby
Single Vision	●	●					
FT-28	●	●					
FT-35	●	●					
7x28	●	●					
8x35	●	●					
Irx Pro™	●	●					

Polycarbonate	Gray-3	Brown-3	Gray-Green	Gray-1	Brown-1	Copper	Ruby
Varilux® Stylistic®	●	●	●	●	●	●	●
Varilux Physio DRx™	●	●	●	●	●	●	●
Varilux Physio Short DRx™	●	●	●	●	●	●	●
Varilux Comfort DRx™	●	●	●	●	●	●	●
Varilux Comfort Short DRx™	●	●	●	●	●	●	●
Single Vision	●	●	●	●	●	●	●
FT-28	●	●	●				
Xplore™	●	●	●				
High Index 1.60	Gray-3	Brown-3	Gray-Green	Gray-1	Brown-1	Copper	Ruby
Single Vision	●	●	●				
High Index 1.67	Gray-3	Brown-3	Gray-Green	Gray-1	Brown-1	Copper	Ruby
Varilux® Stylistic®	●	●	●				
Varilux Physio DRx™	●	●	●				
Varilux Physio Short DRx™	●	●	●				
Varilux Comfort DRx™	●	●	●				
Varilux Comfort Short DRx™	●	●	●				
FT-28	●	●					

(DS) extended base curve range optimized for digital surfacing
Please note that add/base curve/Rx ranges may vary slightly based on color.

†Gray-Green only available in 6 and 8 base.

HC - Hard Coated
UC - Uncoated

For more information,
go to www.kbco.net



S - In stock
M - Made-to-order