

Help your patients bring fashion to function with the new KBco Matte Mirror Polarized Lenses

AVAILABLE NOW IN 3 COLORS



Charcoal Mist

Base Color: Gray
Light Transmittance: 15%
Blue Light Filter: 88%



Sunny Mist

Base Color: Brown
Light Transmittance: 14%
Blue Light Filter: 97%



Sea Mist


Base Color: Gray
Light Transmittance: 15%
Blue Light Filter: 88%

Available in **Semi-Finished Single Vision** lenses (76 mm, base curves 4, 6, 8) and **Finished Plano** (76x55 mm, base curve 6)

ADVANCED POLARIZED TECHNOLOGY:



Blocks UV light coming through the front of the lens, protecting your patients' eyes



Reduces painful and blinding glare



Improves color perception, making colors crisper and more vibrant than tinted sunglass lenses

Recommend Matte Mirror lenses to your patients today and help them protect their eyes from UV with a stylish and trendy solution.



Optimal Base Curve Range for Digital Surfacing

POLYCARBONATE

| Base Color / Mirror Color | Base Curve Range | Power Range | Diameter |
|---|------------------|----------------|----------|
| Brown 3 Base <ul style="list-style-type: none">• Sunny Mist | 4.25, 6.25, 8.25 | -4.00 to +6.00 | 76mm |
| Gray 3 Base <ul style="list-style-type: none">• Charcoal Mist• Sea Mist | 4.25, 6.25, 8.25 | -4.00 to +6.00 | 76mm |

Faster delivery time

Process in-house, eliminating outsourcing

Advanced polarized technology

Provides comprehensive UV protection

Wide power range

Polycarbonate: -4.00D to +6.00D

Optimized base lens color-matching

Stylish Sunny Mist Mirror matched to a Brown Base, while Charcoal and Sea Mist matched to a Cool Gray Base



To learn more, visit www.KBco.net