

Technology Driven



Vision First Design Technology™ is a proprietary technology used to produce **Kodak Unique Lenses**.

Vision First Design Technology provides:

- ▶ A broad field of view in a clear distance area
- ▶ Smooth gradation of power across the surface of the lens to ease patient adaptation
- ▶ Gentle binocular balance for quick, clear object recognition
- ▶ Visual quality often compared to a single vision lens



i-Sync™ Technology elevates the level of optical performance by maximizing clarity in the peripheral areas of the lens and reducing distortion to the lens edge.

i-Sync Technology provides:

- ▶ More consistent optical performance over the range of prescription power
- ▶ Wider near-viewing areas for hyperopes
- ▶ Improved distance area for myopes
- ▶ Improved image quality in principal viewing areas
- ▶ Flatter base curve capabilities

KODAK Unique Lenses.
Because We Are All Unique.

Kodak LENS ids INTELLIGENT DISPENSING SOFTWARE

Kodak Lens Intelligent Dispensing Software

Kodak Lens Intelligent Dispensing Software (IDS) is composed of three separate sections: **FRAME**, **LENS** and **MEASURE**:

FRAME: Compares and shares photos of patient frame options to help their selection.

LENS: Assists the eyecare practitioner in educating the patient about the visual benefits of **Kodak** Progressive Lenses with virtual simulations and augmented reality technology. Comparative screens demonstrate high index, polarized, photochromic and AR coating options.

MEASURE: Easily converts the iPad into a measuring device – providing monocular PD measurements for far and near vision, monocular fitting height measurements, pantoscopic tilt, back vertex distance, face form angle, A, B and DBL measurements.

Check out the **Kodak** Lens IDS informational video on the Signet Armorlite YouTube channel, www.youtube.com/SignetArmorliteUS.



WE ARE ALL UNIQUE.



See the **Colors of Life**®

Kodak
Unique Progressive Lens

Kodak
Unique HD Progressive Lens

For over seventy years, Signet Armorlite has been a leader in the optical industry with quality products and innovative technologies.

Signet Armorlite, Inc.
800.950.5367
www.signetarmorlite.com



Kodak Unique HD Lens customization includes the following additional features:

Easy Dispensing

A digitally-created full backside design, **Kodak Unique™** Progressive Lens delivers increased prescription accuracy with smooth gradation of power across the lens surface for clear, comfortable vision. With its ease in dispensing and high adaptation rate, **Kodak Unique Lens** is a practice favorite.


The eyecare practitioner may request a specific corridor length or allow sophisticated software technology to select it based on the frame size, monocular PD and fitting height for each individual prescription.

See Life in High Definition

Kodak Unique HD™ Lens includes all the features of **Kodak Unique Lens** with added customization.

The HD stands for 'High-Definition' and fulfills this description by including additional measurements in the design to provide a higher level of customization to match the patient's unique viewing needs.

[See additional measurements →](#)



Prescription Compensation

The patient's **Back Vertex Distance (BVD)**, **Pantoscopic Tilt** and **Wrap Angle** are included to provide a customized lens optimized for the individual patient's prescription and frame selection.

Variable Inset

Relying on the **Near Reading Distance**, **Lens Power**, **Back Vertex Distance (BVD)** and **Monocular PD (MPD)**, **Kodak Unique HD Lens** delivers a near viewing experience tailored to the individual patient's needs.


Patient's Viewing Experience

- Wide fields of view offer picture-perfect vision at all distances
- Advanced technologies allow greater visual comfort
- Reduced distortion in the design for smooth, natural viewing and comfortable wear
- Greater color contrast for a more dynamic viewing experience
- Thin lenses resulting in comfortable and more cosmetically-appealing eyeglasses
- Optional customization

*We are all unique.
Our viewing experience
should be as well.*

ECP's Dispensing Experience

- Outstanding patient adaptation
- A trusted brand name
- Easy to dispense and customize
- Advanced technologies at a sensible price
- A wide variety of material options



Feature	KODAK Unique Lens	KODAK Unique HD Lens
• Corridor Lengths	6	6
• Material Availability	60+	60+
• Full Backside Progressive Design	✓	✓
• Vision First Design Technology™	✓	✓
• i-Sync™ Technology	✓	✓
• Variable Decentration	✓	✓
• Wide Distance Vision	✓	✓
• Prescribed Prism	✓	✓
• Variable Inset	✓	✓
• Prescription Compensation		✓