

Regular Lens


KODAK CleAR+ Lens



Engineered for Performance

Both sides of KODAK CleAR+ Lens include:

- Protective layers that provide durability and scratch resistance
- NEW!** • Anti-static properties that resist clinging dust
- Several microscopic layers of minerals that effectively reduce reflections
- Hydrophobic top coat seals the lens to help resist dirt and smudges



The brand is KODAK.
The quality you and
your patients trust.

KODAK CleAR+ Lens carries one of the most respected brands in the world. KODAK CleAR+ Lens was designed and tested by Signet Armorlite, a company with over 66 years of leadership and innovation in ophthalmic lens technology.

See the *Colors* of Life

Kodak LENS



Signet Armorlite
5803 Newton Drive, Suite A
Carlsbad, CA 92008
800-759-0075
www.signetarmorlite.com



Photographic demonstrations simulated.
The Kodak trademark and trade dress are used under license from Kodak by Signet Armorlite, Inc. CleAR is a trademark of Signet Armorlite, Inc. © 2015 Signet Armorlite, Inc.

PN 358-818 JS 2/2015

Kodak CleAR+ Lens



Help Your Patients See Life More Clearly.

The CleAR Choice for Your Patients

Kodak CleAR+ Lens

Easy-to-Clean

Resists oily smudges, dust and dirt, and cleans with ease – making clean lenses easy to maintain.



Scratch Resistant

Avoid irritating scratches so lenses last longer and look better over the life of the prescription.



Reduces Nighttime Reflections

Improves ability to see objects in low light for safer, more comfortable night driving.



Eliminates Irritating Glare

Reduces eyestrain and fatigue from artificial lighting and computer screens.



Help Patients See Clearer and Crisper

With ordinary lenses, approximately 10% of light is reflected away and doesn't reach the eyes. More reflections and less light mean poorer vision. With KODAK CleAR+ Lens, approximately 99% of available light reaches the eyes providing clear and crisp vision.



Regular Lens



KODAK CleAR+ Lens

Easy to Dispense

KODAK CleAR+ Lens allows eyecare professionals to offer patients a recognizable brand known for trust, respect and loyalty.