



Kodak CleAR

Anti-Reflective Coating

Regular Lens

Kodak CleAR
Anti-Reflective Coating

Help Patients See the World More Clearly

Easy to Clean – Resists oily smudges, dust and dirt, and cleans with ease making them easy to maintain.

Scratch Resistance – Protects lenses from irritating scratches so they last longer and look better over the life of the prescription.

Reduces Night-time Reflections – Improves ability to see objects in low light for safer, more comfortable night driving.



Regular Lenses



KODAK CleAR Lenses


Eliminates Irritating Glare – Reduces eyestrain and fatigue from artificial lighting and computer screens.

Easy to Dispense – Kodak is one of the oldest and most respected global brands. It portrays a level of quality that everyone has come to trust.



Engineered for Performance

Both sides of KODAK CleAR™ Lenses include: protective layers that provide durability and scratch resistance; several microscopic layers of minerals that effectively reduce reflections over a broad band of wavelengths. A super-hydrophobic top coat seals the lens to help resist dirt and smudges and makes it easier to clean.



Help Patients See Clearer and Crisper

With ordinary lenses, approximately 10% of light is reflected away and doesn't reach the eyes. More reflections and less light mean poorer vision. With KODAK CleAR Lenses, approximately 99% of available light reaches the eyes providing clearer and crisper vision.



1001 Armorlite Dr., San Marcos, CA 92069 • 800-950-5367 • www.signetarmorlite.com



Photographic demonstrations simulated. Kodak and the Kodak trade dress are trademarks of Kodak, used under license by Signet Armorlite, Inc. CleAR is a trademark of Signet Armorlite, Inc. ©2010 Signet Armorlite, Inc. PN 359-089 01/10