

CLEANING INSTRUCTIONS:

1. Rinse lenses under warm running water to remove any abrasive particles.
2. Clean with mild liquid soap to remove any grease or oily build-up on the lenses.
3. Dry with a soft cotton towel or optical lens cleaning cloth.

* For quick touch-ups, use optical cleaning spray and dry with an optical lens cleaning cloth.



...AND HELP OTHERS SEE YOU MORE CLEARLY TOO!

Your Eyes Will Shine Through

Eliminate distracting reflections and glare visible to others.

Lenses Will Almost Disappear

Less glare makes your lenses appear thinner and more transparent - like no lenses at all!



Regular Lenses



KODAK Clean'N'Clear Lenses

Performance and value you can count on.

KODAK Clean'N'Clear Lenses provide the quality and performance you've come to expect from a Kodak-branded product.

Ask your eyecare professional today if these lenses are right for you.

Signet Armorlite, Inc.
800.950.5367
www.signetarmorlite.com



Photographic demonstrations simulated. Kodak and the Kodak trade dress are trademarks of Kodak, used under license by Signet Armorlite, Inc. Clear is a trademark of Signet Armorlite, Inc. ©2010 Signet Armorlite, Inc.

PN 480-596

PROTECT YOUR LENSES FROM

SCRATCHES

FINGERPRINTS & DIRT



Eliminate Distracting Glare and Reflections.
Protect Your Eyewear.
Clean with Ease.

Kodak Clean'N'Clear
Anti-Reflective Coating

CLEAR VISION...

More Clarity and Contrast

Greatly improve your vision by allowing up to 99% of available light to reach your eyes. More light means sharper vision.

Reduce Eyestrain

Eliminate irritating glare from artificial lighting and computer screens, providing less stress on your eyes.

Safer, Comfortable Night Driving

Significantly reduce night-time reflections, and improve your ability to see objects in low light.



Regular Lenses



KODAK Clean'N'Clear Lenses

...WITH EASY-TO-CLEAN LENSES



Clean with Ease

New optical technology helps resist oily smudges and dirt, and makes your lenses easier to clean and maintain. Anti-static properties repel clingy dust and particles.


Scratch Resistance

Protect your lenses from irritating scratches so they will last longer and look better over the life of your prescription.

...AND INNOVATIVE TECHNOLOGY

Your lens includes:

- A scratch-resistant layer
- Multiple microscopic layers to reduce reflections over a broad band of wavelengths
- A super-hydrophobic top layer that seals the lens and makes it easier to clean
- Anti-static properties



Kodak Clean'N'Clear Anti-Reflective Coating

Your eyewear will last longer and look better if you follow these simple lens care and cleaning instructions.

SPECIAL CARE INSTRUCTIONS:

- Always keep your eyewear in a protective case when not wearing.
- Lay your eyewear face-up rather than face-down.
- Never subject your eyewear to high heat; it may cause your frame to bend.
- Never rub your lenses when they are dry.

