




KODAK Clean'N'Clear Lenses are engineered to enhance the lens for greater all-around performance.

It includes:

- » A scratch-resistant layer.
- » Multiple microscopic layers to reduce reflections over a broad band of wavelengths.
- » A super-hydrophobic top layer that seals the lens and makes it easier to clean.



KODAK Clean'N'Clear Lenses carry one of the most respected brands in the world. It is designed to provide the quality and durability your patients demand, with a brand name they know and trust.

Kodak LENS

See the *Colors* of Life



5803 Newton Drive, Suite A
Carlsbad, CA 92008
800-759-0075
www.signetarmorlite.com



The Kodak trademark and trade dress are used under license from Kodak by Signet Armorlite, Inc. Clear is a trademark of Signet Armorlite, Inc.

© 2014 Signet Armorlite, Inc.

PN 359-907 Rev 12/14

Kodak
Clean'N'Clear Lens



The Difference is **Clear**



Kodak

Clean'N'Clear Lens

KODAK Clean'N'Clear™ Anti-Reflective Coating enhances the performance of prescription lenses by reducing glare for clearer vision, while improving the cosmetic appeal of eyewear. Your patients will not only see better—they will look better, too.

The difference is strikingly clear.



With an untreated lens, approximately 10% of light is reflected away from the eyes, which makes it difficult to see clearly.

With KODAK Clean'N'Clear Lenses, glare is reduced, enabling approximately 99% of available light to reach the eyes, for clear, crisp vision.

KODAK Clean'N'Clear Lenses provide patients with clean, clear protection:



Reduces eye strain and fatigue caused by irritating glare from harsh artificial lighting, tablets and computer screens.




Significantly reduces nighttime reflections, and improves the wearer's ability to see objects in low light.



Oleophobic technology makes lenses easier to care for: Oily smudges and grime wipe away with ease. Anti-static properties resist dust and dirt.



Protects lenses from irritating scratches so they will last longer and look better.



Eliminates unsightly glare on lenses for a more appealing aesthetic.



Provides water-repellent performance
A super-slick hydrophobic top coat seals the lens and repels water.

Prescribe with confidence.

KODAK Clean'N'Clear Lenses are part of a full family of professional-grade coatings that deliver advanced optical technology, with the quality and performance that is synonymous with the trusted Kodak brand.

For more information about our anti-reflective coatings, please visit us at www.signetarmorlite.com