

# PHOTOVIEWS™

PHOTOCHROMIC LENSES

## FADE TIME

Faster fade times are more convenient, changing tints in different light conditions.  
80% Dynamic Range, minutes

**PHOTOVIEWS**  
**5 MINUTES**



**POPULAR IN-MASS  
PHOTOCHROMIC**  
**44 MINUTES**




## UV PROTECTION

More complete direct UV protection promotes better eye health.  
% UVA and UVB frontside light transmission protection

**PHOTOVIEWS**

**100%\*\*  
PROTECTION**



**POPULAR IN-MASS  
PHOTOCHROMIC**

**96%\*\*  
PROTECTION**

**FASTER FADE TIMES.**  
**100% DIRECT UV PROTECTION.**

PhotoViews photochromic lenses rapidly keep up with changing light conditions while providing full direct UV protection.

Introduce more patients to photochromic lenses with the performance, protection and price of PhotoViews lenses.

**FASTER TINT ADJUSTMENT.**  
**FULL DIRECT UV PROTECTION.**  
**SMART PRICING.**

## PHOTOVIEWS AVAILABILITY

### STANDARD RESIN

Index of Refraction: 1.497  
 Abbe Value: 60  
 Specific Gravity (gm/cm<sup>3</sup>): 1.27

### POLYCARBONATE

Index of Refraction: 1.586  
 Abbe Value: 30  
 Specific Gravity (gm/cm<sup>3</sup>): 1.20  
 Recommended for drill mount frames.

## SIGNETEK AVAILABILITY

Lens Design	Colors	STANDARD RESIN		POLYCARBONATE	
		Sphere Range	Maximum Cutout	Sphere Range	Maximum Cutout
KODAK Unique™ HD	gray/brown	+4.00 to -8.00	71mm	+5.00 to -10.00	69/74mm
KODAK Unique	gray/brown	+4.00 to -8.00	71mm	+5.00 to -10.00	69/74mm
KODAK Precise® Digital*	gray/brown	+5.00 to -8.00	71mm	+6.00 to -10.00	73mm
KODAK Precise Short Digital*	gray/brown	+5.00 to -8.00	71mm	+6.00 to -10.00	73mm
KODAK Precise PB	gray/brown	+4.00 to -8.00	71mm	+5.00 to -10.00	69/74mm
KODAK Precise Short PB	gray/brown	+4.00 to -8.00	71mm	+5.00 to -10.00	69/74mm
DiracTek™	gray/brown	+4.00 to -8.00	71mm	+5.00 to -10.00	69/74mm
DiracTek Short	gray/brown	+4.00 to -8.00	71mm	+5.00 to -10.00	69/74mm
KODAK Digital Single Vision	gray/brown	+6.00 to -8.00	71mm	+6.50 to -10.00	69/74mm
Signetek™ Single Vision	gray/brown	+6.00 to -8.00	71mm	+6.50 to -10.00	69/74mm

\*KODAK Precise Digital and Precise Short Digital PhotoViews Poly only available in gray.

## SEMI-FINISHED AVAILABILITY

Lens Design	Colors	STANDARD RESIN			POLYCARBONATE		
		Base Curves	Diameter	Colors	Base Curves	Diameter	
KODAK Precise	gray/brown	1.50,3.25,5.00,7.00	76mm	gray only	1.50,3.25,5.50,7.00	78mm	
KODAK Precise Short	gray/brown	2.00,4.50,5.50,7.50	76mm	gray only	2.00,4.50,5.50,7.50	78mm	
Navigator®	gray/brown	2.00,4.50,5.75,7.75	76mm				
Navigator Short®	gray/brown	1.25,3.25,5.50,7.50	76mm				
Single Vision	gray/brown	0.50,2.25,4.25,5.25,6.25,8.25	76mm	gray/brown	1.50,3.25,4.25,4.75,5.50,7.00,7.50*,8.50*	74*/79mm	
Flat-Top 28	gray/brown	0.50,2.25,4.25,5.25,6.25,8.25	76mm				
Trifocal 7x28	gray only	4.50,6.25,8.75	74mm				
Finished Single Vision	gray/brown	65mm / +0.25 to +4.00 / Cylinder: 0 to -2.00 70mm / Plano to -4.00 / Cylinder: 0 to -2.00					

PhotoViews lenses are available exclusively from Signet Armorlite, a leading innovator in photochromic lenses. Recognizing the need for an array of photochromic options, Signet Armorlite has a long history of supporting the growth of adaptable lenses across the industry.



Signet Armorlite, Inc.

www.signetarmorlite.com