



ELECTRONIC PROGRESSIVE IDENTIFIER CATALOG

PDF Generated: 6/28/2020 4:01 AM

How to Use This Catalog

Each lens in this catalog is labeled with a sequential number (Item Number), shown in the right hand corner of the lens entry. In each index, symbols and lens characteristics are listed along with their relevant lens item numbers. To find a lens in the catalog, locate the item number(s) in one of the indexes below and then simply flip to the section containing the relevant item numbers.

Index by Symbol

Symbol	Item #	Symbol	Item #	Symbol	Item #
*	70, 79, 82, 85-87, 90, 92-94, 96		106		417-421, 426-432
—	172		282, 284		411
==	178-179		275		486
^	39-40, 51-59, 271		492		408
	278		105		282, 298, 303, 307, 311, 314, 391-392
	490		123-124		300, 306, 312, 319
○	1-2, 5-8, 10-11, 15, 23-24, 45, 47, 63-64, 73-74, 76-77, 98-110, 113-118, 120-121, 123-124, 126-158, 160-219, 222-232, 240-252, 257-263, 269-273, 275-281, 285-293, 295-296, 299, 301, 304-305, 309, 313, 315, 317-318, 320-322, 328, 331-366, 368-376, 388-390, 400, 402-404, 407-411, 413-415, 433-435, 464, 469, 471, 473-486, 488		452-453		325
			430		367
			416-417		172, 187, 192
			180		160, 171, 182, 188-190
			378-387	+	133, 436-453, 468, 470
			177		394, 396-397
			444, 446		279, 398
			439-440		253-256
			191		377, 399
			135		115, 145-147
			192		134
			193-195		27
			120		27, 159, 402, 472
			121, 169-170		20-22, 25-26, 28-35, 60, 65-72, 75, 78-96
			160, 171-172		41, 457
			161-162, 164-168, 178-180		464
			181-183		274, 283
			173-175, 177, 184-186		322-326
			187		11-14
			163		489
			419		418
			176		220-221, 439-440, 445
			304		198, 248
			116-117, 151-152, 154		216-217
			118, 156-158		378-380, 382, 393, 404
			472		331-333
			456		337-338
			493		339-340
			409		
			416, 422		
			327		
			407		

Symbol	Item #	Symbol	Item #	Symbol	Item #
2N	31	CC	369	EV	19, 48, 238
2S50	334-336	CD	259-260	EW	49-50
2S59	341-344	CE	18	E <u>W</u>	167, 185
2S60	345-348	CH	65, 370	EWP	230
2S67	349-352	CHD	220	EXC	208
3	240-247, 249-263	CHS	221	F	10, 17-19, 39-42, 44, 46, 48, 50-59, 61-62, 195, 251-258, 454, 456-457
3	20, 327, 353, 412	CHV	40		198, 420
3V60	355-357	C _J	438	FF	426
3V	32	CL	393	FFC	421
4D	194	CN	217	FFO	432
5	391-392	CP	16, 66	FFW	168, 186
6	47, 113-115, 118, 121, 123-124, 127-128, 130-131, 136-137, 141, 143-147, 154-158, 161-165, 169-170, 173-181, 183, 387, 393-394, 397-398, 448, 486	CS	371	F <u>W</u>	360
		C _U	436	G4	248
		CU	43, 454	G	212-213
		C _V	437	H53	331-357
7	9, 383-386	CV	67	H60	
15	416	CW	372	H	5, 36-37, 109, 201-202, 255-256, 383, 459
16	368-376, 422, 465-467, 476-485	CX	394	h	282, 284
50	21, 242, 245, 260, 262, 414	D3	132-133	I2	79
60	22, 26, 28-34, 43, 65-72, 78-96, 98-106, 358-367, 377-382, 388-390, 399, 403-404, 410, 412-413, 415	D	68, 113-114, 130-131, 202	I	78, 251, 323
		DS	69	/	26
		DEM	72	IB	235
67	75, 116-117, 126, 129, 132-134, 138-139, 142, 150-153	DLS	216-217	III6	112
		DNA	204	IS	277
100	85	DSI	70	K ₊	443-444
200	86	E	108, 222, 233, 427	K	447
360°	354	e	231	KB	269-270, 273
400	87	EA	42, 223	Kc	441
A ₃	422	EBA	224, 232	K _{DRO}	449
A	36, 125-126, 129, 215, 240-242, 465	EDR	225, 234	K _{DSII}	442, 450
a	33-34	EH	75	Ks	445-446
A	30	EK	215, 218-219	K _U	448
ABC	105-106	EL	206	L3v	82
AC	38-39	ELC	205	L	252-257, 397
AF	216	EM	44	LD	80-81
AO+	6	EN	239	LX	377
AO	3	EOH	236	M	49, 200, 204-208, 210-214, 258, 385, 447
AOB	7	EP	227	MK	83
AOT	8	EPC	226	MKI	84
AS	243-245, 322, 420	ER	210	N3	320-321
B	21-22, 379, 389	ES	46, 207, 209, 428	N4	392
B	15, 64	ESP	228	N8A	391
bA	1	ET	214	N	361, 451
BIQ	285	EUS	229	Nb	316-317
BPS	196			nc	310-312, 314, 319
C	45, 220-221, 253-254, 400				
c	328				
CA	368				
CAO2	62				
CAOC	61				

Symbol	Item #	Symbol	Item #	Symbol	Item #
NE	148	QS	366	V	285, 373
<i>nm</i>	308-309	QU	365	VCP	374
no symbol	2, 45, 107-110, 120, 125, 159, 222-232	<i>R</i>	377-378, 381-388, 391-398	Ve	480
NS	149, 293	RD	27	VEL	475
NX	150	RNA	211	VG	458
O	324, 453	ROL	358-366	VI	479
ODC	469	S	28, 63, 107, 114, 129, 131, 159, 183, 237, 326, 387, 405-406, 412, 414, 433, 435, 453	Vo1	482
OFF	359		17-19, 38, 41, 455-458, 460-461	Vo5	483
ONM	102	SA	98	Vm	481
OOM	103	SB	99	Vp	484
OPM	104	SC	100, 470	Vw	485
OS	85-87	SD	101	VW	375
p2	88	SE	17	VX	203
P2	89	SH	53	W	241, 244, 246-247, 249, 303, 311-312, 405, 415
P	12, 97, 135, 148-149, 191, 204-211, 214, 219, 293, 329-330, 380, 390, 396	SIG	367		406, 440, 467
<i>p</i>	3-4	SIN	93	<u>W</u>	166, 184
PC	401, 463	SL	413	WG	57
PE	2	SP3	252, 257, 381	WO	58
<u>PE</u>	10	SP	261-262, 276, 306-307	WR	59
PEM	71	SPL	94	WR	139
PF	363	SPP	305, 321	X	11, 13-14, 201
PH	51, 362	SPU	95	X/H	384
PL	90	SSB	96	XL	460
PM	287, 330	<i>STY</i>	192	XLG	461
PN	364	SW	447	X/M	386
PP	320	T	8, 29, 200, 250, 263	XP	376
PRO	4	TLX	291	XPV	212
PS	92, 286, 329	TTX	434	XS	76
PU	91	U	466	XSV	213
PW	478	UH	54	Y	241, 244, 246-247, 263, 289, 491
PX	399, 476	UHC	55	Z	110
Q	468	UHV	56		

Index by Lens Name

Lens Name	Company Name	Item #	Lens Name	Company Name	Item #
Acclaim OEM PAL #1 and #2**	Walman - Soderberg	486	Channel 14 Plastic II	Landon Lens	279
Accolade	Essilor of America	123	Channel 14 Poly II	Landon Lens	280
Accolade Freedom	Essilor of America	124	Channel 14 Variations Gray	Landon Lens	281
Acuity 9/11/13	PFO Global	322	& B		
Adage	Younger Optics	490	Channel 14 Variations Poly	Landon Lens	282
Adaptar	Essilor Canada	111	Pho		
Adaptar	Essilor of America	125	Clarlet Business	Carl Zeiss Vision	15
Adaptar Digital	Essilor of America	126	ClassicLife Free XS. M or L	Rodenstock Canada	377
Adaptar DS	Essilor Canada	112	CLEARi Absolute-W 12	CLEARi	98
Adaptar DS (Canada)	Essilor of America	127	CLEARi DED 14	CLEARi	99
Adaptar DS Short (Canada)	Essilor of America	128	CLEARi DED 16	CLEARi	100
Adaptar Short Digital	Essilor of America	129	CLEARi DED 19	CLEARi	101
Advanced HD	US Optical	465	CLEARi OFFICE BOOK	CLEARi	102
AF2	Seiko Optical Products	400	CLEARi OFFICE	CLEARi	103
Amplitude	Hoya Vision Care - Dallas	240	COMPUTER		
Amplitude BKS	Hoya Vision Care - Dallas	241	CLEARi OFFICE ROOM	CLEARi	104
Amplitude iQ	Hoya Vision Care - Dallas	242	Compact ULTRA HD	SOLA Optical	454
Amplitude Mini	Hoya Vision Care - Dallas	243	Computer	Shamir Insight Inc.	426
Amplitude Mini BKS	Hoya Vision Care - Dallas	244	Constitution HD	FEA Industries, Inc.	220
Amplitude Mini iQ	Hoya Vision Care - Dallas	245	Constitution Short HD	FEA Industries, Inc.	221
AO b'Active		1	Continual Focus Lens (CFL)	Vision Dynamics Lab/Optical	468
AO Easy		2		Dynamics	
AO Easy	American Optical	9	Continuum	Carl Zeiss Vision	16
AO Easy HD	American Optical	10	COZE ADL ADVANCED	Robertson Optical Laboratories, Inc.	368
AO Force 55		3			
AO Pro 15		4	COZE CASATA	Robertson Optical Laboratories, Inc.	369
AO Pro 16		5			
APEX	Classic Optical Laboratories, Inc.	97	COZE HDW	Robertson Optical Laboratories, Inc.	370
Array	Hoya Vision Care - Dallas	246	COZE SOFT	Robertson Optical Laboratories, Inc.	371
Array Wrap	Hoya Vision Care - Dallas	247			
Attitude III Fashion 15	Shamir Insight Inc.	416	COZE WRAP	Robertson Optical Laboratories, Inc.	372
Attitude III Fashion 18	Shamir Insight Inc.	417			
Attitude III Sport	Shamir Insight Inc.	418	Crossbows Custom U	Signet Armorlite, Inc.	436
Augen HD Trinity Spacia SF	Augen Optics	11	Crossbows Custom V	Signet Armorlite, Inc.	437
AugenSol HD Trinity 13/17	Augen Optics	12	Crossbows Junior	Signet Armorlite, Inc.	438
Autograph II -Fixed & Var. Des	Shamir Insight Inc.	419	DC mini	Seiko Optical Products	401
Autograph II Attitude	Shamir Insight Inc.	420	Definity	Essilor Canada	113
Autograph II Office	Shamir Insight Inc.	421	Definity	Essilor of America	130
Autograph III Fixed & VAR.	Shamir Insight Inc.	422	Definity Short	Essilor Canada	114
Autograph Intelligence Fixed & VARI.	Shamir Insight Inc.	423	Definity Short	Essilor of America	131
			Definity 3	Essilor of America	132
Autograph II+ Fixed & VARI.	Shamir Insight Inc.	424	Definity 3 Plus	Essilor of America	133
Autograph Relax	Shamir Insight Inc.	425	DirecTek	Signet Armorlite, Inc.	439
Balance Digital 12	Nikon Optical Canada	294	DirecTek Short	Signet Armorlite, Inc.	440
Balance Digital 14	Nikon Optical Canada	295	DLS Advanced PAL	Eye Kraft Optical, Inc.	215
Blue IQ Variable Progressive	Lenscrafters	285	DLS Anti-Fatigue	Eye Kraft Optical, Inc.	216
Ce-Tru Normal Corridor	L.B.I.	274	DLS Camber Near	Eye Kraft Optical, Inc.	217
Ce-Tru Polarized Progressive	L.B.I.	275	DLS PAL	Eye Kraft Optical, Inc.	218
Ce-Tru Short Corridor	L.B.I.	276	DLS Pro PAL	Eye Kraft Optical, Inc.	219
Ce-Tru Short Corridor Photochr	L.B.I.	277	Eagle	FEA Industries, Inc.	222
Ce-Tru SunSensor	L.B.I.	278	Eagle Apprentice	FEA Industries, Inc.	223
			Eagle Balance	FEA Industries, Inc.	224
			Eagle Driver	FEA Industries, Inc.	225
			Eagle PC	FEA Industries, Inc.	226
			Eagle Precision	FEA Industries, Inc.	227
			Eagle Sport	FEA Industries, Inc.	228
			Eagle Ultra Short	FEA Industries, Inc.	229
			Eagle Workplace	FEA Industries, Inc.	230

Lens Name	Company Name	Item #	Lens Name	Company Name	Item #
Eaglet	FEA Industries, Inc.	231	Geo-FORM 15	Rite-Style Optical Company, Inc.	333
Eaglet Balance	FEA Industries, Inc.	232	Geo-FORM AW 10	Rite-Style Optical Company, Inc.	334
EasyLux	Dehlawi Optical Industries	105	Geo-FORM AW 12	Rite-Style Optical Company, Inc.	335
Elan HD	Carl Zeiss Vision	17	Geo-FORM AW 14	Rite-Style Optical Company, Inc.	336
Elan HD Short	Carl Zeiss Vision	18	Geo-FORM EZ 13	Rite-Style Optical Company, Inc.	337
Elan HDV	Carl Zeiss Vision	19	Geo-FORM EZ 15	Rite-Style Optical Company, Inc.	338
Element	Shamir Insight Inc.	427	Geo-FORM EZ S 11	Rite-Style Optical Company, Inc.	339
Element Short	Shamir Insight Inc.	428	Geo-FORM EZ S 13	Rite-Style Optical Company, Inc.	340
EOS	KBco, Inc.	269	Golf	Shamir Insight Inc.	430
EOS Wrap	KBco, Inc.	270	Gradal 3	Carl Zeiss Vision	20
epik Day	GSRx	233	Gradal Brevis 1.5	Carl Zeiss Vision	21
epik DRV 50	GSRx	234	Gradal Brevis 1.6	Carl Zeiss Vision	22
epik iBifocal	GSRx	235	Gradal Brevity 1.59	Carl Zeiss Vision	23
epik Office HD	GSRx	236	Gradal Brevity 1.67	Carl Zeiss Vision	24
epik Standard	GSRx	237	Gradal HS (All)	Carl Zeiss Vision	25
epik Variable	GSRx	238	Gradal Individual	Carl Zeiss Vision	26
epik Variable Device	GSRx	239	Gradal RD	Carl Zeiss Vision	27
Essilor Anti-Fatigue	Essilor of America	134	Gradal Shorti	Carl Zeiss Vision	28
Essilor Computer	Essilor of America	135	Gradal Top	Carl Zeiss Vision	29
Essilor Ideal	Essilor of America	136	GT2	Carl Zeiss Vision	30
Essilor Ideal Advanced	Essilor of America	137	GT2 3D	Carl Zeiss Vision	31
Essilor Ideal Advanced Fit	Essilor of America	138	GT2 3D 3DV	Carl Zeiss Vision	32
Essilor Ideal Advanced Wrap	Essilor of America	139	GT2 3D Short	Carl Zeiss Vision	33
Essilor Ideal Computer	Essilor of America	140	GT2 Short	Carl Zeiss Vision	34
Essilor Ideal Driving	Essilor of America	141	HD Trinity 13/17 SF	Augen Optics	13
Essilor Ideal Fit	Essilor of America	142	HD Trinity 8/12 SF	Augen Optics	14
Essilor Ideal Short	Essilor of America	143	Hoyalux GP	Hoya Vision Care - Dallas	248
Essilor Ideal Sport	Essilor of America	144	Hoyalux GPW	Hoya Vision Care - Dallas	249
Essilor Natural	Essilor of America	145	Hoyalux Tact (EP40 and 60)	Hoya Vision Care - Dallas	250
Essilor Natural / Natural DS	Essilor Canada	115	HyperPro	PFO Optima LLC	327
Essilor Natural Digital	Essilor of America	146	ID InStyle	Hoya Vision Care - Dallas	251
Executive/ ShoreView ES	Shore Lens Company	433	ID Lifestyle	Hoya Vision Care - Dallas	252
Expert DNA™	Expert Optics	204	iD LifeStyle 2 Clarity	Hoya Vision Care - Dallas	253
Expert Elite C™	Expert Optics	205	ID LifeStyle 2 Clarity cd	Hoya Vision Care - Dallas	254
Expert Elite™	Expert Optics	206	ID LifeStyle 2 Harmony	Hoya Vision Care - Dallas	255
Expert Essential®	Expert Optics	207	ID LifeStyle 2 Harmony cd	Hoya Vision Care - Dallas	256
Expert Extreme C™	Expert Optics	208	ID Lifestyle cd	Hoya Vision Care - Dallas	257
Expert Extreme™	Expert Optics	209	ID Mystyle	Hoya Vision Care - Dallas	258
Expert Relief™ 50 and 75	Expert Optics	210	identity home 9/11/13	PFO Global	323
Expert RNA™	Expert Optics	211	identity play 9/11/13	PFO Global	324
Expert Sport Wrap XPV™	Expert Optics	212	identity pro	PFO Global	325
(progressive)			identity smart 9/11/13	PFO Global	326
Expert Sport Wrap XSV™ (single vision)	Expert Optics	213	Illumina	Vision-Ease Lens, Inc.	472
			Image	Younger Optics	491
Expert Tween™ 1.3; 2.0; 4.0	Expert Optics	214	Image	Younger Optics	492
EZ View Mini	SOMO Optical	462	Impression 2 XS, M, or L	Rodenstock Canada	378
EZ View Std	SOMO Optical	463	Impression Ergo 2 Book	Rodenstock Canada	379
First PAL	Shamir Insight Inc.	429	Impression Ergo 2 PC or Room	Rodenstock Canada	380
Freedom ID	X-CEL Optical	487-489	Impression Ergo FreeSign 2	Rodenstock Canada	381
Futureise	Polycore Optical	329	Impression FreeSign 2	Rodenstock Canada	382
Geo-FORM 11	Rite-Style Optical Company, Inc.	331			
Geo-FORM 13	Rite-Style Optical Company, Inc.	332			

Lens Name	Company Name	Item #	Lens Name	Company Name	Item #
Impression Hyperop	Rodenstock Canada	383	Nikon W	Nikon Optical USA	311
Impression Hyperop XS	Rodenstock Canada	384	Nikon W Short	Nikon Optical USA	312
Impression Myop	Rodenstock Canada	385	Novel	Vision-Ease Lens, Inc.	473
Impression Myop XS	Rodenstock Canada	386	Novella Short Corridor	Vision-Ease Lens, Inc.	474
Impression Sport 2	Rodenstock Canada	387	Outlook	Vision-Ease Lens, Inc.	475
inMotion	IOT	268	Ovation	Essilor Canada	116
Intellego 3	Interstate Optical	264	Ovation	Essilor of America	151
Intellego 3 14	Interstate Optical	265	Ovation ADS	Essilor Canada	117
Intellego 3 16	Interstate Optical	266	Ovation ADS (Canada)	Essilor of America	152
Intellego 3 18	Interstate Optical	267	Ovation ADS Short	Essilor of America	153
iRx Pro	KBco, Inc.	271	(Canada)		
iRx Short	KBco, Inc.	272	Ovation Digital	Essilor of America	154
Kodak Concise	Signet Armorlite, Inc.	441	Ovation DS (Canada)	Essilor of America	155
KODAK DSII	Signet Armorlite, Inc.	442	Overture	Sutherlin Optical, Inc.	464
Kodak Precise	Signet Armorlite, Inc.	443	Paradigm Progressive	Vision Dynamics Lab/Optical	469
KODAK Precise PB	Signet Armorlite, Inc.	444		Dynamics	
Kodak Precise Short	Signet Armorlite, Inc.	445	Paradigm Short Corridor	Vision Dynamics Lab/Optical	470
KODAK Precise Short PB	Signet Armorlite, Inc.	446		Dynamics	
KODAK SoftWear	Signet Armorlite, Inc.	447	PCWide	Seiko Optical Products	402
Kodak Unique	Signet Armorlite, Inc.	448	Percepta	SOLA Optical	455
KODAK Unique DRO	Signet Armorlite, Inc.	449	Perfas Premier II	Seiko Optical Products	403
KODAK Unique DRO HD	Signet Armorlite, Inc.	450	Perfas Prominent	Seiko Optical Products	404
Maui Passport Progressive	Maui Jim	288	Persona	Riverside Opticalab Ltd.	358
Melibrad Progressive	Melibrad	289	Persona Office	Riverside Opticalab Ltd.	359
Micro	Polycore Optical	330	Persona G4	Riverside Opticalab Ltd.	360
MOVE	Nikon Optical Canada	296	Persona Novice	Riverside Opticalab Ltd.	361
Move Digital 13	Nikon Optical USA	308	Persona Phd	Riverside Opticalab Ltd.	362
Move Digital 15	Nikon Optical USA	309	Persona Phd Distance	Riverside Opticalab Ltd.	363
MOVE Short	Nikon Optical Canada	297	Persona Phd Near	Riverside Opticalab Ltd.	364
Multigressiv MyView 2 XS, M, L	Rodenstock Canada	388	Persona Quattro	Riverside Opticalab Ltd.	365
Mvp Platinum Plastic II	Landon Lens	283	Persona Quattro Sport	Riverside Opticalab Ltd.	366
MyView Ergo 2 Book	Rodenstock Canada	389	Polar-Ray Progressive	Melibrad	290-291
MyView Ergo 2 PC or Room	Rodenstock Canada	390	PolySensus	MacLens Optical	286
Nalco Progressive	Nassau Lens Company	292	Presio Digital i Short	Nikon Optical USA	313
Natural DS (Canada)	Essilor of America	147	Presio i	Nikon Optical USA	314
Navigator FBS	Signet Armorlite, Inc.	451	Presio i Digital	Nikon Optical USA	315
Navigator Precision	Signet Armorlite, Inc.	452	Presio i Digital 12	Nikon Optical USA	316
Navigator Short	Signet Armorlite, Inc.	453	Presio i Digital 14	Nikon Optical USA	317
Nexyma Computer (80)	Rodenstock Canada	391	Presio i Digital Short	Nikon Optical USA	318
Nexyma Reading (40)	Rodenstock Canada	392	Presio i Short	Nikon Optical USA	319
Nikon Eyes Advanced	Essilor of America	148	Presio Power 12/13/14	Nikon Optical Canada	304
Nikon Eyes Advanced Short	Essilor of America	149	Presio Power 12/13/14	Nikon Optical USA	320
Nikon Eyes Customized	Essilor of America	150	Privo Advance	Dehlawi Optical Industries	106
Nikon Go	Nikon Optical Canada	298	Pro-Freedom	Excelite, Inc.	196
Nikon Go Digital	Nikon Optical Canada	299	Pro-Freedom Plus	Excelite, Inc.	197
Nikon i	Nikon Optical Canada	300	Pro-Freedom Solace	Excelite, Inc.	198
Nikon i Digital	Nikon Optical Canada	301	Pro-Freedom Sports	Excelite, Inc.	199
Nikon Online	Nikon Optical Canada	302	Proceed II Short	Seiko Optical Products	405
Nikon Online	Nikon Optical USA	310	Proceed III Super Short	Seiko Optical Products	406
Nikon W	Nikon Optical Canada	303	Progressiv ClassicLife	Rodenstock Canada	393
			Progressiv ClassicLife XS	Rodenstock Canada	394
			Progressiv Life 2	Rodenstock Canada	395
			Progressiv Life 2 XS	Rodenstock Canada	396

Lens Name	Company Name	Item #	Lens Name	Company Name	Item #
Progressiv Purelife	Rodenstock Canada	397	synchrony Easy S	Carl Zeiss Vision	45
Progressiv Purelife XS	Rodenstock Canada	398	synchrony Easy S HD	Carl Zeiss Vision	46
ProgressivPureLife Free XS,M,L	Rodenstock Canada	399	synchrony Easy View	Carl Zeiss Vision	47
Promult	MacLens Optical	287	synchrony Easy View HD	Carl Zeiss Vision	48
Proprius DS	Nexus Vision Group	293	synchrony Easy Wear	Carl Zeiss Vision	49
Punktal Gradal HS	Carl Zeiss Vision	35	synchrony Easy Wear HD	Carl Zeiss Vision	50
PureSite Easy Fit	Diversified Ophthalmics, Inc.	107	synchrony Performance HD (13, 15, 17, 19)	Carl Zeiss Vision	51
PureSite Enhanced	Diversified Ophthalmics, Inc.	108	synchrony Performance HDV	Carl Zeiss Vision	52
PureSite Premier	Diversified Ophthalmics, Inc.	109	synchrony Starter HD	Carl Zeiss Vision	53
PureSite Ultima	Diversified Ophthalmics, Inc.	110	synchrony Ultra HD (13, 15, 17, 19, 21)	Carl Zeiss Vision	54
Resolution Response	PFO Optima LLC	328	synchrony Ultra HDC	Carl Zeiss Vision	55
SeeMax Power 12/13/14	Nikon Optical Canada	305	synchrony Ultra HDV	Carl Zeiss Vision	56
SeeMax Power 12/13/14	Nikon Optical USA	321	synchrony Work & Go HD	Carl Zeiss Vision	57
SeeMaxProgressive 12	Nikon Optical Canada	306	synchrony Work & Office HD	Carl Zeiss Vision	58
SeeMaxProgressive 14	Nikon Optical Canada	307	synchrony Work & Read HD	Carl Zeiss Vision	59
ShoreView / Shoreview Mini	Shore Lens Company	434	TACT 40 BKS and TACT 60 BKS	Hoya Vision Care - Dallas	263
Shoreview Mini ES	Shore Lens Company	435	TR1-Far AW 14	Rite-Style Optical Company, Inc.	341
Signature HD	Riverside Opticalab Ltd.	367	TR1-Far AW 16	Rite-Style Optical Company, Inc.	342
Smallfit	Essilor of America	156	TR1-Far AW 18	Rite-Style Optical Company, Inc.	343
SmallFit / SmallFit DS	Essilor Canada	118	TR1-Far AW 20	Rite-Style Optical Company, Inc.	344
SmallFit Digital	Essilor of America	157	TR2-Basic 14	Rite-Style Optical Company, Inc.	345
Smallfit DS (Canada)	Essilor of America	158	TR2-Basic 16	Rite-Style Optical Company, Inc.	346
SOLA HDV	SOLA Optical	456	TR2-Basic 18	Rite-Style Optical Company, Inc.	347
SOLAMAX	Carl Zeiss Vision	36	TR2-Basic 20	Rite-Style Optical Company, Inc.	348
SOLAOne HD	SOLA Optical	457	TR3-Demand AW 14	Rite-Style Optical Company, Inc.	349
Spectrum	Shamir Insight Inc.	431	TR3-Demand AW 16	Rite-Style Optical Company, Inc.	350
Stealth 15	Vision Warehouse	471	TR3-Demand AW 18	Rite-Style Optical Company, Inc.	351
Succeed	Seiko Optical Products	407	TR3-Demand AW 20	Rite-Style Optical Company, Inc.	352
Succeed WS	Seiko Optical Products	408	TR4-Highway AW	Rite-Style Optical Company, Inc.	353
Summit cd	Hoya Vision Care - Dallas	259	TR5-Zone 13	Rite-Style Optical Company, Inc.	354
Summit cd iQ	Hoya Vision Care - Dallas	260	TR8-Active AW 15	Rite-Style Optical Company, Inc.	355
Summit ecp	Hoya Vision Care - Dallas	261	TR8-Active AW 17	Rite-Style Optical Company, Inc.	356
Summit ecp iQ	Hoya Vision Care - Dallas	262	TR8-Active AW 19	Rite-Style Optical Company, Inc.	357
Super No-Line	Essilor Canada	119	TruVision		6
Super No-Line	Essilor of America	159	TruVision Omni		7
Supercede	Seiko Optical Products	409	TruVision Technica		8
Supercede II	Seiko Optical Products	410	Ultimate HD	US Optical	466
Supercede WS	Seiko Optical Products	411	Umbramatic SR Gradal HS	Carl Zeiss Vision	60
Superior	Seiko Optical Products	412	UNITY PLx	VSP	476
Supernal	Seiko Optical Products	413	UNITY PLxtra	VSP	477
Surmount	Seiko Optical Products	414	Unity PLxtreme	VSP	478
Surmount Ws	Seiko Optical Products	415	Unity Via	VSP	479
Synchrony	Carl Zeiss Vision	37	Unity Via Elite	VSP	480
synchrony Access	Carl Zeiss Vision	38	Unity Via Mobile	VSP	481
synchrony Access HD (075, 125, 175)	Carl Zeiss Vision	39	Unity Via OfficePro 10ft	VSP	482
synchrony Curves HDC	Carl Zeiss Vision	40	Unity Via OfficePro 5ft	VSP	483
synchrony Easy Adapt	Carl Zeiss Vision	41	Unity Via Plus	VSP	484
synchrony Easy Adapt HD	Carl Zeiss Vision	42			
synchrony Easy M	Carl Zeiss Vision	43			
synchrony Easy M HD	Carl Zeiss Vision	44			

<u>Lens Name</u>	<u>Company Name</u>	<u>Item #</u>
Unity Via Wrap	VSP	485
Universa HD Office (075, 125, 175)	Carl Zeiss Vision	61
Universa2 HD	Carl Zeiss Vision	62

Lens Name	Company Name	Item #	Lens Name	Company Name	Item #
Variations	Landon Lens	284	VISTA WS	Robertson Optical	375
Varilux Ipseo IV	Essilor of America	160		Laboratories, Inc.	
Varilux Comfort	Essilor of America	161	Visuality	SOLA Optical	459
Varilux Comfort / Comfort DSx	Essilor Canada	120	WorkSpace	Shamir Insight Inc.	432
Varilux Comfort Drx	Essilor of America	162	Wrap HD	US Optical	467
Varilux Comfort Enhanced	Essilor of America	163	X-Pro Discovery	Excelite, Inc.	200
Varilux Comfort Short	Essilor of America	164	X-Pro Discovery Plus	Excelite, Inc.	201
Varilux Comfort Short DRx	Essilor of America	165	X-Pro Discovery Xtra	Excelite, Inc.	202
Varilux Comfort W3+	Essilor of America	166	X-Pro Maxima Plus	Excelite, Inc.	203
Varilux Comfort W3+ eyecode	Essilor of America	167	XL	SOLA Optical	460
Varilux Comfort W3+ Fit	Essilor of America	168	XL Gold	SOLA Optical	461
Varilux Ellipse	Essilor of America	169	XP HD	Robertson Optical	376
Varilux Ellipse / Ellipse 360	Essilor Canada	121		Laboratories, Inc.	
Varilux Ellipse 360°	Essilor of America	170	Xplore	KBco, Inc.	273
Varilux Ipseo IV	Essilor of America	171	Zeal Digital Progressive	Zeal Optics	493
Varilux Ipseo IV eyecode	Essilor of America	172	ZEISS Business	Carl Zeiss Vision	64
Varilux Physio	Essilor of America	173	ZEISS Choice (13, 15, 17, 19)	Carl Zeiss Vision	65
Varilux Physio 360°	Essilor of America	174	ZEISS Choice Plus (13, 15, 17, 19, 21)	Carl Zeiss Vision	66
Varilux Physio DRx	Essilor of America	175	ZEISS Choice Plus V	Carl Zeiss Vision	67
Varilux Physio Enhanced	Essilor of America	176	ZEISS Digital Lens	Carl Zeiss Vision	68
Varilux Physio Enhanced Azio	Essilor of America	177	ZEISS DriveSafe	Carl Zeiss Vision	69
Varilux Physio Enhanced eyecode	Essilor of America	178	ZEISS DriveSafe Individual	Carl Zeiss Vision	70
Varilux Physio Enhanced Fit	Essilor of America	179	ZEISS EnergizeMe (14, 16, 18)	Carl Zeiss Vision	71
Varilux Physio Enhanced India	Essilor of America	180	ZEISS EnergizeMe Digital	Carl Zeiss Vision	72
Varilux Physio Short	Essilor of America	181	ZEISS Experience	Carl Zeiss Vision	73-74
Varilux Physio Short 360°	Essilor of America	182	ZEISS Experience HD2	Carl Zeiss Vision	75
Varilux Physio Short DRx	Essilor of America	183	ZEISS Experience Short	Carl Zeiss Vision	76-77
Varilux Physio W3+	Essilor of America	184	ZEISS Individual	Carl Zeiss Vision	78
Varilux Physio W3+ eyecode	Essilor of America	185	ZEISS Individual 2	Carl Zeiss Vision	79
Varilux Physio W3+ Fit	Essilor of America	186	ZEISS Light D (14, 18)	Carl Zeiss Vision	80
Varilux S 4D	Essilor of America	187	ZEISS Light H (14, 16, 18)	Carl Zeiss Vision	81
Varilux S Design	Essilor of America	188	ZEISS Light V	Carl Zeiss Vision	82
Varilux S Design Short	Essilor of America	189	ZEISS MyoKids Pro	Carl Zeiss Vision	83
Varilux S Fit	Essilor of America	190	ZEISS MyoKids Pro Individual	Carl Zeiss Vision	84
Varilux Sport	Essilor Canada	122	ZEISS Officelens Book	Carl Zeiss Vision	85
Varilux Sport	Essilor of America	191	ZEISS Officelens Desk	Carl Zeiss Vision	86
Varilux Stylistic	Essilor of America	192	ZEISS Officelens Room	Carl Zeiss Vision	87
Varilux X	Essilor of America	193	ZEISS Plus 2 Short Sport	Carl Zeiss Vision	88
Varilux X 4D	Essilor of America	194	ZEISS Plus 2 Sport	Carl Zeiss Vision	89
Varilux X Fit	Essilor of America	195	ZEISS Precision Plus	Carl Zeiss Vision	90
VIP	Carl Zeiss Vision	63	ZEISS Precision Pure (13, 15, 17, 19, 21)	Carl Zeiss Vision	91
VIP Gold	SOLA Optical	458	ZEISS Precision Superb	Carl Zeiss Vision	92
Vista	Robertson Optical Laboratories, Inc.	373	ZEISS SmartLife Individual	Carl Zeiss Vision	93
VISTA CP	Robertson Optical Laboratories, Inc.	374	ZEISS SmartLife Plus	Carl Zeiss Vision	94
			ZEISS SmartLife Pure (S, M, L)	Carl Zeiss Vision	95
			ZEISS SmartLife Superb	Carl Zeiss Vision	96

Extended ADD Range Index

ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #
to 0.50	Acuity 9/11/13	322	to 0.50	Persona Quattro Sport	366	to 0.75	Amplitude BKS	241
	CLEARi Absolute-W 12	98		Privo Advance	106		Amplitude iQ	242
	CLEARi DED 14	99		Proceed II Short	405		Amplitude Mini	243
	CLEARi DED 16	100		Proceed III Super Short	406		Amplitude Mini BKS	244
	CLEARi DED 19	101		Signature HD	367		Amplitude Mini iQ	245
	CLEARi OFFICE BOOK	102		Smallfit DS (Canada)	158		AO Easy	2
	CLEARi OFFICE COMPUTER	103		Succeed	407		AO Easy HD	10
	CLEARi OFFICE ROOM	104		Succeed WS	408		AO Force 55	3
	COZE ADL ADVANCED	368		Summit ecp	261		Array	246
	COZE CASATA	369		Summit ecp iQ	262		Array Wrap	247
	COZE HDW	370		Supercede	409		Attitude III Fashion 15	416
	COZE SOFT	371		Supercede II	410		Attitude III Fashion 18	417
	COZE WRAP	372		Supercede WS	411		Attitude III Sport	418
	DLS Anti-Fatigue	216		Superior	412		Augen HD Trinity Spacia SF	11
	DLS Camber Near	217		Supernal	413		Autograph II -Fixed & Var. Des	419
	Eagle	222		Surmount	414		Autograph II Attitude	420
	Eagle Apprentice	223		Surmount Ws	415		Autograph II Office	421
	Eagle Balance	224		synchrony Starter HD	53		Autograph III Fixed & VAR.	422
	Eagle Driver	225		TR1-Far AW 14	341		Autograph Intelligence Fixed & VARI.	423
	Eagle Precision	227		TR1-Far AW 16	342		Autograph II+ Fixed & VARI.	424
	Eagle Sport	228		TR1-Far AW 18	343		Balance Digital 12	294
	Eagle Ultra Short	229		TR1-Far AW 20	344		Balance Digital 14	295
	Eagle Workplace	230		TR2-Basic 14	345		Ce-Tru Normal Corridor	274
	Eaglet	231		TR2-Basic 16	346		ClassicLife Free XS. M or L	377
	Eaglet Balance	232		TR2-Basic 18	347		Compact ULTRA HD	454
	EasyLux	105		TR2-Basic 20	348		Computer	426
	epik DRV 50	234		TR3-Demand AW 14	349		Constitution HD	220
	Expert Elite™	206		TR3-Demand AW 16	350		Constitution Short HD	221
	Expert Essential®	207		TR3-Demand AW 18	351		Crossbows Custom U	436
	Expert Extreme™	209		TR3-Demand AW 20	352		Crossbows Custom V	437
	Expert Relief™ 50 and 75	210		TR4-Highway AW	353		Crossbows Junior	438
	Expert Tween™ 1.3; 2.0; 4.0	214		TR5-Zone 13	354		DirecTek	439
	Geo-FORM 11	331		TR8-Active AW 15	355		DirecTek Short	440
	Geo-FORM 13	332		TR8-Active AW 17	356		DLS Advanced PAL	215
	Geo-FORM 15	333		TR8-Active AW 19	357		DLS PAL	218
	Geo-FORM AW 10	334		Unity Via	479		DLS Pro PAL	219
	Geo-FORM AW 12	335		Unity Via Elite	480		Elan HD	17
	Geo-FORM AW 14	336		Unity Via Mobile	481		Elan HD Short	18
	Geo-FORM AW 14	336		Unity Via Plus	484		Elan HDV	19
	Geo-FORM EZ 13	337		Unity Via Wrap	485		Element	427
	Geo-FORM EZ 15	338		Vista	373		Element Short	428
	Geo-FORM EZ S 11	339		VISTA CP	374		epik Day	233
	Geo-FORM EZ S 13	340		VISTA WS	375		epik iBifocal	235
	identity home 9/11/13	323		XP HD	376		epik Standard	237
	identity play 9/11/13	324	to 0.75	Accolade	123		epik Variable	238
	identity smart 9/11/13	326		Accolade Freedom	124		epik Variable Device	239
	Perfas Premier II	403		Adage	490		Essilor Ideal	136
	Perfas Prominent	404		Adaptar	125		Essilor Ideal Advanced	137
	Persona	358		Adaptar Digital	126		Essilor Ideal Advanced Fit	138
	Persona Office	359		Adaptar DS	112		Essilor Ideal Advanced Wrap	139
	Persona G4	360		Adaptar DS (Canada)	127		Essilor Ideal Computer	140
	Persona Novice	361		Adaptar DS Short (Canada)	128		Essilor Ideal Driving	141
	Persona Phd	362		Adaptar Short Digital	129		Essilor Ideal Fit	142
	Persona Phd Distance	363		Advanced HD	465		Essilor Ideal Short	143
	Persona Phd Near	364		Amplitude	240		Essilor Ideal Sport	144
	Persona Quattro	365						

ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #
to 0.75	Essilor Natural	145	to 0.75	MyView Ergo 2 Book	389	to 0.75	SOLASOne HD	457
	Essilor Natural / Natural DS	115		MyView Ergo 2 PC or Room	390		Spectrum	431
	Essilor Natural Digital	146		Natural DS (Canada)	147		Summit cd	259
	Executive/ ShoreView ES	433		Navigator FBS	451		Summit cd iQ	260
	Expert DNA™	204		Nikon Eyes Advanced	148		Super No-Line	159
	Expert Elite C™	205		Nikon Eyes Advanced Short	149		synchrony Access HD (075, 125, 175)	39
	Expert Extreme C™	208		Nikon Eyes Customized	150		synchrony Curves HDC	40
	Expert RNA™	211		Nikon Go	298		synchrony Easy Adapt HD	42
	First PAL	429		Nikon Go Digital	299		synchrony Easy M	43
	Golf	430		Nikon i	300		synchrony Easy M HD	44
	Gradal Individual	26		Nikon i Digital	301		synchrony Easy S	45
	Gradal Shorti	28		Nikon W	303		synchrony Easy S HD	46
	Gradal Top	29		Nikon W	311		synchrony Easy View	47
	GT2	30		Nikon W Short	312		synchrony Easy View HD	48
	GT2 3D	31		Novel	473		synchrony Easy Wear HD	50
	GT2 3D Short	33		Outlook	475		synchrony Performance HD (13, 15, 17, 19)	51
	GT2 Short	34		Ovation	116		synchrony Performance HDV	52
	Hoyalux GP	248		Ovation	151		synchrony Ultra HD (13, 15, 17, 19, 21)	54
	Hoyalux GPW	249		Ovation ADS	117		synchrony Ultra HDC	55
	HyperPro	327		Ovation ADS (Canada)	152		synchrony Ultra HDV	56
	ID InStyle	251		Ovation ADS Short (Canada)	153		synchrony Work & Go HD	57
	ID Lifestyle	252		Ovation Digital	154		synchrony Work & Office HD	58
	iD LifeStyle 2 Clarity	253		Ovation DS (Canada)	155		synchrony Work & Read HD	59
	ID LifeStyle 2 Clarity cd	254		Percepta	455		TACT 40 BKS and TACT 60 BKS	263
	ID LifeStyle 2 Harmony	255		Presio Digital i Short	313		Ultimate HD	466
	ID LifeStyle 2 Harmony cd	256		Presio i	314		Umbramatic SR Gradal HS	60
	ID Lifestyle cd	257		Presio i Digital	315		UNITY PLx	476
	ID Mystyle	258		Presio i Digital 12	316		UNITY PLxtra	477
	identity pro	325		Presio i Digital 14	317		Unity PLxtreme	478
	Illumina	472		Presio i Digital Short	318		Universa HD Office (075, 125, 175)	61
	Impression 2 XS, M, or L	378		Presio i Short	319		Universa2 HD	62
	Impression Ergo 2 Book	379		Presio Power 12/13/14	304		Varilux Ipseo IV	160
	Impression Ergo 2 PC or Room	380		Presio Power 12/13/14	320		Varilux Comfort	161
	Impression Ergo FreeSign 2	381		Pro-Freedom Solace	198		Varilux Comfort Drx	162
	Impression FreeSign 2	382		Progressiv ClassicLife	393		Varilux Comfort Enhanced	163
	Impression Hyperop	383		Progressiv ClassicLife XS	394		Varilux Comfort Short	164
	Impression Hyperop XS	384		Progressiv Life 2	395		Varilux Comfort Short DRx	165
	Impression Myop	385		Progressiv Life 2 XS	396		Varilux Comfort W3+	166
	Impression Myop XS	386		Progressiv Purelife	397		Varilux Comfort W3+ eyecode	167
	Intellego 3	264		Progressiv Purelife XS	398		Varilux Comfort W3+ Fit	168
	Intellego 3 14	265		ProgressivPureLife Free	399		Varilux Ellipse	169
	Intellego 3 16	266		XS,M,L			Varilux Ellipse / Ellipse 360	121
	Intellego 3 18	267		Proprius DS	293		Varilux Ellipse 360°	170
	KODAK DSII	442		PureSite Easy Fit	107		Varilux Ipseo IV	171
	KODAK Precise PB	444		PureSite Enhanced	108		Varilux Ipseo IV eyecode	172
	KODAK Precise Short PB	446		PureSite Premier	109		Varilux Physio	173
	KODAK SoftWear	447		PureSite Ultima	110		Varilux Physio 360°	174
	Kodak Unique	448		SeeMax Power 12/13/14	305		Varilux Physio DRx	175
	KODAK Unique DRO	449		SeeMax Power 12/13/14	321		Varilux Physio Enhanced	176
	KODAK Unique DRO HD	450		SeeMaxProgressive 12	306		Varilux Physio Enhanced Azio	177
	MOVE	296		SeeMaxProgressive 14	307		Varilux Physio Enhanced eyecode	178
	Move Digital 13	308		Shoreview Mini ES	435		Varilux Physio Enhanced Fit	179
	Move Digital 15	309		Smallfit	156		Varilux Physio Enhanced India	180
	MOVE Short	297		SmallFit / SmallFit DS	118		Varilux Physio Short	181
	Multigressiv MyView 2 XS, M, L	388		SmallFit Digital	157		Varilux Physio Short 360°	182
				SOLA HDV	456			

ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #
to 0.75	Varilux Physio Short DRx	183	to 3.50	Advanced HD	465	to 3.50	Essilor Ideal	136
	Varilux Physio W3+	184		Amplitude	240		Essilor Ideal Advanced	137
	Varilux Physio W3+ eyecode	185		Amplitude BKS	241		Essilor Ideal Advanced Fit	138
	Varilux Physio W3+ Fit	186		Amplitude iQ	242		Essilor Ideal Computer	140
	Varilux S 4D	187		Amplitude Mini	243		Essilor Ideal Driving	141
	Varilux S Design	188		Amplitude Mini BKS	244		Essilor Ideal Fit	142
	Varilux S Design Short	189		Amplitude Mini iQ	245		Essilor Ideal Sport	144
	Varilux S Fit	190		AO b'Active	1		Essilor Natural	145
	Varilux Sport	122		AO Easy	2		Essilor Natural / Natural DS	115
	Varilux Sport	191		AO Easy HD	10		Essilor Natural Digital	146
	Varilux Stylistic	192		AO Force 55	3		Executive/ ShoreView ES	433
	Varilux X	193		APEX	97		Expert DNA™	204
	Varilux X 4D	194		Array	246		Expert Elite C™	205
	Varilux X Fit	195		Array Wrap	247		Expert Extreme C™	208
	VIP	63		Attitude III Fashion 15	416		Expert RNA™	211
	VIP Gold	458		Attitude III Fashion 18	417		Expert Tween™ 1.3; 2.0; 4.0	214
	Visuality	459		Attitude III Sport	418		First PAL	429
	WorkSpace	432		Augen HD Trinity Spacia SF	11		Futureise	329
	Wrap HD	467		Autograph II Attitude	420		Geo-FORM 11	331
	XL	460		Autograph II Office	421		Geo-FORM 13	332
	XL Gold	461		Autograph Intelligence Fixed & VARI.	423		Geo-FORM 15	333
	ZEISS Choice (13, 15, 17, 19)	65		Autograph II+ Fixed & VARI.	424		Geo-FORM AW 10	334
	ZEISS Choice Plus (13, 15, 17, 19, 21)	66		Ce-Tru Normal Corridor	274		Geo-FORM AW 12	335
	ZEISS Choice Plus V	67		CLEARi Absolute-W 12	98		Geo-FORM AW 14	336
	ZEISS Digital Lens	68		CLEARi DED 14	99		Geo-FORM EZ 13	337
	ZEISS DriveSafe	69		CLEARi DED 16	100		Geo-FORM EZ 15	338
	ZEISS DriveSafe Individual	70		CLEARi DED 19	101		Geo-FORM EZ S 11	339
	ZEISS EnergizeMe (14, 16, 18)	71		CLEARi OFFICE BOOK	102		Geo-FORM EZ S 13	340
	ZEISS EnergizeMe Digital	72		CLEARi OFFICE COMPUTER	103		Golf	430
	ZEISS Experience HD2	75		CLEARi OFFICE ROOM	104		Gradal Individual	26
	ZEISS Individual	78		Compact ULTRA HD	454		Gradal Shorti	28
	ZEISS Individual 2	79		Constitution HD	220		Gradal Top	29
	ZEISS Light D (14, 18)	80		Constitution Short HD	221		GT2	30
	ZEISS Light H (14, 16, 18)	81		COZE HDW	370		GT2 3D	31
	ZEISS Light V	82		COZE SOFT	371		GT2 3D 3DV	32
	ZEISS Officelens Book	85		COZE WRAP	372		GT2 3D Short	33
	ZEISS Officelens Desk	86		Crossbows Custom U	436		GT2 Short	34
	ZEISS Officelens Room	87		Crossbows Custom V	437		Hoyalux GP	248
	ZEISS Plus 2 Short Sport	88		Crossbows Junior	438		Hoyalux GPW	249
	ZEISS Plus 2 Sport	89		DirecTek	439		HyperPro	327
	ZEISS Precision Plus	90		DirecTek Short	440		ID InStyle	251
	ZEISS Precision Pure (13, 15, 17, 19, 21)	91		DLS Advanced PAL	215		ID Lifestyle	252
	ZEISS Precision Superb	92		DLS PAL	218		iD LifeStyle 2 Clarity	253
	ZEISS SmartLife Individual	93		DLS Pro PAL	219		ID LifeStyle 2 Clarity cd	254
	ZEISS SmartLife Plus	94		Elan HD	17		ID LifeStyle 2 Harmony	255
	ZEISS SmartLife Pure (S, M, L)	95		Elan HD Short	18		ID LifeStyle 2 Harmony cd	256
	ZEISS SmartLife Superb	96		Elan HDV	19		ID Lifestyle cd	257
to 3.50	Accolade	123		Element	427		ID Mystyle	258
	Accolade Freedom	124		Element Short	428		Illumina	472
	Acuity 9/11/13	322		epik Day	233		Impression 2 XS, M, or L	378
	Adage	490		epik iBifocal	235		Impression Ergo 2 Book	379
	Adaptar	125		epik Office HD	236		Impression Ergo 2 PC or Room	380
	Adaptar DS	112		epik Standard	237		Impression Ergo FreeSign 2	381
	Adaptar DS (Canada)	127		epik Variable	238		Impression FreeSign 2	382
	Adaptar DS Short (Canada)	128		epik Variable Device	239		Intellego 3	264

ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #
to 3.50	Intellego 3 14	265	to 3.50	Proceed II Short	405	to 3.50	UNITY PLxtra	477
	Intellego 3 16	266		Proceed III Super Short	406		Unity PLxtreme	478
	Intellego 3 18	267		Progressiv Life 2	395		Unity Via OfficePro 10ft	482
	iRx Pro	271		Progressiv Purelife	397		Unity Via OfficePro 5ft	483
	KODAK DSII	442		ProgressivPureLife Free XS,M,L	399		Universa HD Office (075, 125, 175)	61
	KODAK Precise PB	444		Promult	287		Universa2 HD	62
	KODAK Precise Short PB	446		Proprius DS	293		Varilux Ipseo IV	160
	KODAK SoftWear	447		PureSite Easy Fit	107		Varilux Comfort	161
	Kodak Unique	448		PureSite Enhanced	108		Varilux Comfort Drx	162
	KODAK Unique DRO	449		PureSite Premier	109		Varilux Comfort Enhanced	163
	KODAK Unique DRO HD	450		PureSite Ultima	110		Varilux Ipseo IV	171
	Micro	330		SeeMax Power 12/13/14	305		Varilux Ipseo IV eyecode	172
	MOVE	296		SeeMax Power 12/13/14	321		Varilux Physio Enhanced	176
	Move Digital 13	308		SeeMaxProgressive 12	306		Varilux Sport	122
	Move Digital 15	309		SeeMaxProgressive 14	307		Varilux Sport	191
	MOVE Short	297		Smallfit DS (Canada)	158		Varilux Stylistic	192
	Multigressiv MyView 2 XS, M, L	388		SOLA HDV	456		VIP	63
	Mvp Platinum Plastic II	283		SOLAOne HD	457		VIP Gold	458
	MyView Ergo 2 Book	389		Spectrum	431		Vista	373
	MyView Ergo 2 PC or Room	390		Succeed	407		VISTA CP	374
	Natural DS (Canada)	147		Succeed WS	408		VISTA WS	375
	Navigator FBS	451		Summit cd	259		Visuality	459
	Nikon Eyes Advanced	148		Summit cd iQ	260		WorkSpace	432
	Nikon Eyes Customized	150		Super No-Line	159		Wrap HD	467
	Nikon Go	298		Supercede	409		X-Pro Discovery	200
	Nikon Go Digital	299		Supercede II	410		X-Pro Discovery Plus	201
	Nikon i	300		Supercede WS	411		X-Pro Discovery Xtra	202
	Nikon i Digital	301		Supernal	413		XL	460
	Nikon W	303		Surmount	414		XL Gold	461
	Nikon W	311		Surmount Ws	415		ZEISS Choice (13, 15, 17, 19)	65
	Nikon W Short	312		Synchrony	37		ZEISS Choice Plus (13, 15, 17, 19, 21)	66
	Novel	473		synchrony Access HD (075, 125, 175)	39		ZEISS Choice Plus V	67
	Novella Short Corridor	474		synchrony Curves HDC	40		ZEISS Digital Lens	68
	Outlook	475		synchrony Easy Adapt	41		ZEISS DriveSafe	69
	Ovation	116		synchrony Easy Adapt HD	42		ZEISS DriveSafe Individual	70
	Ovation	151		synchrony Easy M	43		ZEISS EnergizeMe (14, 16, 18)	71
	Ovation ADS	117		synchrony Easy M HD	44		ZEISS EnergizeMe Digital	72
	Ovation ADS (Canada)	152		synchrony Easy S	45		ZEISS Experience HD2	75
	Ovation Digital	154		synchrony Easy S HD	46		ZEISS Individual 2	79
	Ovation DS (Canada)	155		synchrony Easy View	47		ZEISS Light D (14, 18)	80
	Percepta	455		synchrony Easy View HD	48		ZEISS Light H (14, 16, 18)	81
	Perfas Premier II	403		synchrony Easy Wear	49		ZEISS Light V	82
	Perfas Prominent	404		synchrony Easy Wear HD	50		ZEISS Officelens Book	85
	PolySensus	286		synchrony Performance HD (13, 15, 17, 19)	51		ZEISS Officelens Desk	86
	Presio Digital i Short	313		synchrony Performance HDV	52		ZEISS Officelens Room	87
	Presio i	314		synchrony Ultra HD (13, 15, 17, 19, 21)	54		ZEISS Plus 2 Short Sport	88
	Presio i Digital	315		synchrony Ultra HDC	55		ZEISS Plus 2 Sport	89
	Presio i Digital Short	318		synchrony Ultra HDV	56		ZEISS Precision Plus	90
	Presio i Short	319		synchrony Work & Go HD	57		ZEISS Precision Pure (13, 15, 17, 19, 21)	91
	Presio Power 12/13/14	304		synchrony Work & Office HD	58		ZEISS Precision Superb	92
	Presio Power 12/13/14	320		synchrony Work & Read HD	59		ZEISS SmartLife Individual	93
	Pro-Freedom	196		TACT 40 BKS and TACT 60 BKS	263		ZEISS SmartLife Plus	94
	Pro-Freedom Plus	197		Ultimate HD	466		ZEISS SmartLife Pure (S, M, L)	95
	Pro-Freedom Solace	198		Umbromatic SR Gradal HS	60		ZEISS SmartLife Superb	96
	Pro-Freedom Sports	199		UNITY PLx	476	to 4.00	Autograph II -Fixed & Var. Des	419

ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #	ADD Range	Lens Name	Item #
to 4.00	Autograph III Fixed & VAR.	422	to	Varilux Physio Enhanced	178	to 4.50	Expert Elite™	206
	Balance Digital 12	294	4.00	eyecode			Expert Essential®	207
	Balance Digital 14	295		Varilux Physio Enhanced Fit	179		Expert Extreme™	209
	Definity 3	132		Varilux Physio Enhanced India	180		Persona Quattro	365
	Definity 3 Plus	133		Varilux Physio W3+	184		Persona Quattro Sport	366
	EasyLux	105		Varilux Physio W3+ eyecode	185		TR1-Far AW 14	341
	identity home 9/11/13	323		Varilux Physio W3+ Fit	186		TR1-Far AW 16	342
	identity play 9/11/13	324		Varilux S 4D	187		TR1-Far AW 18	343
	identity smart 9/11/13	326		Varilux S Design	188		TR1-Far AW 20	344
	Persona	358		Varilux S Design Short	189		TR2-Basic 14	345
	Persona G4	360		Varilux S Fit	190		TR2-Basic 16	346
	Persona Phd	362		Varilux X	193		TR2-Basic 18	347
	Persona Phd Distance	363		Varilux X 4D	194		TR2-Basic 20	348
	Persona Phd Near	364		Varilux X Fit	195		TR3-Demand AW 14	349
	Presio i Digital 12	316		XP HD	376		TR3-Demand AW 16	350
	Presio i Digital 14	317	to	COZE ADL ADVANCED	368		TR3-Demand AW 18	351
	Privo Advance	106	4.50				TR3-Demand AW 20	352
	Signature HD	367		COZE CASATA	369		TR4-Highway AW	353
	Summit ecp	261		DLS Camber Near	217		TR5-Zone 13	354
	Summit ecp iQ	262		Eagle	222		TR8-Active AW 15	355
	Superior	412		Eagle Apprentice	223		TR8-Active AW 17	356
	Varilux Comfort W3+	166		Eagle Balance	224		TR8-Active AW 19	357
	Varilux Comfort W3+ eyecode	167		Eagle Driver	225		Unity Via	479
	Varilux Comfort W3+ Fit	168		Eagle PC	226		Unity Via Elite	480
	Varilux Physio	173		Eagle Precision	227		Unity Via Mobile	481
	Varilux Physio 360°	174		Eagle Sport	228		Unity Via Plus	484
	Varilux Physio DRx	175		Eagle Ultra Short	229		Unity Via Wrap	485
	Varilux Physio Enhanced Azio	177		Eagle Workplace	230			
				Eaglet	231			
				Eaglet Balance	232			

Index of Companies

Company Name	Item #	Company Name	Item #
American Optical	9	Nikon Optical Canada	294
Augen Optics	11	Nikon Optical USA	308
Carl Zeiss Vision	15	PFO Global	322
Classic Optical Laboratories, Inc.	97	PFO Optima LLC	327
CLEARi	98	Polycore Optical	329
Dehlawi Optical Industries	105	Rite-Style Optical Company, Inc.	331
Diversified Ophthalmics, Inc.	107	Riverside Opticalab Ltd.	358
Essilor Canada	111	Robertson Optical Laboratories, Inc.	368
Essilor of America	123	Rodenstock Canada	377
Excelite, Inc.	196	Seiko Optical Products	400
Expert Optics	204	Shamir Insight Inc.	416
Eye Kraft Optical, Inc.	215	Shore Lens Company	433
FEA Industries, Inc.	220	Signet Armorlite, Inc.	436
GSRx	233	SOLA Optical	454
Hoya Vision Care - Dallas	240	SOMO Optical	462
Interstate Optical	264	Sutherlin Optical, Inc.	464
IOT	268	US Optical	465
KBco, Inc.	269	Vision Dynamics Lab/Optical Dynamics	468
L.B.I.	274	Vision Warehouse	471
Landon Lens	279	Vision-Ease Lens, Inc.	472
Lenscrafters	285	VSP	476
MacLens Optical	286	Walman - Soderberg	486
Maui Jim	288	X-CEL Optical	487
Melibrad	289	Younger Optics	490
Nassau Lens Company	292	Zeal Optics	493
Nexus Vision Group	293		

Index by Recommended Minimum Height

Min Height	Lens Name	Item #	Min Height	Lens Name	Item #	Min Height	Lens Name	Item #
10mm	Eagle Ultra Short	229	13mm	synchrony Performance HD (13, 15, 17, 19)	51	14mm	Expert Essential®	207
11mm	AO Easy	9		synchrony Performance HDV	52		Expert Extreme C™	208
	Autograph II -Fixed & Var. Des	419		synchrony Performance	54		Expert Extreme™	209
	Autograph III Fixed & VAR.	422		synchrony Ultra HD (13, 15, 17, 19, 21)	55		Expert RNA™	211
	Autograph Intelligence Fixed & VARI.	423		synchrony Ultra HDC	56		Expert Tween™ 1.3; 2.0; 4.0	214
	Autograph II+ Fixed & VARI.	424		Universa2 HD	62		Geo-FORM AW 10	334
	Blue IQ Variable Progressive	285		Varilux Ellipse 360°	170		ID InStyle	251
12mm	Unity Via	479		ZEISS Choice (13, 15, 17, 19)	65		ID LifeStyle 2 Clarity cd	254
	Unity Via Elite	480		ZEISS Choice Plus (13, 15, 17, 19, 21)	66		ID LifeStyle 2 Harmony cd	256
	Unity Via Mobile	481		ZEISS Choice Plus V	67		ID Lifestyle cd	257
	Unity Via Plus	484		ZEISS Digital Lens	68		Impression 2 XS, M, or L	378
	Unity Via Wrap	485		ZEISS DriveSafe	69		Impression Ergo 2 Book	379
13mm	Adage	490		ZEISS DriveSafe Individual	70		Impression Ergo FreeSign 2	381
	Autograph III Fixed & VAR.	422		ZEISS Experience HD2	75		Impression Hyperop XS	384
	Autograph Intelligence Fixed & VARI.	423		ZEISS Experience Short	76,77		Impression Myop XS	386
	Autograph II+ Fixed & VARI.	424		ZEISS Individual	78		iRx Short	272
	Compact ULTRA HD	454		ZEISS Individual 2	79		Kodak Concise	441
	Crossbows Custom U	436		ZEISS Light V	82		Multigressiv MyView 2 XS, M, L	388
	Crossbows Custom V	437		ZEISS Officelens Book	85		MyView Ergo 2 Book	389
	Crossbows Junior	438		ZEISS Officelens Desk	86		Navigator Short	453
	DLS Advanced PAL	215		ZEISS Officelens Room	87		Nikon Eyes Advanced Short	149
	DLS PAL	218		ZEISS Precision Plus	90		Novella Short Corridor	474
	DLS Pro PAL	219		ZEISS Precision Pure (13, 15, 17, 19, 21)	91		Perfas Premier II	403
	Elan HD Short	18		ZEISS Precision Superb	92		Perfas Prominent	404
	Elan HDV	19		ZEISS SmartLife Individual	93		Progressiv ClassicLife XS	394
	epik Day	233		ZEISS SmartLife Plus	94		Progressiv Life 2 XS	396
	Essilor Anti-Fatigue	134		ZEISS SmartLife Superb	96		Progressiv Purelife XS	398
	Essilor Ideal Computer	140	13.5mm	DirecTek Short	440		ProgressivPureLife Free XS,M,L	399
	Expert DNA™	204	14mm	Adaptar Short Digital	129		PureSite Enhanced	108
	Expert RNA™	211		Advanced HD	465		PureSite Premier	109
	GT2 3D 3DV	32		Amplitude Mini	243		PureSite Ultima	110
	GT2 3D Short	33		Amplitude Mini BKS	244		Signature HD	367
	GT2 Short	34		Amplitude Mini iQ	245		Smallfit	156
	Impression FreeSign 2	382		APEX	97		SmallFit Digital	157
	Intellego 3	264		Channel 14 Plastic II	279		Spectrum	431
	Intellego 3 14	265		Channel 14 Poly II	280		Summit cd	259
	Intellego 3 16	266		Channel 14 Variations Poly	282		Summit cd iQ	260
	Intellego 3 18	267		Pho			Supercede II	410
	KODAK DSII	442		ClassicLife Free XS. M or L	377		Supercede WS	411
	Kodak Precise Short	445		CLEARi Absolute-W 12	98		Supernal	413
	KODAK Precise Short PB	446		Computer	426		Surmount	414
	Kodak Unique	448		DLS Anti-Fatigue	216		Surmount Ws	415
	KODAK Unique DRO	449		DLS Camber Near	217		synchrony Starter HD	53
	KODAK Unique DRO HD	450		Eagle	222		TR1-Far AW 14	341
	Nexyma Reading (40)	392		Eagle Balance	224		TR2-Basic 14	345
	Persona G4	360		Eagle PC	226		TR3-Demand AW 14	349
	Persona Novice	361		Eagle Precision	227		Ultimate HD	466
	Persona Phd	362		Eaglet	231		Unity PLxtreme	478
	Persona Phd Distance	363		Eaglet Balance	232		Unity Via OfficePro 10ft	482
	Persona Phd Near	364		epik Office HD	236		Unity Via OfficePro 5ft	483
	Persona Quattro	365		epik Standard	237		Varilux Ipseo IV	160
	SOLA HDV	456		Essilor Ideal Fit	142		Varilux Comfort Enhanced	163
	synchrony Curves HDC	40		Essilor Ideal Short	143		Varilux Comfort Short	164
	synchrony Easy M	43		Essilor Natural Digital	146		Varilux Comfort Short DRx	165
	synchrony Easy S HD	46		Expert Elite C™	205		Varilux Comfort W3+	166
				Expert Elite™	206			


Min Height	Lens Name	Item #	Min Height	Lens Name	Item #	Min Height	Lens Name	Item #
14mm	Varilux Comfort W3+	167	15mm	Gradal Brevis 1.5	21	16mm	identity play 9/11/13	324
	eyecode			Gradal Brevis 1.6	22		identity smart 9/11/13	326
	Varilux Comfort W3+ Fit	168		Gradal Shorti	28		inMotion	268
	Varilux Ellipse	169		ID Mystyle	258		Micro	330
	Varilux Ellipse / Ellipse 360	121		identity pro	325		Move Digital 13	308
	Varilux Ipseo IV	171		Nikon Eyes Customized	150		MOVE Short	297
	Varilux Ipseo IV eyecode	172		Nikon Online	302,310		Nikon i	300
	Varilux Physio Enhanced	176		Nikon W Short	312		Nikon i Digital	301
	Varilux Physio Enhanced	177		Ovation ADS Short (Canada)	153		Novel	473
	Azio			PCWide	402		Paradigm Short	470
	Varilux Physio Enhanced	178		Persona Office	359		Corridor	
	eyecode			Persona Quattro Sport	366		Polar-Ray Progressive	291
	Varilux Physio Enhanced Fit	179		Presio i Digital 12	316		PolySensus	286
	Varilux Physio Enhanced	180		Pro-Freedom	196		Presio Digital i Short	313
	India			Pro-Freedom Plus	197		Presio i Digital Short	318
	Varilux Physio Short	181		Pro-Freedom Sports	199		Presio i Short	319
	Varilux Physio Short 360°	182		PureSite Easy Fit	107		Proceed III Super Short	406
	Varilux Physio Short DRx	183		SeeMaxProgressive 12	306		Proprius DS	293
	Varilux Physio W3+	184		ShoreView / Shoreview Mini	434		Spectrum	431
	Varilux Physio W3+ eyecode	185		Shoreview Mini ES	435		synchrony Work &	58
	Varilux Physio W3+ Fit	186		SmallFit / SmallFit DS	118		Office HD	
	Varilux S 4D	187		Smallfit DS (Canada)	158		synchrony Work &	59
	Varilux S Design Short	189		Succeed WS	408		Read HD	
	Varilux S Fit	190		synchrony Access	38		TR1-Far AW 16	342
	WorkSpace	432		synchrony Access HD (075, 125, 175)	39		TR2-Basic 16	346
	Wrap HD	467		synchrony Easy M HD	44		TR3-Demand AW 16	350
	ZEISS EnergizeMe (14, 16, 18)	71		synchrony Easy S	45		ZEISS Business	64
	ZEISS EnergizeMe Digital	72		TR8-Active AW 15	355		ZEISS EnergizeMe (14, 16, 18)	71
	ZEISS Light D (14, 18)	80		Universa HD Office (075, 125, 175)	61		ZEISS Experience	73,74
	ZEISS Light H (14, 16, 18)	81		ZEISS Choice (13, 15, 17, 19)	65		ZEISS Light H (14, 16, 18)	81
	ZEISS MyoKids Pro	83		ZEISS Precision Pure (13, 15, 17, 19, 21)	91		ZEISS SmartLife Pure (S, M, L)	95
	ZEISS MyoKids Pro Individual	84	16mm	Acuity 9/11/13	322	16.5mm	DirectTek	439
	ZEISS Plus 2 Short Sport	88		Autograph II Office	421		Accolade	123
	ZEISS SmartLife Pure (S, M, L)	95		Autograph Relax	425		AF2	400
15mm	Accolade Freedom	124		Ce-Tru Short Corridor	276		AugenSol HD Trinity 13/17	12
	Adaptar DS Short (Canada)	128		Ce-Tru Short Corridor	277		Balance Digital 14	295
	Attitude III Fashion 15	416		Photochr			CLEARi OFFICE BOOK	102
	Autograph II Attitude	420		Channel 14 Variations Gray & B	281		Constitution HD	220
	Autograph III Fixed & VAR.	422		Clarlet Business	15		Elan HD	17
	Autograph Intelligence Fixed & VARI.	423		CLEARi DED 14	99		EOS	269
	Autograph II+ Fixed & VARI.	424		Eagle Apprentice	223		EOS Wrap	270
	Balance Digital 12	294		Eagle Sport	228		Essilor Ideal	136
	Constitution Short HD	221		Element Short	428		Essilor Ideal Advanced	137
	Continuum	16		Expert Elite C™	205		Essilor Ideal Advanced	138
	Definity Short	114,131		Expert Elite™	206		Fit	
	Definity 3	132		Expert Essential®	207		Essilor Ideal Advanced	139
	Definity 3 Plus	133		Expert Extreme C™	208		Wrap	
	epik Variable	238		Expert Extreme™	209		Expert DNA™	204
	epik Variable Device	239		Expert RNA™	211		Expert Elite C™	205
	Essilor Computer	135		Expert Sport Wrap XPV™ (progressive)	212		Expert Elite™	206
	Expert DNA™	204		EZ View Mini	462		Expert Essential®	207
	Expert Elite C™	205		Geo-FORM 11	331		Expert Extreme C™	208
	Expert Elite™	206		Geo-FORM AW 12	335		Expert Extreme™	209
	Expert Essential®	207		Gradal Brevity 1.59	23		Expert RNA™	211
	Expert Extreme C™	208		Gradal Brevity 1.67	24		Freedom ID	487,488,489
	Expert Extreme™	209		identity home 9/11/13	323		Geo-FORM EZ 13	337
	Expert RNA™	211					Geo-FORM EZ S 13	340
	First PAL	429					GT2	30
	Geo-FORM EZ S 11	339					GT2 3D	31
							HD Trinity 13/17 SF	13

Min Height	Lens Name	Item #	Min Height	Lens Name	Item #	Min Height	Lens Name	Item #		
17mm	HyperPro	327	18mm	Eagle Driver	225	18mm	SOLAOne HD	457		
	Illumina	472		Eagle Workplace	230		Stealth 15	471		
	Kodak Precise	443		Essilor Ideal Driving	141		Summit ecp	261		
	KODAK Precise PB	444		Essilor Ideal Sport	144		Summit ecp iQ	262		
	KODAK SoftWear	447		Essilor Natural	145		synchrony Easy Adapt	41		
	Nikon Eyes Advanced	148		Essilor Natural / Natural DS	115		synchrony Easy Adapt HD	42		
	Nikon W	303,311		Executive/ ShoreView ES	433		synchrony Easy Wear	49		
	Ovation	116,151		Expert Elite C™	205		synchrony Easy Wear HD	50		
	Ovation ADS	117		Expert Elite™	206		TACT 40 BKS and TACT 60	263		
	Ovation ADS (Canada)	152		Expert Essential®	207		BKS			
	Ovation Digital	154		Expert Extreme C™	208		TR1-Far AW 18	343		
	Overture	464		Expert Extreme™	209		TR2-Basic 18	347		
	Presio i Digital 14	317		Expert Relief™ 50 and 75	210		TR3-Demand AW 18	351		
	Pro-Freedom Solace	198		Expert RNA™	211		TR4-Highway AW	353		
	SeeMaxProgressive 14	307		Expert Sport Wrap XPV™	212		TR5-Zone 13	354		
	SOLAMAX	36		(progressive)			TruVision Technica	8		
	Succeed	407		Expert Tween™ 1.3; 2.0; 4.0	214		Variations	284		
	Synchrony	37		Futureise	329		Varilux Comfort / Comfort	120		
	synchrony Easy View	47		Geo-FORM 13	332		DSx			
	synchrony Easy View HD	48		Geo-FORM AW 14	336		Varilux Sport	122,191		
	synchrony Work & Go HD	57		Golf	430		Varilux Stylistic	192		
	TR8-Active AW 17	356		Gradal Individual	26		VIP	63		
	Varilux Comfort	161		Gradal Top	29		VIP Gold	458		
	Varilux Comfort Drx	162		Hoyalux GP	248		Visuality	459		
	Varilux Physio	173		Hoyalux GPW	249		X-Pro Discovery Xtra	202		
	Varilux Physio 360°	174		Hoyalux Tact (EP40 and 60)	250		X-Pro Maxima Plus	203		
	Varilux Physio DRx	175		ID Lifestyle	252		XL	460		
	Varilux S Design	188		iD LifeStyle 2 Clarity	253		XL Gold	461		
	Varilux X	193		ID LifeStyle 2 Harmony	255		Zeal Digital Progressive	493		
	Varilux X 4D	194		Image	491		ZEISS EnergizeMe (14, 16, 18)	71		
	Varilux X Fit	195		Image	492		ZEISS Light D (14, 18)	80		
	Xplore	273		Impression Ergo 2 PC or Room	380		ZEISS Light H (14, 16, 18)	81		
	ZEISS Choice (13, 15, 17, 19)	65		Impression Hyperop	383		ZEISS Plus 2 Sport	89		
	ZEISS Precision Pure (13, 15, 17, 19, 21)	91		Impression Myop	385		ZEISS SmartLife Pure (S, M, L)	95		
	18mm	Adaptar		125	Impression Sport 2		387	19mm	Amplitude	240
		Adaptar Digital		126	iRx Pro		271		Amplitude BKS	241
		Adaptar DS		112	Melibrad Progressive		289		Amplitude iQ	242
		Adaptar DS (Canada)		127	MOVE		296		DC mini	401
		AO b'Active		1	Move Digital 15		309		Element	427
		AO Easy		2	Mvp Platinum Plastic II		283		Expert DNA™	204
AO Easy HD		10	MyView Ergo 2 PC or Room	390	Expert Elite C™	205				
AO Force 55		3	Natural DS (Canada)	147	Expert Elite™	206				
AO Pro 15		4	Navigator FBS	451	Expert Essential®	207				
AO Pro 16		5	Nikon Go	298	Expert Extreme C™	208				
Attitude III Fashion 18		417	Nikon Go Digital	299	Expert Extreme™	209				
Attitude III Sport		418	Outlook	475	Geo-FORM EZ 15	338				
Autograph II Attitude		420	Ovation DS (Canada)	155	Nalco Progressive	292				
Autograph III Fixed & VAR.		422	Percepta	455	Spectrum	431				
Autograph Intelligence Fixed & VARI.		423	Presio i	314	TR8-Active AW 19	357				
Autograph II+ Fixed & VARI.		424	Presio i Digital	315	ZEISS Choice (13, 15, 17, 19)	65				
CLEARi DED 16		100	Proceed II Short	405	ZEISS Precision Pure (13, 15, 17, 19, 21)	91				
CLEARi OFFICE COMPUTER		103	Progressiv ClassicLife	393	20mm	Expert Elite C™	205			
CLEARi OFFICE ROOM		104	Progressiv Life 2	395		Expert Elite™	206			
Continual Focus Lens (CFL)		468	Progressiv Purelife	397		Expert Essential®	207			
Definity		113,130	Promult	287		Expert Extreme C™	208			
			Resolution Response	328		Expert Extreme™	209			
						EZ View Std	463			

Min Height	Lens Name	Item #	Min Height	Lens Name	Item #	Min Height	Lens Name	Item #
20mm	Geo-FORM 15	333	22mm	Ce-Tru SunSensor	278	Several	Expert Essential®	207
	Gradal 3	20		Gradal HS (All)	25		Expert Extreme™	209
	Navigator Precision	452		Polar-Ray Progressive	290		Expert RNA™	211
	Nexyma Computer (80)	391		Punktal Gradal HS	35		Maui Passport Progressive	288
	Paradigm Progressive	469		TruVision	6		Persona	358
	TR1-Far AW 20	344		TruVision Omni	7		Presio Power 12/13/14	304,320
	TR2-Basic 20	348		Umbramatic SR Gradal HS	60		Privo Advance	106
	TR3-Demand AW 20	352	Several	Acclaim OEM PAL #1 and #2**	486		SeeMax Power 12/13/14	305,321
	X-Pro Discovery	200		Array	246		Supercede	409
	X-Pro Discovery Plus	201		Array Wrap	247		UNITY PLx	476
21mm	CLEARi DED 19	101		COZE ADL ADVANCED	368		UNITY PLxtra	477
	Gradal RD	27		COZE CASATA	369		Vista	373
	Super No-Line	159		COZE HDW	370		VISTA CP	374
	ZEISS Precision Pure (13, 15, 17, 19, 21)	91		COZE SOFT	371		VISTA WS	375
				COZE WRAP	372		XP HD	376
22mm	Augen HD Trinity Spacia SF	11		EasyLux	105			
	Ce-Tru Normal Corridor	274		Expert Elite™	206			
	Ce-Tru Polarized Progressive	275						

AO b'Active

1




Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada

AO Easy

2



no symbol

P

Materials:
Plastic (CR39 and 1.50)
Polycarbonate



bA

Filters:
Photochromic
Polarized

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

AO Force 55

3



Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada

AO Pro 15

4



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada

AO Pro 16

5




H

Materials:
Plastic High Index (1.60)

PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada

TruVision
6



M

Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 22mm
Available as: Clear
Available in: US and Canada

TruVision Omni
7



M

Materials:
Plastic (CR39 and 1.50)
Glass (1.523)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 22mm
Available as: Clear; Photochromic;
Available in: US and Canada

TruVision Technica
8



M

Materials:
Plastic (CR39 and 1.50)

PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear
Available in: US and Canada

American Optical | AO Easy
9



M

7

Materials:
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height: 11mm
Available as: Clear; Photochromic;
Available in: US and Canada

American Optical | AO Easy HD
10



M

H

T

P

F

M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Augen Optics | Augen HD Trinity Spacia SF

11



Materials:
Trivex

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
22mm
Available as: Clear
Available in: US and Canada

Augen Optics | AugenSol HD Trinity 13/17

12



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic;
Available in: US and Canada

Augen Optics | HD Trinity 13/17 SF

13



no symbol

no symbol

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic;
Available in: US and Canada

Augen Optics | HD Trinity 8/12 SF

14



Materials:
Plastic (CR39 and 1.50)
Trivex

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
12
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | Clarlet Business

15



Materials:
Plastic (CR39 and 1.50)

Near Variable Focus
PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height:
16mm
Available as: Clear
Available in: Canada Only

Carl Zeiss Vision | Continuum

16



M

CP

Materials:
Polycarbonate

Near Variable Focus
PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height:
15mm
Available as: Clear
Available in: US and Canada

Carl Zeiss Vision | Elan HD

17



M

H

T

P

F

M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | Elan HD Short

18



M

H

T

P

F

M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | Elan HDV

19



M

H

T

P

F

M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | Gradal 3

20




M

Materials:
Plastic (CR39 and 1.50)

PAL Design on Front side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height:
20mm
Available as: Clear
Available in: Canada Only

Carl Zeiss Vision | Gradal Brevis 1.5

21



M


50

Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear
Available in: Canada Only

Carl Zeiss Vision | Gradal Brevis 1.6

22



M

60

Materials:
Plastic High Index (1.60)

Standard Progressive
PAL Design on Front side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear
Available in: Canada Only

Carl Zeiss Vision | Gradal Brevity 1.59

23



M

Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 3mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | Gradal Brevity 1.67

24




M

Materials:
Plastic High Index (1.66 - 1.67)

PAL Design on Front side
Fitting cross location: 3mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | Gradal HS (All)

25




M

Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height:
22mm
Available as: Clear
Available in: US and Canada

Carl Zeiss Vision | Gradal Individual

26



M

60

66

Materials:

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height: 18mm

Available as: Clear

Available in: US and Canada

Carl Zeiss Vision | Gradal RD

27



M

Materials:

Plastic (CR39 and 1.50)

Standard Progressive

PAL Design on Front side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height: 21mm

Available as: Clear

Available in: US and Canada

Carl Zeiss Vision | Gradal Shorti

28



M

60

Materials:

Plastic High Index (1.60)

Standard Progressive

PAL Design on Back side

Fitting cross location: 6mm above 180 line


Rec. minimum fitting height: 15mm

Available as: Clear

Available in: US and Canada

Carl Zeiss Vision | Gradal Top

29



M

50

58

60

66

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Front side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height: 18mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | GT2

30



M

50

59

60

67

60

80

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic Ultra High Index (> 1.67)

Glass High Index (> 1.60)

Glass Ultra High Index (> 1.70)

Standard Progressive

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 17mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | GT2 3D

31



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | GT2 3D 3DV

32



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | GT2 3D Short

33



M	Materials:
50	Plastic (CR39 and 1.50)
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | GT2 Short

34




M	Materials:
50	Plastic (CR39 and 1.50)
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | Punktal Gradal HS

35



M

Materials:
Glass High Index (> 1.60)

PAL Design on Front side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height:
22mm
Available as: Clear
Available in: US and Canada

Carl Zeiss Vision | SOLAMAX

36



M

H
P

Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | Synchrony

37



M

H
P

Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | synchrony Access

38



M

AC
ACP

Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Near Variable Focus
PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height:
15mm
Available as: Clear
Available in: US and Canada

Carl Zeiss Vision | synchrony Access HD (075, 125, 175)

39



M


H
T
P
F
M
J

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Curves HDC

40



M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy Adapt

41




M	Materials:
H	Plastic (CR39 and 1.50)
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy Adapt HD

42



M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy M

43




M	Materials:
C	Plastic (CR39 and 1.50)
P	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy M HD

44




M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy S

45




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy S HD

46



M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy View

47




M	Materials:
C	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
F	Filters:
	Photochromic

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy View HD

48




M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy Wear

49




M	Materials:
H	Plastic (CR39 and 1.50)
P	Polycarbonate
M	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Easy Wear HD

50




M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Performance HD (13, 15, 17, 19)

51




M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Performance HDV

52



M

H
T
P
F
M
J

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Starter HD

53



M


H
T
P
F
M
J

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Ultra HD (13, 15, 17, 19, 21)

54



M

H
T
P
F
M
M
J

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Ultra HDC

55



M


H
T
P
F
M
J

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Ultra HDV

56



M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Work & Go HD

57



M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Work & Office HD

58




M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | synchrony Work & Read HD

59




M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic High Index (1.66 - 1.67)
J	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | Umbramatic SR Gradal HS

60



M

Materials:
Glass (1.523)

PAL Design on Front side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height:
22mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | Universa HD Office (075, 125, 175)

61

□

M


H
T
P
F
M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
15mm
Available as: Clear
Available in: Canada Only

Carl Zeiss Vision | Universa2 HD

62



M

H
T
P
F
M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Carl Zeiss Vision | VIP

63



M

S

Materials:
Glass (1.523)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Carl Zeiss Vision | ZEISS Business

64



M

Materials:
Plastic (CR39 and 1.50)

Near Variable Focus
PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height:
16mm
Available as: Clear
Available in: US and Canada

Carl Zeiss Vision | ZEISS Choice (13, 15, 17, 19)

65



M	Materials:	Standard Progressive
50	Plastic (CR39 and 1.50)	PAL Design on Back side
53	Trivex	Fitting cross location: 6mm above 180 line
59	Polycarbonate	Rec. minimum fitting height:
60	Plastic High Index (1.60)	13mm
67	Plastic High Index (1.66 - 1.67)	15mm
		17mm
		19mm
		Available as: Clear; Photochromic; Polarized
		Available in: US and Canada

Carl Zeiss Vision | ZEISS Choice Plus (13, 15, 17, 19, 21)

66



M	Materials:	Standard Progressive
50	Plastic (CR39 and 1.50)	PAL Design on Back side
53	Trivex	Fitting cross location: 6mm above 180 line
59	Polycarbonate	Rec. minimum fitting height:
60	Plastic High Index (1.60)	13mm
67	Plastic High Index (1.66 - 1.67)	Available as: Clear; Photochromic; Polarized
74	Plastic Ultra High Index (> 1.67)	Available in: US and Canada

Carl Zeiss Vision | ZEISS Choice Plus V

67



M	Materials:	Standard Progressive
50	Plastic (CR39 and 1.50)	PAL Design on Back side
53	Trivex	Fitting cross location: 6mm above 180 line
59	Polycarbonate	Rec. minimum fitting height:
60	Plastic High Index (1.60)	13mm
67	Plastic High Index (1.66 - 1.67)	Available as: Clear; Photochromic; Polarized
74	Plastic Ultra High Index (> 1.67)	Available in: US and Canada

Carl Zeiss Vision | ZEISS Digital Lens

68



M	Materials:	Standard Progressive
50	Plastic (CR39 and 1.50)	PAL Design on Back side
53	Trivex	Fitting cross location: 6mm above 180 line
59	Polycarbonate	Rec. minimum fitting height:
60	Plastic High Index (1.60)	13mm
67	Plastic High Index (1.66 - 1.67)	Available as: Clear; Photochromic; Polarized
74	Plastic Ultra High Index (> 1.67)	Available in: US and Canada

Carl Zeiss Vision | ZEISS DriveSafe

69



M

50

53

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height:


13mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | ZEISS DriveSafe Individual

70



M

50

53

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height:

13mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | ZEISS EnergizeMe (14, 16, 18)

71



M

50

53

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height:

14mm

16mm

18mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | ZEISS EnergizeMe Digital

72



M

50

53

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height:


14mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | ZEISS Experience

73




Materials:
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 3mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US Only

Carl Zeiss Vision | ZEISS Experience

74




Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 3mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US Only

Carl Zeiss Vision | ZEISS Experience HD2

75



59
67

Materials:
Polycarbonate
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | ZEISS Experience Short

76




Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic;
Available in: US Only

Carl Zeiss Vision | ZEISS Experience Short

77




Materials:
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic;
Available in: US Only

Carl Zeiss Vision | ZEISS Individual

78



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | ZEISS Individual 2

79



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | ZEISS Light D (14, 18)

80



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 14mm
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | ZEISS Light H (14, 16, 18)

81




M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 14mm
16mm
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | ZEISS Light V

82



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | ZEISS MyoKids Pro

83



M	Materials:
50	Plastic (CR39 and 1.50)
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: Canada Only

Carl Zeiss Vision | ZEISS MyoKids Pro Individual

84



M	Materials:
50	Plastic (CR39 and 1.50)
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: Canada Only

Carl Zeiss Vision | ZEISS Officelens Book

85



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear
Available in: US and Canada

Carl Zeiss Vision | ZEISS Officelens Desk

86

	M	Materials:	Standard Progressive
50		Plastic (CR39 and 1.50)	PAL Design on Back side
53		Trivex	Fitting cross location: 6mm above 180 line
59		Polycarbonate	Rec. minimum fitting height:
60		Plastic High Index (1.60)	13mm
67		Plastic High Index (1.66 - 1.67)	Available as: Clear
74		Plastic Ultra High Index (> 1.67)	Available in: US and Canada

Carl Zeiss Vision | ZEISS Officelens Room

87

	M	Materials:	Standard Progressive
50		Plastic (CR39 and 1.50)	PAL Design on Back side
53		Trivex	Fitting cross location: 6mm above 180 line
59		Polycarbonate	Rec. minimum fitting height:
60		Plastic High Index (1.60)	13mm
67		Plastic High Index (1.66 - 1.67)	Available as: Clear
74		Plastic Ultra High Index (> 1.67)	Available in: US and Canada

Carl Zeiss Vision | ZEISS Plus 2 Short Sport

88

ADD p2 M	M	Materials:	Standard Progressive
50		Plastic (CR39 and 1.50)	PAL Design on Back side
53		Trivex	Fitting cross location: 6mm above 180 line
59		Polycarbonate	Rec. minimum fitting height:
60		Plastic High Index (1.60)	14mm
67		Plastic High Index (1.66 - 1.67)	Available as: Clear; Photochromic; Polarized
			Available in: US and Canada

Carl Zeiss Vision | ZEISS Plus 2 Sport

89

ADD P2 M	M	Materials:	Standard Progressive
50		Plastic (CR39 and 1.50)	PAL Design on Back side
53		Trivex	Fitting cross location: 6mm above 180 line
59		Polycarbonate	Rec. minimum fitting height:
60		Plastic High Index (1.60)	18mm
67		Plastic High Index (1.66 - 1.67)	Available as: Clear; Photochromic; Polarized
			Available in: US and Canada

Carl Zeiss Vision | ZEISS Precision Plus

90



M

50

53

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height:

13mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | ZEISS Precision Pure (13, 15, 17, 19, 21)

91



M

50

53

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height:

13mm

15mm

17mm

19mm

21mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | ZEISS Precision Superb

92



M

50

53

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height:

13mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | ZEISS SmartLife Individual

93



M

50

53

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 6mm above 180 line

Rec. minimum fitting height:


13mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Carl Zeiss Vision | ZEISS SmartLife Plus

94



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | ZEISS SmartLife Pure (S, M, L)

95



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 14mm
16mm
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Carl Zeiss Vision | ZEISS SmartLife Superb

96



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 6mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Classic Optical Laboratories, Inc. | APEX

97



M	Materials:
P	Plastic (CR39 and 1.50)
T	Trivex
C	Polycarbonate

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

CLEARi | CLEARi Absolute-W 12

98



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

CLEARi | CLEARi DED 14

99



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

CLEARi | CLEARi DED 16

100



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

CLEARi | CLEARi DED 19

101



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 21mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

CLEARi | CLEARi OFFICE BOOK

102



M

50

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: On 180 line

Rec. minimum fitting height:

17mm

Available as: Clear

Available in: US and Canada

CLEARi | CLEARi OFFICE COMPUTER

103



M

50

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 2mm above 180 line

Rec. minimum fitting height:

18mm

Available as: Clear

Available in: US and Canada

CLEARi | CLEARi OFFICE ROOM

104



M

50

59

60

67

74

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

18mm

Available as: Clear

Available in: US and Canada

Dehlawi Optical Industries | EasyLux

105



M

50

53

56

58

60

67

74

Materials:

Plastic (CR39 and 1.50)

Plastic Mid-Index (1.52 - 1.58)

Plastic Mid-Index (1.52 - 1.58)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

Several

Available as: Clear; Photochromic; Polarized

Available in: US and Canada



M	Materials:
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
no symbol	Plastic Mid-Index (1.52 - 1.58)
no symbol	Trivex
no symbol	Polycarbonate
no symbol	Plastic High Index (1.60)
no symbol	Plastic High Index (1.66 - 1.67)
no symbol	Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
no symbol	Plastic Mid-Index (1.52 - 1.58)
no symbol	Trivex
no symbol	Polycarbonate
no symbol	Plastic High Index (1.60)
no symbol	Plastic High Index (1.66 - 1.67)
no symbol	Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Diversified Ophthalmics, Inc. | PureSite Premier109



- M

no symbol

Materials:

Plastic (CR39 and 1.50)
- no symbol

Plastic Mid-Index (1.52 - 1.58)
- no symbol

Trivex
- no symbol

Polycarbonate
- no symbol

Plastic High Index (1.60)
- no symbol

Plastic High Index (1.66 - 1.67)
- no symbol

Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Diversified Ophthalmics, Inc. | PureSite Ultima110



- M

no symbol

Materials:

Plastic (CR39 and 1.50)
- no symbol

Plastic Mid-Index (1.52 - 1.58)
- no symbol

Trivex
- no symbol

Polycarbonate
- no symbol

Plastic High Index (1.60)
- no symbol

Plastic High Index (1.66 - 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor Canada | Adaptar111




- M

Materials:

Plastic High Index (1.60)

Fitting cross location:
Rec. minimum fitting height:
Available as: Clear
Available in:

Essilor Canada | Adaptar DS112



- M

III

Materials:

Plastic (CR39 and 1.50)
- IIIP

Polycarbonate
- III6

Plastic High Index (1.60)
- III67

Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: Canada Only

Essilor Canada | Definity

113



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor Canada | Definity Short

114



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor Canada | Essilor Natural / Natural DS

115



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Essilor Canada | Ovation

116



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Essilor Canada

Ovation ADS

117



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Essilor Canada

SmallFit / SmallFit DS

118



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor Canada

Super No-Line

119



M	Materials:
	Plastic (CR39 and 1.50)

Fitting cross location:
Rec. minimum fitting height:
Available as: Clear
Available in:

Essilor Canada

Varilux Comfort / Comfort DSx

120



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
no symbol	Glass High Index (> 1.60)
8	Glass Ultra High Index (> 1.70)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Essilor Canada

Varilux Ellipse / Ellipse 360

121



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Essilor Canada | Varilux Sport

122




Materials:
Polycarbonate

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear Polarized
Available in: US and Canada

Essilor of America | Accolade

123




*
P
6
67
74

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Accolade Freedom

124



no symbol
*
P
6
74

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Adaptar

125




no symbol
P

Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Essilor of America | Adaptar Digital

126



no symbol
*
P
67

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Adaptar DS (Canada)

127



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Essilor of America

Adaptar DS Short (Canada)

128




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.66 - 1.67)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Essilor of America

Adaptar Short Digital

129



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Definity

130




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Definity Short

131




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Definity 3

132




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Definity 3 Plus

133




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Essilor Anti-Fatigue

134



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height: 13mm
Available as: Clear
Available in: US and Canada

Essilor of America

Essilor Computer

135



Materials:
Polycarbonate

Near Variable Focus
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear
Available in: US Only

Essilor of America

Essilor Ideal

136



no symbol



P

6

67


Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Essilor Ideal Advanced

137



no symbol



P

6

67


Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Essilor Ideal Advanced Fit

138



no symbol



P

60

67

74


Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Essilor Ideal Advanced Wrap

139



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Essilor Ideal Computer

140



M	Materials:
*	Plastic (CR39 and 1.50)
P	Trivex
6	Polycarbonate
67	Plastic High Index (1.60)
74	Plastic High Index (1.66 - 1.67)
	Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Essilor Ideal Driving

141



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Essilor Ideal Fit

142




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Essilor Ideal Short

143



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Essilor Ideal Sport

144




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Essilor Natural

145




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Essilor of America | Essilor Natural Digital

146




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Natural DS (Canada)

147



M

no symbol

P

6

67

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 18mm

Available as: Clear; Photochromic; Polarized

Available in: Canada Only

Essilor of America

Nikon Eyes Advanced

148



M

P

Materials:

Polycarbonate

Standard Progressive

PAL Design on Front and Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 17mm

Available as: Clear; Photochromic;

Available in: US Only

Essilor of America

Nikon Eyes Advanced Short

149



M

P

Materials:

Polycarbonate

Standard Progressive

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 14mm

Available as: Clear; Photochromic;

Available in: US Only

Essilor of America

Nikon Eyes Customized

150



M

67

Materials:

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Front and Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 15mm


Available as: Clear; Photochromic;

Available in: US Only

Essilor of America

Ovation

151



M

no symbol

P

67

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 17mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Essilor of America

Ovation ADS (Canada)

152



M

no symbol

P

67

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Front and Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 17mm

Available as: Clear; Photochromic; Polarized

Available in: Canada Only

Essilor of America

Ovation ADS Short (Canada)

153



M

no symbol

67

Materials:

Plastic (CR39 and 1.50)

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Front and Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 15mm

Available as: Clear; Photochromic;

Available in: Canada Only

Essilor of America

Ovation Digital

154



M

no symbol

*

P

6

67

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 17mm

Available as: Clear; Photochromic; Polarized

Available in: US Only

Essilor of America

Ovation DS (Canada)

155



M

no symbol

P

6

67

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 18mm

Available as: Clear; Photochromic; Polarized

Available in: Canada Only

Essilor of America

Smallfit

156



M

no symbol

P

6

67

Materials:

Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 14mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Essilor of America | SmallFit Digital

157



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Smallfit DS (Canada)

158



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.66 - 1.67)
67	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Essilor of America | Super No-Line

159



M	Materials:
no symbol	Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
21mm
Available as: Clear
Available in: US Only

Essilor of America | Varilux Ipseo IV

160




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
Π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic;
Available in: US Only

Essilor of America | Varilux Comfort

161




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
Π	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America | Varilux Comfort Drx

162



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America | Varilux Comfort Enhanced

163



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Varilux Comfort Short

164



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Varilux Comfort Short DRx

165




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux Comfort W3+

166



M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

*

P

6


π

74

Essilor of America

Varilux Comfort W3+ eyecode

167



M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

*

P

6

π

74

Essilor of America

Varilux Comfort W3+ Fit

168



M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

*

P

6


π

74

Essilor of America

Varilux Ellipse

169



M

no symbol
Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada


P

6

π

Essilor of America | Varilux Ellipse 360°

170




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America | Varilux Ipseo IV

171



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Essilor of America | Varilux Ipseo IV eyecode

172




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Varilux Physio

173




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America | Varilux Physio 360°

174




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America | Varilux Physio DRx

175




M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Varilux Physio Enhanced

176




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America | Varilux Physio Enhanced Azio

177




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Varilux Physio Enhanced eyecode

178



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
Π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Varilux Physio Enhanced Fit

179




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
Π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Varilux Physio Enhanced India

180




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
Π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Varilux Physio Short

181




M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
Π	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US Only

Essilor of America

Varilux Physio Short 360°

182



M

P

π

Materials:
Polycarbonate


Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Essilor of America

Varilux Physio Short DRx

183



M

no symbol

P

6

π

74


Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Essilor of America

Varilux Physio W3+

184



M

*

P

6

π

74


Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux Physio W3+ eyecode

185



M

*

P

6

π

74


Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux Physio W3+ Fit

186



M	Materials:
	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux S 4D

187



M	Materials:
	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux S Design

188



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux S Design Short

189



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux S Fit

190



M	Materials:
no symbol	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux Sport

191



M	Materials:
P	Polycarbonate

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Polarized
Available in: US and Canada

Essilor of America

Varilux Stylistic

192



M	Materials:
	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux X

193



M	Materials:
	Plastic (CR39 and 1.50)
*	Trivex
P	Polycarbonate
π	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: ; Photochromic; Polarized
Available in: US and Canada

Essilor of America

Varilux X 4D

194



M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: ; Photochromic; Polarized
Available in: US and Canada

★

P

π

74

Essilor of America

Varilux X Fit

195



M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: ; Photochromic; Polarized
Available in: US and Canada

★

P

π

74

Excelite, Inc.

Pro-Freedom

196



M

Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 15mm
Available as: Clear
Available in: US and Canada

Excelite, Inc.

Pro-Freedom Plus

197



M


Materials:
Plastic High Index (1.60)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear
Available in: US and Canada

Excelite, Inc.

Pro-Freedom Solace

198



M

1

Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Excelite, Inc. | Pro-Freedom Sports

199



M

Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear
Available in: US Only

Excelite, Inc. | X-Pro Discovery

200



M

G

V

M

Materials:
Plastic (CR39 and 1.50)

Trivex

Plastic High Index (1.60)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
20mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Excelite, Inc. | X-Pro Discovery Plus

201



M

C

H

Materials:
Plastic (CR39 and 1.50)

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
20mm
Available as: Clear; Photochromic;
Available in: US and Canada

Excelite, Inc. | X-Pro Discovery Xtra

202



M

V

H

Materials:
Plastic (CR39 and 1.50)
Trivex

Plastic High Index (1.60)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Excelite, Inc. | X-Pro Maxima Plus

203



M

Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)

Near Variable Focus
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada



M	
B50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
PRT	Trivex
TPT	Trivex
B53	Trivex
H53	Trivex
T53	Trivex
PTX	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
BLY	Polycarbonate
ZPY	Polycarbonate
B60	Plastic High Index (1.60)
HTR	Plastic High Index (1.60)
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
B67	Plastic High Index (1.66 - 1.67)


Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
15mm
17mm
19mm
P Digitally balanced design for an economically conscious patient with 4 fitting heights. Also available with Blue Blocking filter.
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
P	
TRN	
PRY	
H53	
T53	
PLY	
TPY	
PLP	
TPP	
H60	
T60	
P60	
HTR	
H67	
T67	
P67	
H74	
T74	

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
15mm
16mm
17mm
18mm
19mm
20mm
P Advanced two-sided digitally well-balanced design that is customized and optimized using Position of Wear values (actual or standard)
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
B50	Plastic (CR39 and 1.50)
H53	Trivex
B53	Trivex
T53	Trivex
PTX	Trivex
PRT	Trivex
TPT	Trivex
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
ZPY	Polycarbonate
BLY	Polycarbonate
PLY	Polycarbonate
HTR	Plastic High Index (1.60)
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
B60	Plastic High Index (1.60)
P67	Plastic High Index (1.66 - 1.67)
H67	Plastic High Index (1.66 - 1.67)
B67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H70	Plastic Ultra High Index (> 1.67)
T70	Plastic Ultra High Index (> 1.67)
H74	Plastic Ultra High Index (> 1.67)
T74	Plastic Ultra High Index (> 1.67)
B74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
15mm
16mm
17mm
18mm
19mm
20mm
Several
P Digitally Personalized and Optimized with Position of Wear (actual or average) for Everyday Use. Also available with Blue Blocking filter.
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

**M****P****TRN****PRY****TPR****B50****H53****B53****T53****PTX****PRT****TPT****PLP****TPP****ZPY****BLY****PLY****TPY****HTR****H60****T60****P60****B60****H67****B67****T67****P67****T70****H74****T74****B74****H70****Materials:**

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Trivex

Trivex

Trivex

Trivex

Trivex

Trivex

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Plastic Ultra High Index (> 1.67)

Plastic Ultra High Index (> 1.67)

Plastic Ultra High Index (> 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design**PAL Design on Back side****Fitting cross location:** 4mm above 180 line**Rec. minimum fitting height:**

14mm

15mm

16mm

17mm

18mm

19mm

20mm


Several

P Digitally Balanced Design for Everyday Use. Also available with Blue Blocking filter.**Available as:** Clear; Photochromic; Polarized**Available in:** US and Canada




M	
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
H53	Trivex
T53	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
HTR	Plastic High Index (1.60)
H60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)
T74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
15mm
16mm
17mm
18mm
19mm
20mm
P Advanced two-sided digitally well-balanced design that is customized and optimized using Position of Wear values (actual or standard)
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
B50	Plastic (CR39 and 1.50)
H53	Trivex
B53	Trivex
T53	Trivex
PTX	Trivex
PRT	Trivex
TPT	Trivex
PLP	Polycarbonate
TPP	Polycarbonate
ZPY	Polycarbonate
BLY	Polycarbonate
PLY	Polycarbonate
TPY	Polycarbonate
HTR	Plastic High Index (1.60)
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
B60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
B67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T70	Plastic Ultra High Index (> 1.67)
H74	Plastic Ultra High Index (> 1.67)
T74	Plastic Ultra High Index (> 1.67)
B74	Plastic Ultra High Index (> 1.67)
H70	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
15mm
16mm
17mm
18mm
19mm
20mm
Several
P Design Highlighting Distance Vision for More Active Patients. Individually compensated and optimized using Position of Wear Measurements (actual or average). Also available with Blue Blocking filter.
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
B50	Plastic (CR39 and 1.50)
H53	Trivex
B53	Trivex
T53	Trivex
PTX	Trivex
PRT	Trivex
TPT	Trivex
PLP	Polycarbonate
TPP	Polycarbonate
ZPY	Polycarbonate
BLY	Polycarbonate
PLY	Polycarbonate
TPY	Polycarbonate
HTR	Plastic High Index (1.60)
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
B60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
B67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T70	Plastic Ultra High Index (> 1.67)
H74	Plastic Ultra High Index (> 1.67)
T74	Plastic Ultra High Index (> 1.67)
B74	Plastic Ultra High Index (> 1.67)
H70	Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
P Designs are +0.50 and +0.75 adds.
Relieves visual fatigue symptoms for wearers ages 18 to 45. Also available with Blue Blocking filter.
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
B50	Plastic (CR39 and 1.50)
B53	Trivex
H53	Trivex
T53	Trivex
PTX	Trivex
PRT	Trivex
TPT	Trivex
PLP	Polycarbonate
TPP	Polycarbonate
BLY	Polycarbonate
ZPY	Polycarbonate
PLY	Polycarbonate
TPY	Polycarbonate
B60	Plastic High Index (1.60)
HTR	Plastic High Index (1.60)
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
B67	Plastic High Index (1.66 - 1.67)
H67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
13mm
14mm
15mm
16mm
17mm
18mm
Several
P Digitally customized and optimized for any frame and fine-tuned compensation with Position of Wear values (standard or actual). Also available with Blue Blocking filter.
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
H53	Trivex
T53	Trivex
PRT	Trivex
TPT	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
ZPY	Polycarbonate
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
18mm
P Progressive specifically designed for Wrap Frames. Individually compensated and optimized using Position of Wear values (actual or average).
Available as: Clear; Photochromic; Polarized
Available in: US Only



M	Materials:
H53	Trivex
T53	Trivex
PRT	Trivex
PRT	Trivex
TPT	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
ZPY	Polycarbonate
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)

PAL Design on Back side
Fitting cross location: None
Rec. minimum fitting height:
P SV F.H. within 5 mm of O.C. Single Vision specifically designed for Wrap Frames. Individually compensated and optimized using Position of Wear values (actual or average).
Available as: Clear; Photochromic; Polarized
Available in: US Only




M	Materials:
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
B50	Plastic (CR39 and 1.50)
H53	Trivex
B53	Trivex
T53	Trivex
PTX	Trivex
PRT	Trivex
TPT	Trivex
PLP	Polycarbonate
TPP	Polycarbonate
ZPY	Polycarbonate
BLY	Polycarbonate
PLY	Polycarbonate
TPY	Polycarbonate
HTR	Plastic High Index (1.60)
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
B60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
B67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T70	Plastic Ultra High Index (> 1.67)
H74	Plastic Ultra High Index (> 1.67)
T74	Plastic Ultra High Index (> 1.67)
B74	Plastic Ultra High Index (> 1.67)
H70	Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
14mm
18mm
P Digitally design for reducing eye fatigue during extended computer/task specific working conditions. Also available with Blue Blocking filter.
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Eye Kraft Optical, Inc. | DLS Advanced PAL

215




M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Eye Kraft Optical, Inc. | DLS Anti-Fatigue

216



M

1.50
1.53
1.58
1.60
1.66
1.74

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only
Manufacturer Tech Specs: [Click Here](#)

Eye Kraft Optical, Inc. | DLS Camber Near

217



M


1.50
1.53
1.58
1.60
1.66
1.74

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only
Manufacturer Tech Specs: [Click Here](#)

Eye Kraft Optical, Inc. | DLS PAL

218



M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Eye Kraft Optical, Inc. | DLS Pro PAL

219



M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

FEA Industries, Inc. | Constitution HD

220



M

no symbol
I
P
L

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Polycarbonate
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

FEA Industries, Inc. | Constitution Short HD

221



M

no symbol
I
P
L

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Polycarbonate
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

FEA Industries, Inc. | Eagle

222



M

no symbol
I
Δ
P
6
L
7
no symbol
6
7

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)
Glass (1.523)
Glass Mid-Index (> 1.53)
Glass High Index (> 1.60)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass High Index (> 1.60)
7	Glass Ultra High Index (> 1.70)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass High Index (> 1.60)
7	Glass Ultra High Index (> 1.70)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass High Index (> 1.60)
7	Glass Ultra High Index (> 1.70)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic Ultra High Index (> 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass High Index (> 1.60)
7	Glass Ultra High Index (> 1.70)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass High Index (> 1.60)
7	Glass Ultra High Index (> 1.70)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass High Index (> 1.60)
7	Glass Ultra High Index (> 1.70)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada




M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass Mid-Index (> 1.53)
7	Glass High Index (> 1.60)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 10mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass Mid-Index (> 1.53)
7	Glass High Index (> 1.60)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 8mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada




[M]	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)
no symbol	Glass (1.523)
6	Glass High Index (> 1.60)
7	Glass Ultra High Index (> 1.70)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



[M]	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic High Index (1.66 - 1.67)
no symbol	Glass (1.523)
6	Glass High Index (> 1.60)
7	Glass Ultra High Index (> 1.70)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada




[M]	Materials:
	Plastic (CR39 and 1.50)
	Trivex
	Polycarbonate
	Plastic High Index (1.60)
	Plastic High Index (1.66 - 1.67)
	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm 
Available as: Clear; Photochromic; Polarized
Available in: US Only

GSRx | epik DRV 50

234




Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Available as: Clear; Photochromic; Polarized
Available in: US Only

GSRx | epik iBifocal

235



Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Available as: Clear; Photochromic; Polarized
Available in: US Only

GSRx | epik Office HD

236




Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm 
Available as: Clear; Photochromic;
Available in: US Only

GSRx | epik Standard

237




Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm 
Available as: Clear; Photochromic; Polarized
Available in: US Only

GSRx | epik Variable

238




Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm 
Available as: Clear; Photochromic; Polarized
Available in: US Only

GSRx | epik Variable Device

239



M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm 
Available as: Clear; Photochromic; Polarized
Available in: US Only

Hoya Vision Care - Dallas | Amplitude

240



M


1
X
P
3
5

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 19mm
Available as: Clear; Photochromic;
Available in: US and Canada

Hoya Vision Care - Dallas | Amplitude BKS

241



M

1
X
P
3
5
4

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 19mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Hoya Vision Care - Dallas | Amplitude iQ

242



M

1
X
P
3
5

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 19mm
Available as: Clear; Photochromic;
Available in: US and Canada

Hoya Vision Care - Dallas

Amplitude Mini

243



M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Hoya Vision Care - Dallas

Amplitude Mini BKS

244




M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)
4	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Hoya Vision Care - Dallas

Amplitude Mini iQ

245




M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Hoya Vision Care - Dallas

Array

246




M	Materials:
1	Plastic (CR39 and 1.50)
U	Plastic Mid-Index (1.52 - 1.58)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)
4	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: Several
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Hoya Vision Care - Dallas

Array Wrap

247



M	Materials:
1	Plastic (CR39 and 1.50)
U	Plastic Mid-Index (1.52 - 1.58)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US Only

Hoya Vision Care - Dallas

Hoyalux GP

248



M	Materials:
1	Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US Only

Hoya Vision Care - Dallas

Hoyalux GPW

249




M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
3	Plastic High Index (1.60)
4	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US Only

Hoya Vision Care - Dallas

Hoyalux Tact (EP40 and 60)

250




M	Materials:
1	Plastic (CR39 and 1.50)
3	Plastic High Index (1.60)

Near Variable Focus
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada

Hoya Vision Care - Dallas

ID InStyle

251




M	Materials:
X	Trivex
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)
4	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Hoya Vision Care - Dallas

ID Lifestyle

252



M

1

X

P

3

5

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Custom Progressive Design

PAL Design on Front and Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 18mm


Available as: Clear; Photochromic; Polarized

Available in: US Only

Hoya Vision Care - Dallas

iD LifeStyle 2 Clarity

253



M

1

X

P

3

5

7

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Front and Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 18mm


Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Hoya Vision Care - Dallas

ID LifeStyle 2 Clarity cd

254



M

1

X

P

3

5

7

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Front and Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 14mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Hoya Vision Care - Dallas

ID LifeStyle 2 Harmony

255



M

1

X

P

3

5

7

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Front and Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 18mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada

Hoya Vision Care - Dallas

ID LifeStyle 2 Harmony cd

256



M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Hoya Vision Care - Dallas

ID Lifestyle cd

257



M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Hoya Vision Care - Dallas

ID Mystyle

258



M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic;
Available in: US and Canada

Hoya Vision Care - Dallas

Summit cd

259



M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)
4	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Hoya Vision Care - Dallas

Summit cd iQ

260



M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US Only

Hoya Vision Care - Dallas

Summit ecp

261



M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)
4	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Hoya Vision Care - Dallas

Summit ecp iQ

262



M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Hoya Vision Care - Dallas

TACT 40 BKS and TACT 60 BKS

263




M	Materials:
1	Plastic (CR39 and 1.50)
X	Trivex
P	Polycarbonate
3	Plastic High Index (1.60)
5	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Interstate Optical | Intellego 3

264




Materials:

PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height:
13mm
Available as:
Available in: US and Canada

Interstate Optical | Intellego 3 14

265



Materials:

PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height:
13mm
Available as:
Available in: US and Canada

Interstate Optical | Intellego 3 16

266




Materials:

PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height:
13mm
Available as:
Available in: US and Canada

Interstate Optical | Intellego 3 18

267



Materials:

PAL Design on Front side
Fitting cross location: On 180 line
Rec. minimum fitting height:
13mm
Available as:
Available in: US and Canada

IOT | inMotion

268



Materials:

PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
Available as:
Available in: US and Canada

KBco, Inc. | EOS

269




Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear Polarized
Available in: US and Canada

KBco, Inc. | EOS Wrap

270



Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear Polarized
Available in: US and Canada

KBco, Inc. | iRx Pro

271



Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

KBco, Inc. | iRx Short

272



Materials:
Plastic Mid-Index (1.52 - 1.58)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear Polarized
Available in: US and Canada

KBco, Inc. | Xplore

273



Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Filters:
Polarized

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Polarized
Available in: US and Canada

L.B.I. | Ce-Tru Normal Corridor

274



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
22mm
Available as: Clear
Available in: US and Canada

L.B.I. | Ce-Tru Polarized Progressive275



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
22mm
Available as: Polarized
Available in: US Only

L.B.I. | Ce-Tru Short Corridor276



Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear
Available in: US and Canada

L.B.I. | Ce-Tru Short Corridor Photochr277



Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US and Canada


L.B.I. | Ce-Tru SunSensor278



Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
22mm
Available as: Clear; Photochromic;
Available in: US and Canada

Landon Lens | Channel 14 Plastic II279




Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear
Available in: US and Canada

Landon Lens | Channel 14 Poly II

280




Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Landon Lens | Channel 14 Variations Gray & B

281



Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 3mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US and Canada

Landon Lens | Channel 14 Variations Poly Pho

282




Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
14mm
Available as: ; Photochromic;
Available in: US and Canada

Landon Lens | Mvp Platinum Plastic II

283



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada

Landon Lens | Variations


284



Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: ; Photochromic;
Available in: US and Canada


Lenscrafters | Blue IQ Variable Progressive285



Materials:

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
11mm
Available as:
Available in: US Only

MacLens Optical | PolySensus286



Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US and Canada


MacLens Optical | Promult287



Materials:
Plastic (CR39 and 1.50)
Glass (1.523)
Glass Ultra High Index (> 1.70)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada


Maui Jim | Maui Passport Progressive288



Materials:
Polycarbonate
Plastic High Index (1.60)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Several
Available as: Polarized
Available in: US and Canada

Melibrad | Melibrad Progressive289



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Polarized
Available in: US and Canada

Melibrad | Polar-Ray Progressive290



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
22mm
Available as: ; Photochromic;
Available in: US and Canada

Melibrad | Polar-Ray Progressive291



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Polarized
Available in: US and Canada

Nassau Lens Company | Nalco Progressive292



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
19mm
Available as: Clear; Photochromic;
Available in: US Only

Nexus Vision Group | Proprius DS293



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear
Available in: US Only

Nikon Optical Canada | Balance Digital 12294



Materials:
Plastic (CR39 and 1.50)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic;
Available in: Canada Only

Nikon Optical Canada | Balance Digital 14

295



Materials:
Plastic (CR39 and 1.50)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic;
Available in: Canada Only

Nikon Optical Canada | MOVE

296




Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: Canada Only

Nikon Optical Canada | MOVE Short

297



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: Canada Only

Nikon Optical Canada | Nikon Go

298



Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: Canada Only

Nikon Optical Canada | Nikon Go Digital

299



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: Canada Only

Nikon Optical Canada | Nikon i300



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear
Available in: Canada Only

Nikon Optical Canada | Nikon i Digital301



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear
Available in: Canada Only

Nikon Optical Canada | Nikon Online302



Materials:
Plastic (CR39 and 1.50)

Near Variable Focus
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear
Available in: US and Canada

Nikon Optical Canada | Nikon W303



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear
Available in: US and Canada

Nikon Optical Canada | Presio Power 12/13/14304



Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic;
Available in: US and Canada

Nikon Optical Canada | SeeMax Power 12/13/14

305



Materials:
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic;
Available in: US and Canada

Nikon Optical Canada | SeeMaxProgressive 12

306



Materials:
Plastic High Index (1.60)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic;
Available in: Canada Only

Nikon Optical Canada | SeeMaxProgressive 14

307




Materials:
Plastic High Index (1.60)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic;
Available in: Canada Only

Nikon Optical USA | Move Digital 13

308




Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Nikon Optical USA | Move Digital 15

309



Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Nikon Optical USA | Nikon Online

310



M

PK

NC

NC

Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Near Variable Focus
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear
Available in: US and Canada

Nikon Optical USA | Nikon W

311



M

PC

NC

NC

N5

Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear
Available in: US and Canada

Nikon Optical USA | Nikon W Short

312



M

PC

NC

NC

N5

Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic;
Available in: US Only

Nikon Optical USA | Presio Digital i Short

313



M

Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic;
Available in: US Only

Nikon Optical USA | Presio i

314



M

PC

NC

NC

N5

Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear
Available in: US Only

Nikon Optical USA | Presio i Digital

315



Materials:
Plastic (CR39 and 1.50)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US Only

Nikon Optical USA | Presio i Digital 12

316



Nb
NB
Nb
NB
N5B

Materials:
Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Nikon Optical USA | Presio i Digital 14

317



Nb
NB
Nb
NB
N5B

Materials:
Plastic (CR39 and 1.50)

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Nikon Optical USA | Presio i Digital Short

318



Materials:
Plastic High Index (1.60)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Nikon Optical USA | Presio i Short

319



PC
NC
NC
N5

Materials:
Plastic (CR39 and 1.50)

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US Only

Nikon Optical USA | Presio Power 12/13/14

320



M	Materials:
N1	Plastic (CR39 and 1.50)
N3	Plastic High Index (1.60)
N4	Plastic High Index (1.66 - 1.67)
N5	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic;
Available in: US and Canada

Nikon Optical USA | SeeMax Power 12/13/14

321



M	Materials:
N3	Plastic High Index (1.60)
N4	Plastic High Index (1.66 - 1.67)
N5	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic;
Available in: US and Canada

PFO Global | Acuity 9/11/13

322



M	Materials:
○	Plastic (CR39 and 1.50)
◇	Polycarbonate
△	Plastic High Index (1.60)
□	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

PFO Global | identity home 9/11/13

323



M	Materials:
○	Plastic (CR39 and 1.50)
◇	Polycarbonate
△	Plastic High Index (1.60)
□	Plastic High Index (1.66 - 1.67)
⬠	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

PFO Global | identity play 9/11/13

324



M




Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

PFO Global | identity pro

325



M



Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 15mm
Available as: Clear
Available in: US and Canada

PFO Global | identity smart 9/11/13

326



M




Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

PFO Optima LLC | HyperPro

327



M



Materials:
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear
Available in: US and Canada

PFO Optima LLC | Resolution Response

328



Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 5mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Polycore Optical | Futureise

329



Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Polycore Optical | Micro

330



Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic;
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
20mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	
P	
PRY	
TPR	
TRN	
H53	
PRT	
T53	
PLP	
PLY	
TPP	
TPY	
H60	
P60	
T60	
H67	
P67	
T67	
H74	

Materials:	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Trivex	
Trivex	
Trivex	
Polycarbonate	
Polycarbonate	
Polycarbonate	
Polycarbonate	
Plastic High Index (1.60)	
Plastic High Index (1.60)	
Plastic High Index (1.60)	
Plastic High Index (1.66 - 1.67)	
Plastic High Index (1.66 - 1.67)	
Plastic High Index (1.66 - 1.67)	
Plastic Ultra High Index (> 1.67)	

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
19mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M

F50

P

PRY

TPR

TRN

H53

PRT

T53

PLP

PLY

TPP

TPY

H60

P60

T60

H67

P67

T67

H74

Materials:

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Trivex

Trivex

Trivex

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: On 180 line

Rec. minimum fitting height:

15mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada



M

F50

P

PRY

TPR

TRN

H53

PRT

T53

PLP

PLY

TPP

TPY

H60

P60

T60

H67

P67

T67

H74

Materials:

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Trivex

Trivex

Trivex

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side

Fitting cross location: On 180 line

Rec. minimum fitting height:

17mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
TRN	Plastic (CR39 and 1.50)
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
TRN	Plastic (CR39 and 1.50)
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



TRN
F50
P
PRY
TPR
H53
T53
PRT
PLP
PLY
TPP
TPY
H60
P60
T60
H67
P67
T67
H74

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 20mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
TRN	Plastic (CR39 and 1.50)
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLP	Polycarbonate
PLY	Polycarbonate
TPP	Polycarbonate
TPY	Polycarbonate
H60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
H67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	
P	
TRN	
PRY	
TPR	
H53	
PRT	
T53	
PLY	
TPY	
PLP	
TPP	
H60	
T60	
P60	
T67	
P67	
H67	
H74	

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
H67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
H67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
20mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	Plastic (CR39 and 1.50)
T53	Plastic (CR39 and 1.50)
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
H67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	
P	
TRN	
PRY	
TPR	
H53	
PRT	
T53	
PLY	
TPY	
PLP	
TPP	
H60	
T60	
P60	
T67	
P67	
H67	
H74	

Materials:	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Trivex	
Trivex	
Trivex	
Polycarbonate	
Polycarbonate	
Polycarbonate	
Polycarbonate	
Plastic High Index (1.60)	
Plastic High Index (1.60)	
Plastic High Index (1.60)	
Plastic High Index (1.66 - 1.67)	
Plastic High Index (1.66 - 1.67)	
Plastic High Index (1.66 - 1.67)	
Plastic Ultra High Index (> 1.67)	

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
F50	
P	
TRN	
PRY	
TPR	
H53	
PRT	
T53	
PLY	
TPY	
PLP	
TPP	
H60	
T60	
P60	
T67	
P67	
H67	
H74	

Materials:	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Trivex	
Trivex	
Trivex	
Polycarbonate	
Polycarbonate	
Polycarbonate	
Polycarbonate	
Plastic High Index (1.60)	
Plastic High Index (1.60)	
Plastic High Index (1.60)	
Plastic High Index (1.66 - 1.67)	
Plastic High Index (1.66 - 1.67)	
Plastic High Index (1.66 - 1.67)	
Plastic Ultra High Index (> 1.67)	

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
P	
TRN	
PRY	
TPR	
F50	
H53	
PRT	
T53	
PLY	
TPY	
PLP	
TPP	
H60	
T60	
P60	
T67	
P67	
H67	
H74	

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
20mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
P	
TRN	
PRY	
TPR	
F50	
H53	
PRT	
T53	
PLY	
TPY	
PLP	
TPP	
H60	
T60	
P60	
T67	
P67	
H67	
H74	

Materials:	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Trivex	
Trivex	
Trivex	
Polycarbonate	
Polycarbonate	
Polycarbonate	
Polycarbonate	
Plastic High Index (1.60)	
Plastic High Index (1.60)	
Plastic High Index (1.60)	
Plastic High Index (1.66 - 1.67)	
Plastic High Index (1.66 - 1.67)	
Plastic High Index (1.66 - 1.67)	
Plastic Ultra High Index (> 1.67)	

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
F50	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
H67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
F50	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
H67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
P	Plastic (CR39 and 1.50)
TRN	Plastic (CR39 and 1.50)
PRY	Plastic (CR39 and 1.50)
TPR	Plastic (CR39 and 1.50)
F50	Plastic (CR39 and 1.50)
H53	Trivex
PRT	Trivex
T53	Trivex
PLY	Polycarbonate
TPY	Polycarbonate
PLP	Polycarbonate
TPP	Polycarbonate
H60	Plastic High Index (1.60)
T60	Plastic High Index (1.60)
P60	Plastic High Index (1.60)
T67	Plastic High Index (1.66 - 1.67)
P67	Plastic High Index (1.66 - 1.67)
H67	Plastic High Index (1.66 - 1.67)
H74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M

P

TRN

PRY

TPR

F50

H53

PRT

T53

PLY

TPY

PLP

TPP

H60

T60

P60

T67

P67

H67

H74

Materials:

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Trivex

Trivex

Trivex

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side


Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

19mm

Available as: Clear; Photochromic; Polarized

Available in: US and Canada



M

50

56

53

58

60

67

74

Materials:

Plastic (CR39 and 1.50)

Plastic Mid-Index (1.52 - 1.58)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive

PAL Design on Back side


Fitting cross location: On 180 line

Rec. minimum fitting height:

Several

Available as: Clear; Photochromic; Polarized

Available in: Canada Only




M	Materials:
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
58	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only




M	Materials:
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
58	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only




M	Materials:
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
58	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic;
Available in: Canada Only



M	
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
58	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only



M	
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
58	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only



M	
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
58	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Riverside Opticalab Ltd. | Persona Quattro

365



M	Materials:
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Riverside Opticalab Ltd. | Persona Quattro Sport

366



M	Materials:
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
58	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Riverside Opticalab Ltd. | Signature HD

367



M	Materials:
50	Plastic (CR39 and 1.50)
56	Plastic Mid-Index (1.52 - 1.58)
53	Trivex
58	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

Robertson Optical Laboratories, Inc.

COZE ADL ADVANCED

368



M

50

54

58

16

67

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

Several


Available as: Clear; Photochromic; Polarized

Available in: US Only

Robertson Optical Laboratories, Inc.

COZE CASATA

369



M

50

54

58

16

67

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

Several


Available as: Clear; Photochromic; Polarized

Available in: US Only

Robertson Optical Laboratories, Inc.

COZE HDW

370



M

50

54

58

16

67

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

Several


Available as: Clear; Photochromic; Polarized

Available in: US Only

Robertson Optical Laboratories, Inc.

COZE SOFT

371



M

50

54

58

16

67

Materials:

Plastic (CR39 and 1.50)

Trivex

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

Several

Available as: Clear; Photochromic; Polarized

Available in: US Only

Robertson Optical Laboratories, Inc. | COZE WRAP

372



M	Materials:
50	Plastic (CR39 and 1.50)
54	Trivex
58	Polycarbonate
16	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US Only

Robertson Optical Laboratories, Inc. | Vista

373



M	Materials:
50	Plastic (CR39 and 1.50)
54	Trivex
58	Polycarbonate
16	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US Only

Robertson Optical Laboratories, Inc. | VISTA CP

374



M	Materials:
50	Plastic (CR39 and 1.50)
54	Trivex
58	Polycarbonate
16	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US Only

Robertson Optical Laboratories, Inc. | VISTA WS

375




M	Materials:
50	Plastic (CR39 and 1.50)
54	Trivex
58	Polycarbonate
16	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US Only

Robertson Optical Laboratories, Inc. | XP HD

376



M	Materials:
50	Plastic (CR39 and 1.50)
54	Trivex
58	Polycarbonate
16	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US Only

Rodenstock Canada | ClassicLife Free XS. M or L

377




M	Materials:
50	Plastic (CR39 and 1.50)
54	Plastic Mid-Index (1.52 - 1.58)
60	Plastic High Index (1.60)
67	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | Impression 2 XS, M, or L

378




M	Materials:
5	Plastic (CR39 and 1.50)
54	Plastic Mid-Index (1.52 - 1.58)
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)
60	Glass High Index (> 1.60)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | Impression Ergo 2 Book

379



M	Materials:
60	Plastic High Index (1.60)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
14mm
Available as: Clear
Available in: US and Canada

Rodenstock Canada | Impression Ergo 2 PC or Room

380



M


60

Materials:
Plastic High Index (1.60)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada

Rodenstock Canada | Impression Ergo FreeSign 2

381



M


60

Materials:
Plastic High Index (1.60)

Near Variable Focus
PAL Design on Front and Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
14mm
Available as: Clear
Available in: US and Canada

Rodenstock Canada | Impression FreeSign 2

382



M

60

67

74

Materials:
Plastic High Index (1.60)


Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | Impression Hyperop

383



M


7

Materials:
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | Impression Hyperop XS

384



M

7

Materials:
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | Impression Myop

385



M

7

Materials:
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | Impression Myop XS

386



M


7

Materials:
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | Impression Sport 2

387



M

54

6

Materials:
Plastic Mid-Index (1.52 - 1.58)
Plastic High Index (1.60)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | Multigressiv MyView 2 XS, M, L

388



M

50

54

60

67

74

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Rodenstock Canada | MyView Ergo 2 Book

389



M

50

60

Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear
Available in: US and Canada

Rodenstock Canada

MyView Ergo 2 PC or Room

390



M

50

60

Materials:

Plastic (CR39 and 1.50)

Plastic High Index (1.60)

Near Variable Focus

PAL Design on Back side

Fitting cross location: On 180 line

Rec. minimum fitting height: 18mm

Available as: Clear

Available in: US and Canada

Rodenstock Canada

Nexyma Computer (80)

391



M

5

Materials:

Plastic (CR39 and 1.50)

Near Variable Focus

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 20mm


Available as: Clear

Available in: US and Canada

Rodenstock Canada

Nexyma Reading (40)

392



M

5

Materials:

Plastic (CR39 and 1.50)

Near Variable Focus

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 13mm

Available as: Clear

Available in: US and Canada

Rodenstock Canada

Progressiv ClassicLife

393



M

5

6

Materials:

Plastic (CR39 and 1.50)

Plastic High Index (1.60)

Standard Progressive

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 18mm

Available as: Clear

Available in: US and Canada

Rodenstock Canada

Progressiv ClassicLife XS

394



M

5

6

Materials:

Plastic (CR39 and 1.50)

Plastic High Index (1.60)

Standard Progressive

PAL Design on Front side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height: 14mm

Available as: Clear

Available in: US and Canada

Rodenstock Canada

Progressiv Life 2

395



M

Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear
Available in: US and Canada

Rodenstock Canada

Progressiv Life 2 XS

396



M

P

Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear
Available in: US and Canada

Rodenstock Canada

Progressiv Purelife

397



M

5

Materials:
Plastic (CR39 and 1.50)

5

Plastic Mid-Index (1.52 - 1.58)

6

Plastic High Index (1.60)

7

Plastic High Index (1.66 - 1.67)

6

Glass High Index (> 1.60)

8

Glass Ultra High Index (> 1.70)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada

Progressiv Purelife XS

398



M

5

Materials:
Plastic (CR39 and 1.50)

5

Plastic Mid-Index (1.52 - 1.58)

6

Plastic High Index (1.60)

7

Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Rodenstock Canada | ProgressivPureLife Free XS,M,L

399



M

50

54

1.53

60

67

Materials:

- Plastic (CR39 and 1.50)
- Plastic Mid-Index (1.52 - 1.58)
- Trivex
- Plastic High Index (1.60)
- Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Seiko Optical Products | AF2

400



M

A

C

S

Materials:

- Plastic (CR39 and 1.50)
- Plastic High Index (1.60)
- Plastic High Index (1.66 - 1.67)

Standard Progressive

PAL Design on Front and Back side
Fitting cross location: 3mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear
Available in: US and Canada

Seiko Optical Products | DC mini

401



M

Materials:

- Polycarbonate

Standard Progressive

PAL Design on Front and Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 19mm
Available as: Clear
Available in: US Only

Seiko Optical Products | PCWide

402



M

□

□

▽

□

△

□

Materials:


- Plastic (CR39 and 1.50)
- Trivex
- Polycarbonate
- Plastic High Index (1.60)
- Plastic High Index (1.66 - 1.67)
- Plastic Ultra High Index (> 1.67)

Near Variable Focus

PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic;
Available in: US Only

Seiko Optical Products | Perfas Premier II

403



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Seiko Optical Products | Perfas Prominent

404



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Seiko Optical Products | Proceed II Short

405




M	Materials:
	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: US Only

Seiko Optical Products | Proceed III Super Short

406



M	Materials:
	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front and Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic;
Available in: US Only

Seiko Optical Products | Succeed

407



M

Materials:
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 17mm
Available as: Clear
Available in: US and Canada

Seiko Optical Products | Succeed WS

408



M

Materials:
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 15mm
Available as: Clear
Available in: US and Canada

Seiko Optical Products | Supercede

409



M

Materials:
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: Several
Available as: Clear
Available in: US and Canada

Seiko Optical Products | Supercede II

410



M

50

53

59

60

67

74

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Seiko Optical Products | Supercede WS

411



M

Materials:
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US and Canada

Seiko Optical Products | Superior

412



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 12
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Seiko Optical Products | Supernal

413



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Seiko Optical Products | Surmount

414



M	Materials:
50	Plastic (CR39 and 1.50)
59	Polycarbonate
	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only


Seiko Optical Products | Surmount Ws

415



M	Materials:
50	Plastic (CR39 and 1.50)
53	Trivex
59	Polycarbonate
60	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only



M


15
TR
X5
PZ
56
TX
TTX
XX
VX

PT
PX
PV
PD
PC
16
T6
X6
PZ6
PZ67
67
T67
X67

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Trivex
Trivex
Trivex
Trivex

Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as:
Available in: US Only



M

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Shamir Insight Inc. | Attitude III Sport

418



Materials:

Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

18mm

Available as: Clear; Photochromic; Polarized

Available in: US Only

Shamir Insight Inc. | Autograph II -Fixed & Var. Des

419



Materials:

Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

11mm

Available as: Clear

Available in: US Only

Shamir Insight Inc. | Autograph II Attitude

420



Materials:

Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Custom Progressive Design

PAL Design on Back side

Fitting cross location: 4mm above 180 line

Rec. minimum fitting height:

15mm

18mm

Available as: Clear; Photochromic; Polarized

Available in: US Only

Shamir Insight Inc. | Autograph II Office

421



Materials:

Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.60)

Near Variable Focus

PAL Design on Back side

Fitting cross location: On 180 line

Rec. minimum fitting height:


16mm

Available as: Clear

Available in: US Only

Shamir Insight Inc. | Autograph III Fixed & VAR.

422



M

PZ
15
TR
X5
56
TX
PC
16
74

Materials:	Custom Progressive Design PAL Design on Back side Fitting cross location: 4mm above 180 line Rec. minimum fitting height: 11mm 13mm 15mm 18mm Available as: Clear; Photochromic; Polarized Available in: US Only
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic (CR39 and 1.50)	
Plastic Mid-Index (1.52 - 1.58)	
Trivex	
Polycarbonate	
Plastic High Index (1.60)	
Plastic Ultra High Index (> 1.67)	

Shamir Insight Inc. | Autograph Intelligence Fixed & VARI.

423



M

Materials:	Custom Progressive Design PAL Design on Back side Fitting cross location: 4mm above 180 line Rec. minimum fitting height: 11mm 13mm 15mm 18mm Available as: Clear; Photochromic; Polarized Available in: US and Canada
Plastic (CR39 and 1.50)	
Trivex	
Polycarbonate	
Plastic High Index (1.60)	
Plastic High Index (1.66 - 1.67)	
Plastic Ultra High Index (> 1.67)	

Shamir Insight Inc. | Autograph II+ Fixed & VARI.

424



M

Materials:	Custom Progressive Design PAL Design on Back side Fitting cross location: 4mm above 180 line Rec. minimum fitting height: 11mm 13mm 15mm 18mm Available as: Clear; Photochromic; Polarized Available in: US Only
Plastic (CR39 and 1.50)	
Trivex	
Polycarbonate	
Plastic High Index (1.60)	
Plastic High Index (1.66 - 1.67)	
Plastic Ultra High Index (> 1.67)	

Shamir Insight Inc. | Autograph Relax

425




M

Materials:	PAL Design on Back side Fitting cross location: On 180 line Rec. minimum fitting height: 16mm Available as: Clear; Photochromic; Polarized Available in: US Only
Plastic (CR39 and 1.50)	
Trivex	
Polycarbonate	
Plastic High Index (1.60)	
Plastic High Index (1.66 - 1.67)	
Plastic Ultra High Index (> 1.67)	

Shamir Insight Inc. | Computer

426




Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Polycarbonate
Plastic High Index (1.60)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Shamir Insight Inc. | Element

427



Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 19mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Shamir Insight Inc. | Element Short

428




Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Shamir Insight Inc. | First PAL

429



Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.60)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Shamir Insight Inc. | Golf

430



Materials:
Trivex
Polycarbonate
Plastic High Index (1.60)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Shamir Insight Inc. | Spectrum

431



M

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
16mm
19mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Shamir Insight Inc. | WorkSpace

432



M

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Polycarbonate
Plastic High Index (1.60)

Near Variable Focus
PAL Design on Back side
Fitting cross location: On 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

Shore Lens Company | Executive/ ShoreView ES

433



M

S

π

S

F



Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.60)

Filters:
Photochromic

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Shore Lens Company | ShoreView / Shoreview Mini

434



M

TTX

TTX

TV

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Trivex
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
15mm
Available as: Clear; Photochromic;
Available in: US and Canada

Shore Lens Company | Shoreview Mini ES

435




- M
- S
- S
- S
- S
- F
-
- Materials:
- Plastic (CR39 and 1.50)
- Polycarbonate
- Plastic High Index (1.60)
- Plastic High Index (1.66 - 1.67)
- Filters:
- Photochromic

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 15mm
Available as: Clear; Photochromic;
Available in: US Only

Signet Armorlite, Inc. | Crossbows Custom U

436




- M
- I
- Δ
- P
- ◇
- ⊥
- Materials:
- Plastic (CR39 and 1.50)
- Plastic Mid-Index (1.52 - 1.58)
- Trivex
- Polycarbonate
- Plastic High Index (1.60)
- Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Signet Armorlite, Inc. | Crossbows Custom V

437




- M
- I
- Δ
- P
- ◇
- ⊥
- Materials:
- Plastic (CR39 and 1.50)
- Plastic Mid-Index (1.52 - 1.58)
- Trivex
- Polycarbonate
- Plastic High Index (1.60)
- Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada


Signet Armorlite, Inc. | Crossbows Junior

438




- M
- I
- P
- ⊥
- Materials:
- Plastic (CR39 and 1.50)
- Plastic Mid-Index (1.52 - 1.58)
- Polycarbonate
- Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada




M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
L	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 16.5mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada




M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
L	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13.5mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada




M	Materials:
no symbol	Plastic (CR39 and 1.50)
Δ	Trivex
P	Polycarbonate
◇	Plastic High Index (1.60)
L	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear
Available in: US and Canada



M	Materials:
	Plastic (CR39 and 1.50)
Δ	Trivex
P	Polycarbonate
◇	Plastic High Index (1.60)
L	Plastic High Index (1.66 - 1.67)

Custom Progressive Design
PAL Design on Front and Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada




M	Materials:
no symbol	Plastic (CR39 and 1.50)
	Plastic (CR39 and 1.50)
Δ	Trivex
P	Polycarbonate
◇	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
I	Plastic (CR39 and 1.50)
	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
◇	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



M	Materials:
no symbol	Plastic (CR39 and 1.50)
P	Polycarbonate
⊥	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic;
Available in: US and Canada




M	Materials:
I	Plastic (CR39 and 1.50)
	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
◇	Plastic High Index (1.60)
⊥	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Signet Armorlite, Inc. | KODAK SoftWear

447




M	Materials:
M	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
P	Polycarbonate
L	Plastic High Index (1.66 - 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear
Available in: US and Canada

Signet Armorlite, Inc. | Kodak Unique

448




M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
6	Plastic High Index (1.60)
L	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Signet Armorlite, Inc. | KODAK Unique DRO

449




M	Materials:
I	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
◇	Plastic High Index (1.60)
L	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Signet Armorlite, Inc. | KODAK Unique DRO HD

450




M	Materials:
I	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
Δ	Trivex
P	Polycarbonate
◇	Plastic High Index (1.60)
L	Plastic High Index (1.66 - 1.67)
7	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Signet Armorlite, Inc. | Navigator FBS

451




M	Materials:
I	Plastic (CR39 and 1.50)
Δ	Plastic Mid-Index (1.52 - 1.58)
P	Trivex
⊥	Polycarbonate
7	Plastic High Index (1.66 - 1.67)
	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Signet Armorlite, Inc. | Navigator Precision

452




M	Materials:
no symbol	Plastic (CR39 and 1.50)
I	Plastic Mid-Index (1.52 - 1.58)
P	Polycarbonate
◇	Plastic High Index (1.60)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 20mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Signet Armorlite, Inc. | Navigator Short

453



M	Materials:
O	Plastic (CR39 and 1.50)
sP	Polycarbonate
⊥	Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic;
Available in: US and Canada

SOLA Optical | Compact ULTRA HD

454




M	Materials:
H	Plastic (CR39 and 1.50)
T	Trivex
P	Polycarbonate
F	Plastic High Index (1.60)
M	Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

SOLA Optical | Percepta

455



M
H
S
P
6

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Polycarbonate
Plastic High Index (1.60)
Glass High Index (> 1.60)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

SOLA Optical | SOLA HDV

456



M
H
T
P
F
M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 13mm
Available as: Clear; Photochromic; Polarized
Available in: Canada Only

SOLA Optical | SOLAOne HD

457




M
H
T
P
F
M

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

SOLA Optical | VIP Gold

458



M

Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: US and Canada

SOLA Optical | Visuality

459




M
H
S
P

Materials:
Plastic (CR39 and 1.50)
Plastic Mid-Index (1.52 - 1.58)
Polycarbonate

PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: US and Canada

SOLA Optical | XL

460



M

Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Glass (1.523)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

SOLA Optical | XL Gold

461




M

Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

SOMO Optical | EZ View Mini

462



M

Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
16mm
Available as: Clear
Available in: US and Canada

SOMO Optical | EZ View Std

463



M

Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
20mm
Available as: Clear
Available in: US and Canada

Sutherlin Optical, Inc. | Overture

464



M


*
P
6
67

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

US Optical | Advanced HD

465



M


15
TR
TX
PC
16
67
74

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

US Optical | Ultimate HD

466



M


15
TR
TX
PC
16
74

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

US Optical | Wrap HD

467



M


15
TR
TX
PC
16
74

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Vision Dynamics Lab/Optical Dynamics | Continual Focus Lens (CFL)

468



M

Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Vision Dynamics Lab/Optical Dynamics | Paradigm Progressive 469



Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 20mm
Available as: Clear; Photochromic;
Available in: US and Canada

Vision Dynamics Lab/Optical Dynamics | Paradigm Short Corridor 470



Materials:
Plastic Mid-Index (1.52 - 1.58)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic;
Available in: US and Canada

Vision Warehouse | Stealth 15 471



Materials:
Trivex

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 18mm
Available as: Clear
Available in: US and Canada

Vision-Ease Lens, Inc. | Illumina 472



Materials:
Plastic (CR39 and 1.50)
Polycarbonate


Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Vision-Ease Lens, Inc. | Novel 473



Materials:
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height: 16mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada



Materials:
Plastic (CR39 and 1.50)
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
18mm
Available as:
Available in: US and Canada



- X5

Dw

15

TR

PZ

TX

TTX

XX

PC

PPZ

PX

PD

16

T6

PZ6

X6

67

T67

PZ67

X67

74
- Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Trivex

Trivex

Trivex

Trivex

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US Only



M	
X5	Plastic (CR39 and 1.50)
Dw	Plastic (CR39 and 1.50)
15	Plastic (CR39 and 1.50)
TR	Plastic (CR39 and 1.50)
PZ	Plastic (CR39 and 1.50)
TX	Trivex
TTX	Trivex
XX	Trivex
PC	Polycarbonate
PPZ	Polycarbonate
PX	Polycarbonate
PD	Polycarbonate
16	Plastic High Index (1.60)
T6	Plastic High Index (1.60)
PZ6	Plastic High Index (1.60)
X6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
PZ67	Plastic High Index (1.66 - 1.67)
X67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Standard Progressive
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
Several
Available as: Clear; Photochromic; Polarized
Available in: US Only



M	
X5	Plastic (CR39 and 1.50)
Dw	Plastic (CR39 and 1.50)
15	Plastic (CR39 and 1.50)
TR	Plastic (CR39 and 1.50)
PZ	Plastic (CR39 and 1.50)
TX	Trivex
TTX	Trivex
XX	Trivex
PC	Polycarbonate
PPZ	Polycarbonate
PX	Polycarbonate
PD	Polycarbonate
	Polycarbonate
16	Plastic High Index (1.60)
T6	Plastic High Index (1.60)
PZ6	Plastic High Index (1.60)
X6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
PZ67	Plastic High Index (1.66 - 1.67)
X67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)


Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only



M	
X5	Plastic (CR39 and 1.50)
Dw	Plastic (CR39 and 1.50)
15	Plastic (CR39 and 1.50)
TR	Plastic (CR39 and 1.50)
PZ	Plastic (CR39 and 1.50)
TX	Trivex
TTX	Trivex
XX	Trivex
PC	Polycarbonate
PPZ	Polycarbonate
PX	Polycarbonate
PD	Polycarbonate
16	Plastic High Index (1.60)
T6	Plastic High Index (1.60)
PZ6	Plastic High Index (1.60)
X6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
PZ67	Plastic High Index (1.66 - 1.67)
X67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 12mm
Available as: Clear; Photochromic; Polarized
Available in: US Only



M	
X5	Plastic (CR39 and 1.50)
Dw	Plastic (CR39 and 1.50)
15	Plastic (CR39 and 1.50)
TR	Plastic (CR39 and 1.50)
PZ	Plastic (CR39 and 1.50)
TX	Trivex
TTX	Trivex
XX	Trivex
PC	Polycarbonate
PPZ	Polycarbonate
PX	Polycarbonate
PD	Polycarbonate
	Polycarbonate
16	Plastic High Index (1.60)
T6	Plastic High Index (1.60)
PZ6	Plastic High Index (1.60)
X6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
PZ67	Plastic High Index (1.66 - 1.67)
X67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 12mm
Available as: Clear; Photochromic; Polarized
Available in: US Only



M	
X5	Plastic (CR39 and 1.50)
Dw	Plastic (CR39 and 1.50)
15	Plastic (CR39 and 1.50)
TR	Plastic (CR39 and 1.50)
PZ	Plastic (CR39 and 1.50)
TX	Trivex
TTX	Trivex
XX	Trivex
PC	Polycarbonate
PPZ	Polycarbonate
PX	Polycarbonate
PD	Polycarbonate
	Polycarbonate
16	Plastic High Index (1.60)
T6	Plastic High Index (1.60)
PZ6	Plastic High Index (1.60)
X6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
PZ67	Plastic High Index (1.66 - 1.67)
X67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
12mm
Available as: Clear; Photochromic; Polarized
Available in: US Only

□

M

PZ

X5

Dw

15

TR

TTX

XX

TX

PC

PPZ

PX

PD

16

T6

PZ6

X6

67

T67

PZ67

X67

74

Materials:

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Plastic (CR39 and 1.50)

Trivex

Trivex

Trivex

Trivex

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Polycarbonate

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.60)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic High Index (1.66 - 1.67)

Plastic Ultra High Index (> 1.67)

Near Variable Focus

PAL Design on Back side

Fitting cross location: 2mm above 180 line

Rec. minimum fitting height:

14mm

Available as: Clear; Photochromic; Polarized

Available in: US Only

□

M

PZ	Plastic (CR39 and 1.50)
X5	Plastic (CR39 and 1.50)
Dw	Plastic (CR39 and 1.50)
15	Plastic (CR39 and 1.50)
TR	Plastic (CR39 and 1.50)
TTX	Trivex
XX	Trivex
TX	Trivex
PC	Polycarbonate
PPZ	Polycarbonate
PX	Polycarbonate
PD	Polycarbonate
	Polycarbonate
16	Plastic High Index (1.60)
T6	Plastic High Index (1.60)
PZ6	Plastic High Index (1.60)
X6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
PZ67	Plastic High Index (1.66 - 1.67)
X67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Materials:


Near Variable Focus
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
14mm
Available as: Clear; Photochromic; Polarized
Available in: US Only



M	
PZ	Plastic (CR39 and 1.50)
X5	Plastic (CR39 and 1.50)
Dw	Plastic (CR39 and 1.50)
15	Plastic (CR39 and 1.50)
TR	Plastic (CR39 and 1.50)
TTX	Trivex
XX	Trivex
TX	Trivex
PC	Polycarbonate
PPZ	Polycarbonate
PX	Polycarbonate
PD	Polycarbonate
	Polycarbonate
16	Plastic High Index (1.60)
T6	Plastic High Index (1.60)
PZ6	Plastic High Index (1.60)
X6	Plastic High Index (1.60)
67	Plastic High Index (1.66 - 1.67)
T67	Plastic High Index (1.66 - 1.67)
PZ67	Plastic High Index (1.66 - 1.67)
X67	Plastic High Index (1.66 - 1.67)
74	Plastic Ultra High Index (> 1.67)

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 12mm
Available as: Clear; Photochromic; Polarized
Available in: US Only



M

PZ
X5
Dw
15
TR
TTX
XX


TX
PC
PPZ
PX
PD

16
T6
PZ6
X6
67
T67
PZ67
X67
74

Materials:
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Plastic (CR39 and 1.50)
Trivex
Trivex
Trivex
Trivex
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate
Polycarbonate

Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: 12mm
Available as: Clear; Photochromic; Polarized
Available in: US Only



M


no symbol
TX
P
6
67
74

Materials:
Plastic (CR39 and 1.50)
Trivex
Polycarbonate
Plastic High Index (1.60)
Plastic High Index (1.66 - 1.67)
Plastic Ultra High Index (> 1.67)

Near Variable Focus
PAL Design on Back side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height: Several
Available as: Clear; Photochromic; Polarized
Available in: US Only

X-CEL Optical | Freedom ID

487



Materials:
Trivex
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

X-CEL Optical | Freedom ID

488



Materials:
Plastic (CR39 and 1.50)
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

X-CEL Optical | Freedom ID

489




Materials:
Glass (1.523)
Glass High Index (> 1.60)

Standard Progressive
PAL Design on Front side
Fitting cross location: 4mm above 180 line
Rec. minimum fitting height:
17mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Younger Optics | Adage

490



Materials:
Trivex
Polycarbonate

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
13mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Younger Optics | Image

491




Materials:
Plastic (CR39 and 1.50)
Polycarbonate
Plastic High Index (1.66 - 1.67)

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic; Polarized
Available in: US and Canada

Younger Optics | Image

492




Materials:
Trivex

Standard Progressive
PAL Design on Front side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Clear; Photochromic;
Available in: US and Canada

Zeal Optics | Zeal Digital Progressive

493



Materials:
Polycarbonate

Custom Progressive Design
PAL Design on Back side
Fitting cross location: 2mm above 180 line
Rec. minimum fitting height:
18mm
Available as: Polarized
Available in: US and Canada