

Shamir Golf™

Focusing on what's important on the green.


Every golfer has his day... you know the day we're talking about. The sun was shining, the birds were chirping, your swing was powerful and exact, and the ball did everything you wanted it to... you played the best game of your life. Now with Shamir Golf™, every day out on the green can be your best day!

Designed for your visual sweet spots.

Shamir Golf™ is a Freeform® lens designed to ensure all your patients experience their perfect game. Unlike standard progressives, the three vision zones are specifically tailored to the needs of golfers to maximize their viewing: the scorecard in their hand, the ball at their feet and the green in the distance. Shamir Golf™ provides sharp focus without distortion in these three crucial areas. And before we forget, it's the first Golf lens designed with a stable mid-range zone to ensure the ball on the ground can be seen clearly so your patients will have an accurate drive or putt.

One? Two? THREE enhanced vision zones!

Vertical Cross-Section along the Convergence Path




Near
(scorecard)



Mid-range
(golf ball)




Far
(green)



The Secret? It's all in the hips.

Well, for your lenses, the secret is in the peripheral vision. A golfer's peripheral vision is extremely important when they are lining up to take a swing. They need to be able to focus on the ball at their feet and at the same time see where they want the ball to go in the distance. This is where they utilize their peripheral zone and need it to be crystal clear. Lucky for them, Shamir Golf™ is designed specifically for wraparound frames and has a clear peripheral vision zone so your patients can keep on improving their game.

The Secret is in the Peripheral Vision.



Shamir Golf™ offers clear peripheral vision just where golfers need it most.

Take it from the experts who designed the first distortion-free solution for wrap frames, when it comes to providing those patients who enjoy spending time on the green with a progressive solution, Shamir Golf™ will give them the vision they need.

Material Availability

Material	Power Range	Add Power	Cylinder
DLC™ Trivex™ [C,TV,T7,P,TX,BZ]	-5.00 to +4.00	1.00 to 3.50	to -6.50
1.56 [C,P]	-5.50 to +5.25	1.00 to 3.50	to -6.75
1.6 SuperLite™ [C,T7,P,TX]	-6.00 to +5.25	1.00 to 3.50	to -7.00
Polycarbonate [C,TV,T7,P,TX,D,BZ]	-6.00 to +5.25	1.00 to 3.50	to -7.00
NXT [PH]	-5.00 to +4.00	1.00 to 3.50	to -6.50

POWER RANGE IS DETERMINED BY SHAMIR CERTIFIED BLANK RANGE PRISM TO 10 DIOPTERS IN EACH QUADRANT

LEGEND:

- C-Clear
- T7-Transitions® Signature7™
- PH-Photochromic
- P-Polarized
- BZ-Blue Zero™
- TV-Transitions® Vantage™
- TX-XTRActive™
- B-BluTech (Indoor & Outdoor)
- D-Drivewear

The low-down on Shamir Golf™

The benefits

- Specialty lens designed specifically for the golf enthusiast
- Designed to ensure sharp vision at 3 crucial focal points
- Utilizes As-Worn Technology™ and EyePoint Technology®

Why patients need it...

- Provides the player with a more comfortable golfing experience
- Clear peripheral viewing
- Personalized to the patient
- Designed specifically for wrap frames

Minimum Fitting Height

- 19mm

Scan for material availability:



Technical Inquiries: 888.370.0736



www.facebook.com/ShamirInsight



www.twitter.com/ShamirInsight