




# Shamir Super Centration Chart




Workspace/Computer - MFH 14mm/Top 10mm  
 Autograph II Office - MFH 16mm/Top 10mm  
 All Shamir Single Vision Lenses  
 Relax 50 • 65 • 80

Autograph Series 11 13 15 18  
 InTouch 15 18  
 Spectrum/Spectrum+ 14 16 18  
 Element 16 19  
 Genesis HD 16  
 Attitude Fashion/Attitude II 15 18  
 Attitude Sport 18  
 Golf 19  
 FirstPAL 15

11mm, 13mm, 14mm, 15mm, 16mm, 18mm, 19mm & Variable minimum fitting heights




## Shamir Super Centration Chart



Duo - MFH 15mm

### Glacier PLUS™ FSV

- 75mm ●
- 73mm ●
- 70mm ●
- 65mm ●



Contact Shamir Technical Services at 888-370-0736

