

Exhibit A: Transitions® XTRActive® New Generation Product Availability and Power Range

Design	Plastic/1.50	Trexa®***	Airwear®	Thin & Lite® 1.67	Add Power	Cyl
Varilux® X 4D	-10.00 to +6.00	x	-10.00 to +6.00	-12.00 to +9.00	+0.75 to +4.00	-6.00
Varilux® X Design Fit	-10.00 to +6.00	x	-10.00 to +6.00	-12.00 to +9.00	+0.75 to +4.00	-6.00
Varilux® X Design	-10.00 to +6.00	x	-10.00 to +6.00	-12.00 to +9.00	+0.75 to +4.00	-6.00
Varilux® Physio® W3+ eyecode™	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-14.00 to +9.00	+0.75 to +3.50, +4.00	-6.00
Varilux® Physio® W3+ Fit	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-14.00 to +9.00	+0.75 to +3.50, +4.00	-6.00
Varilux® Physio® W3+	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-14.00 to +9.00	+0.75 to +3.50, +4.00	-6.00
Varilux® Comfort® Max Fit	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Varilux® Comfort® Max	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Varilux Comfort® W2+ eyecode™	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Varilux Comfort® W2+ Fit	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Varilux Comfort® W2+	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
TruClear® XD	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-12.00 to +9.00	+0.75 to +3.50	-6.00
TruClear® XD2	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-12.00 to +9.00	+0.75 to +3.50	-6.00
TruClear® SD	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-12.00 to +9.00	+0.75 to +3.50	-6.00
TruClear® SD2	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-12.00 to +9.00	+0.75 to +3.50	-6.00
Varilux Comfort®*	-10.00 to +6.00	x	-10.00 to +6.00	x	+0.75 to +3.50	-6.00
Varilux Comfort Short™*	-10.00 to +6.00	x	-10.00 to +6.00	x	+0.75 to +3.00	-6.00
Accolade®**	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Accolade Freedom™**	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Varilux® Stylistic®	-10.00 to +7.50	x	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Varilux Physio DRx™	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50, +4.00	-6.00
Varilux Physio Short DRx™	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.00	-6.00
Varilux Comfort DRx™	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Varilux Comfort Short DRx™	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.00	-6.00
Essilor Ideal®	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Essilor Ideal® Advanced	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Essilor Ideal Advanced™ FIT	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Essilor Ideal® FIT	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Essilor Ideal® Advanced Wrap	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Essilor Ideal® Short	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.00	-6.00
Essilor Ideal® Sport	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Essilor Ideal® Computer	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Essilor Ideal® Driving	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Essilor Natural® Digital	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Adaptar® Digital	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Adaptar® Short Digital	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.00	-6.00
Brillenz®	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
Brillenz® Short	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.00	-6.00
Ovation® Digital	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.50	-6.00
SmallFit® Digital	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	+0.75 to +3.00	-6.00
Essilor® Stylistic® SV	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Eyezen™ Start	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Eyezen™+ 0	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Eyezen™+ 1	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Eyezen™+ 2	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Eyezen™+ 3	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Eyezen™+ 4	-10.00 to +7.50	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Single Vision	-10.50 to +9.00	-9.50 to +4.00	-10.00 to +8.00	-14.00 to +9.00	x	-4.00
Single Vision AS	x	x	-10.00 to +8.00	x	x	-4.00
Essilor Eyecode SV	x	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Essilor Fit™ SV	x	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Essilor Azio™ SV	x	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Essilor Single Vision 360°	x	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
TruClear® XD Single Vision	x	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
TruClear® SD Single Vision	x	-10.00 to +7.50	-10.00 to +8.00	-14.00 to +9.00	x	-6.00
Finished SV 1.50 HC (Gray)	-4.00 to +4.00	x	x	x	x	-2.00

*Use HD02 blanks, **Use SFSV blanks, ***Available in Gray and Brown