



If what you
SEE is what
you **GET**

Get MORE with **TREXA**TM



TrexaTM is a Trivex[®] based material from Essilor

5 REASONS WHY YOU WILL LOVE TREXATM OVER STANDARD PLASTIC



THINNER & 19% LIGHTER

Compared to Standard Plastic, *TREXA* has a higher refractive index (1.53 vs. 1.50) and lower density (1.11g/cm³ vs 1.32g/cm³), making it thinner and much lighter.



6X MORE IMPACT RESISTANT

Thanks to its special flexible chemical structure, *TREXA* has great impact resistance. Standard Plastic, rigid in structure, is more likely to break.



EASY TO DRILL MOUNT

TREXA's tensile strength is almost 4 times Standard Plastic's (61.2kg/f vs. 15.6kg/f), thus, a lot easier to drill and suitable for rimless frames.



BLOCKS 100% FRONT-SIDE UV TRANSMISSION

Dangerous UV light (UVA & UVB) has a wavelength from 280nm to 380nm. While Standard Plastic only blocks UV light with a wavelength up to 355nm, *TREXA* offers full protection from front-side UV transmission.



KID-FRIENDLY

Thanks to its durability and great optical quality, *TREXA* is suitable for children.

Available with
TOP BRANDS
You Trust

Transitions[®]
Adaptive Lenses[™]

Crizal[®]
No-Glare Lenses

VARILUX[®]
Progressive Lenses

Eyezen[™] +
Enhanced Single Vision Lenses



TREXATM

TREXA™ Availability with Essilor Premium Lenses

TREXA™ Availability: Varilux® & Eyezen™+ Clear Lenses

| PRODUCT | POWER RANGE | BASE CURVES | ADD POWER | MAX CYLINDER |
|----------------------------------|-----------------|------------------------------|---------------------|--------------|
| *Varilux S 4D™ | -10.00 to +6.00 | | 0.75 to 4.00D | -6.00 |
| *Varilux S Fit™ | -10.00 to +6.00 | | 0.75 to 4.00D | -6.00 |
| *Varilux S Design™ | -10.00 to +6.00 | | 0.75 to 4.00D | -6.00 |
| Varilux® Physio® W3+ | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| Varilux® Physio® W3+ Fit | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| Varilux® Physio® W3+ Eyecode™ | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| Varilux Physio Enhanced | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| Varilux Physio Enhanced Fit™ | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| Varilux Physio Enhanced Eyecode™ | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| Varilux Physio DRx™ | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| Varilux Physio Short DRx™ | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.00 | -6.00 |
| Varilux® Physio® | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -4.00 |
| Varilux Physio® Short | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -4.00 |
| Varilux Comfort DRx™ | -10.00 to +6.00 | 2.00, 3.50, 4.50, 6.25, 7.25 | 0.75 to 3.50 | -6.00 |
| Varilux Comfort Short DRx™ | -10.00 to +6.00 | 2.00, 3.50, 4.50, 6.25, 7.25 | 0.75 to 3.00 | -6.00 |
| *Eyezen™+ | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | | -6.00 |

TREXA™ Availability: Varilux® & Eyezen™+ with Transitions®

| | Transitions® Signature™ VII | | | Transitions® XTRActive® | | |
|----------------------------------|-----------------------------|---|----|-------------------------|---|----|
| | G | B | GG | G | B | GG |
| *Varilux S 4D™ | • | • | • | • | • | • |
| *Varilux S Fit™ | • | • | • | • | • | • |
| *Varilux S Design™ | • | • | • | • | • | • |
| Varilux® Physio® W3+ | • | • | • | • | • | • |
| Varilux® Physio® W3+ Fit | • | • | • | • | • | • |
| Varilux® Physio® W3+ Eyecode™ | • | • | • | • | • | • |
| Varilux Physio Enhanced | • | • | • | • | • | • |
| Varilux Physio Enhanced Fit™ | • | • | • | • | • | • |
| Varilux Physio Enhanced Eyecode™ | • | • | • | • | • | • |
| Varilux Physio DRx™ | • | • | • | • | • | • |
| Varilux Physio Short DRx™ | • | • | • | • | • | • |
| Varilux® Physio® | • | • | • | • | • | • |
| Varilux Comfort DRx™ | • | • | • | • | • | • |
| Varilux Comfort Short DRx™ | • | • | • | • | • | • |
| *Eyezen™+ | • | • | • | • | • | • |

*Available Q1 2017

Trexa™ is available with Crizal® coating. Please refer to Crizal.com for Trexa Crizal availability information. Trexa & Trexa with Transitions® are also available on other Essilor Progressive and Single Vision lenses. For further availability information, please refer to EssilorLenses.com.

G: Gray - **B:** Brown - **GG:** Graphite Green

©2016 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries. Transitions and the swirl are registered trademarks of Transitions Optical, Inc., used under license by Transitions Optical Limited. Photochromic performance is influenced by temperature, UV exposure, and lens material. Trivex® is a registered trademark of PPG Industries Ohio, Inc. LMAT000035 SHK/ESCB 8/16

