




Accolade® and Accolade Freedom™ lenses...
The PALs that adapt to your patients in an instant


- **Harmonix™ Technology** – the progressive lens design customized to individual eye shape
- **9 out of 10** wearers adapted instantly
- **Preferred 3:1** over Ovation® lenses by wearers
- **FrameOptimization™ Technology** for the ultimate in frame choice
- **DualOptix™** digital format, for wider fields of clear vision



If your customers like Essilor Natural lenses or Ovation lenses, they'll love **Accolade** lenses—upgrade to **Accolade** today!

Just add **Accolade®** progressive lenses for...

Instant adaptation and a richer viewing experience



If your presbyopic patients like Essilor Natural® or Ovation® lenses, they'll love **Accolade** lenses!

Accolade®

Accolade Freedom™



Learn more at www.AccoladeLenses.com

Accolade®

Accolade Freedom™


Want your customers to see better in an instant?

Upgrade to **Accolade**[®] or **Accolade Freedom**[™] lenses

Specifically engineered for instant wearer adaptation, **Accolade** and **Accolade Freedom** lenses are the high-performance PALs your patients will love. In wearer tests:

- **9** out of **10** patients adapted instantly (≤ 1 hour) compared to 6 out of 10 with Essilor Natural[®] lenses
- In a **100**-patient study, **Accolade** lenses outperformed Ovation[®] lenses in all vision parameters
- **3:1** preference over Ovation lenses in near, intermediate and distance viewing

Accolade lenses utilize Harmonix[™] Technology, Essilor's proprietary system that takes individual eye shape into account:



Accolade Freedom lenses give your patients the freedom to choose any frame

Accolade Freedom lenses offer all of the benefits of **Accolade** lenses, plus:

- FrameOptimization Technology, which allows patients to pick the frame that's right for them
 - Lens design is customized to the patient Rx, frame dimensions and fitting measurements
 - Fits any frame that's at least 15 mm mfh*
 - No need for additional measurements
- DualOptix[™] Technology
 - Places designs on both sides of the lens, utilizing the entire lens surface
 - 100% digital format
 - Noticeably wider fields of clear vision

*Accolade minimum fitting height (mfh) is 17 mm.



Give your patients a more robust visual experience—upgrade them to **Accolade** or **Accolade Freedom** lenses today.