

Varilux®

Available with Transitions®
Light Intelligent Lenses™

Varilux®

#1 progressive lens brand

worldwide¹



VARILUX® LENSES
AVAILABILITY
for Eyecare Professionals
July 2026

(1) Euromonitor International, 2025 data; Retail value amongst spectacle lenses category, brands representing progressive lenses.

(2) Crizal® coating is an Essilor® lens treatment that acts as an invisible shield of protection.



essilor

evolving
vision

VARILUX® XR DESIGN			LENS ID		VARILUX® XR TRACK										LENS ID		POWER RANGE	ADD POWER*	MAX CYLINDER
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	ESC								
PLASTIC	★ Clear	1.50												-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Essilor Polarized	1.50	✓	✓										-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Xperio UV™ Polarized	1.50	✓	✓										-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® GEN S™	1.50	✓	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive® new generation	1.50	✓	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
TREXA® (TRIVEX®)	★ Clear	1.53												-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® GEN S™	1.53	✓	✓										-10.00 to +6.00	0.75 to 4.00D	-6.00			
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59												-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Essilor Polarized	1.59	✓	✓										-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Xperio UV™ Polarized	1.59	✓	✓										-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓				-10.00 to +6.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive® new generation	1.59	✓	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00			
THIN & LITE®	★ Clear	1.67												-12.00 to +9.00	0.75 to 4.00D	-6.00			
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓				-12.00 to +9.00	0.75 to 4.00D	-6.00			
	Transitions® XTRActive® new generation	1.67	✓	✓	✓									-12.00 to +9.00	0.75 to 4.00D	-6.00			
	★ Clear**	1.74												-14.00 to +13.00	0.75 to 4.00D	-6.00			
	Transitions® GEN S™**	1.74	✓	✓										-14.00 to +13.00	0.75 to 4.00D	-6.00			

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**THIN & LITE 1.74 lenses come systematic with Crizal.

★ Available in Essential Blue Series® lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ Currently Available

VARILUX® X DESIGN			LENS ID		VARILUX® X FIT								LENS ID		VARILUX® X 4D			LENS ID	
																			
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	ESC	POWER RANGE		ADD POWER*	MAX CYLINDER				
PLASTIC	★ Clear	1.50										-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Essilor Polarized	1.50	✓	✓								-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Xperio UV™ Polarized	1.50	✓	✓								-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Transitions® GEN S™	1.50	✓	✓	✓							-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Transitions® XTRActive® new generation	1.50	✓	✓	✓							-10.00 to +6.00		0.75 to 4.00D	-6.00				
TREXA® (TRIVEX®)	★ Clear	1.53										-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Transitions® GEN S™	1.53	✓	✓								-10.00 to +6.00		0.75 to 4.00D	-6.00				
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59										-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Essilor Polarized	1.59	✓	✓								-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Xperio UV™ Polarized	1.59	✓	✓								-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00		0.75 to 4.00D	-6.00				
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓							-10.00 to +6.00		0.75 to 4.00D	-6.00				
THIN & LITE®	★ Clear	1.67										-12.00 to +9.00		0.75 to 4.00D	-6.00				
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓		-12.00 to +9.00		0.75 to 4.00D	-6.00				
	Transitions® XTRActive® new generation	1.67	✓SM	✓	✓							-12.00 to +9.00		0.75 to 4.00D	-6.00				
	★ Clear**	1.74										-14.00 to +13.00		0.75 to 4.00D	-6.00				
	Transitions® GEN S™**	1.74	✓	✓								-14.00 to +13.00		0.75 to 4.00D	-6.00				

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal® Sapphire™ HR, Crizal® Rock™ or Crizal® Easy Pro.

★ Available in Essential Blue Series® lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ **Currently Available**

SM: Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

VARILUX® PHYSIO® EXTENSEE™			LENS ID		VARILUX® PHYSIO® EXTENSEE™ TRACK										LENS ID	
																
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER	
PLASTIC	★ Clear	1.50										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00	
	Transitions® GEN S™	1.50	✓	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.50	✓	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00	
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59										-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00	
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.59	✓	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00	
THIN & LITE®	★ Clear	1.67										-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00	
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓		-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00	
	Transitions® XTRActive® new generation	1.67	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00	
	★ Clear**	1.74										-16.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00	
	Transitions® GEN S™***	1.74	✓	✓								-16.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00	

VARILUX® PHYSIO® EXTENSEE™ CLASSIC EDITION												LENS ID			
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50										-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓							☼	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓							☼	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® GEN S™	1.50	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓							-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® Polarized	1.50	✓										-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D
TREXA® (TRIVEX®)	★ Clear	1.53										-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® GEN S™	1.53	✓	✓	✓							-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.53	✓	✓								-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59										-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓							☼	-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓							☼	-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓	✓	✓							-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® Polarized	1.59	✓										-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50, 4.00D
THIN & LITE®	★ Clear	1.60										-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® GEN S™	1.60	✓	✓	✓							-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	★ Clear	1.67										-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.67	✓	✓							☼	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓							☼	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® Polarized	1.67	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	★ Clear**	1.74										-18.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® GEN S™***	1.74	✓	✓								-18.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal® Sapphire™ HR, Crizal® Rock™ or Crizal® Easy Pro.

★ Available in Essential Blue Series® lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ Currently Available

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue).

VARILUX® COMFORT MAX			LENS ID										VARILUX® COMFORT MAX FIT			LENS ID									
			€													Fe									
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER										
PLASTIC	★ Clear	1.50										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00										
	Essilor Polarized	1.50	✓	✓							☼	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00										
	Xperio UV™ Polarized	1.50	✓	✓							☼	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00										
	Transitions® GEN S™	1.50	✓	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00										
	Transitions® XTRActive® new generation	1.50	✓	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00										
TREXA® (TRIVEX®)	★ Clear	1.53										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00										
	Transitions® GEN S™	1.53	✓	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00										
	Transitions® XTRActive® new generation	1.53	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00										
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59										-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00										
	Essilor Polarized	1.59	✓	✓							☼	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00										
	Xperio UV™ Polarized	1.59	✓	✓							☼	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00										
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00										
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00										
THIN & LITE®	★ Clear	1.60										-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00										
	Transitions® GEN S™	1.60	✓	✓	✓							-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00										
	★ Clear	1.67										-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00										
	Essilor Polarized	1.67	✓	✓							☼	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00										
	Xperio UV™ Polarized	1.67	✓	✓							☼	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00										
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓		-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00										
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-6.00										

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

***THIN & LITE* 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal® Sapphire™ HR, Crizal® Rock™ or Crizal® Easy Pro.

★ Available in [Essential Blue Series®](#) lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ **Currently Available**

SM: Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ **Essilor Sun Colors** (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue).

VARILUX COMFORT DRx™			LENS ID 														
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER		
PLASTIC	★ Clear	1.50										-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D [†]	-6.00		
	Essilor Polarized	1.50	✓	✓							☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D [†]	-6.00		
	Xperio UV™ Polarized	1.50	✓	✓							☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® GEN S™	1.50	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® XTRActive® new generation	1.50	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® XTRActive® Polarized™	1.50	✓										-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D [†]	-6.00	
TREXA® (TRIVEX®)	★ Clear	1.53										-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® GEN S™	1.53	✓	✓	✓							-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® XTRActive® new generation	1.53	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D [†]	-6.00		
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D [†]	-6.00		
	Essilor Polarized	1.59	✓	✓							☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D [†]	-6.00		
	Xperio UV™ Polarized	1.59	✓	✓							☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® XTRActive® Polarized™	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D [†]	-6.00	
THIN & LITE®	★ Clear	1.60										-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® GEN S™	1.60	✓	✓								-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D [†]	-6.00		
	★ Clear	1.67										-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D [†]	-6.00		
	Essilor Polarized	1.67	✓	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D [†]	-6.00		
	Xperio UV™ Polarized	1.67	✓	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® XTRActive® Polarized™	1.67	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D [†]	-6.00		
	★ Clear**	1.74										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D [†]	-6.00		
	Transitions® GEN S™**	1.74	✓	✓								-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D [†]	-6.00		

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevenir®, Crizal® Sapphire™ HR, Crizal® Rock™ or Crizal® Easy Pro.

† 4.00D add only available in Varilux Physio DRX.

★ Available in Essential Blue Series® lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ **Currently Available**

SM: Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ **Essilor Sun Colors** (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue).

VARILUX COMFORT SHORT DRX™

VARILUX COMFORT SHORT DRX™

RX AVAILABILITY FOR VARILUX COMFORT SHORT DRX MATCHES VARILUX COMFORT DRX EXCEPT ADD POWER RANGE IS 0.75 TO 3.00D ACROSS ALL MATERIALS.

MFH=14mm

VARILUX® SHIFT™

LENS ID



MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50										-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D†	-6.00
	Essilor Polarized	1.50	✓	✓							☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D†	-6.00
	Xperio UV™ Polarized	1.50	✓	✓							☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D†	-6.00
	Transitions® GEN S™	1.50	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D†	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D†	-6.00
	Transitions® XTRActive® Polarized™	1.50	✓									-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D†	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53										-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D†	-6.00
	Transitions® GEN S™	1.53	✓	✓	✓							-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50, 4.00D†	-6.00
	Transitions® XTRActive® new generation	1.53	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D†	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D†	-6.00
	Essilor Polarized	1.59	✓	✓							☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D†	-6.00
	Xperio UV™ Polarized	1.59	✓	✓							☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D†	-6.00
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓		-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50, 4.00D†	-6.00
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D†	-6.00
	Transitions® XTRActive® Polarized™	1.59	✓									-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D†	-6.00
THIN & LITE®	★ Clear	1.60										-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D†	-6.00
	Transitions® GEN S™	1.60	✓	✓								-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D†	-6.00
	★ Clear	1.67										-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D†	-6.00
	Essilor Polarized	1.67	✓	✓							☼	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D†	-6.00
	Xperio UV™ Polarized	1.67	✓	✓							☼	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D†	-6.00
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D†	-6.00
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D†	-6.00
	Transitions® XTRActive® Polarized™	1.67	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D†	-6.00
	★ Clear**	1.74										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D†	-6.00
	Transitions® GEN S™**	1.74	✓	✓								-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D†	-6.00

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal® Sapphire™ HR, Crizal® Rock™ or Crizal® Easy Pro.

† 4.00D add only available in Varilux Physio DRx.

★ Available in Essential Blue Series® lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ **Currently Available**

SM: Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ **Essilor Sun Colors** (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue).

VARILUX COMFORT®		LENS ID				POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER
MATERIAL	PRODUCT	INDEX	G	B	GG					
PLASTIC	Clear	1.50				-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® GEN S™	1.50	✓	✓		-10.00 to +6.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓	-10.00 to +6.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59				-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Essilor Polarized	1.59	✓	✓		-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓		-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Transitions® GEN S™	1.59	✓	✓	✓	-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Transitions® XTRActive® new generation	1.59	✓	✓	✓	-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
THIN & LITE®	Clear	1.60				-12.00 to +8.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Transitions® GEN S™	1.60	✓			-12.00 to +8.00	73/78	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Clear	1.67				-14.00 to +9.00	78/83	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Transitions® GEN S™	1.67	✓	✓		-14.00 to +9.00	73/78	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00

MFH=17mm

NOTES: *ADD Power in steps of 0.25D.

G: Gray, B: Brown, GG: Graphite Green

✓ Currently Available

VARILUX COMFORT SHORT™		LENS ID			

RX AVAILABILITY FOR **VARILUX COMFORT SHORT** MATCHES **VARILUX COMFORT**, EXCEPT ADD POWER RANGE IS 0.75 TO 3.00D ACROSS ALL MATERIALS.

MFH=14mm

OCCUPATIONAL LENSES:

VARILUX® IMMERSIA™ MID			LENS ID				VARILUX® IMMERSIA™ ROOM				LENS ID			
			I0								I01			
MATERIAL	PRODUCT	INDEX	G	B	GG	ESC	POWER RANGE		BASE CURVES		ADD POWER*		CYLINDER	
PLASTIC	★ Clear	1.50					-10.00 to +6.00		2.00, 3.50, 4.25, 6.25, 7.25		0.75 to 3.50D		-6.00	
	Transitions® GEN S™	1.50	✓	✓	✓		-10.00 to +6.00		2.00, 3.50, 4.25, 6.25, 7.25		0.75 to 3.50D		-6.00	
TREXA® (TRIVEX®)	★ Clear	1.53					-10.00 to +6.00		2.00, 3.50, 4.25, 6.25, 7.25		0.75 to 3.50D		-6.00	
	Transitions® GEN S™	1.53	✓	✓	✓		-10.00 to +6.00		2.00, 3.50, 4.25, 6.25, 7.25		0.75 to 3.50D		-6.00	
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59					-10.00 to +6.00		1.50, 3.25, 4.25, 5.50, 7.00		0.75 to 3.50D		-6.00	
	Transitions® GEN S™	1.59	✓	✓	✓		-10.00 to +6.00		1.50, 3.25, 4.25, 5.50, 7.00		0.75 to 3.50D		-6.00	
THIN & LITE®	★ Clear	1.67					-14.00 to +9.00		1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D		-6.00	
	Transitions® GEN S™	1.67	✓	✓	✓		-14.00 to +9.00		1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D		-6.00	
	★ Clear	1.74					-20.00 to +13.00		0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00		0.75 to 3.50D		-6.00	
	Transitions® GEN S™	1.74	✓	✓			-20.00 to +13.00		0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00		0.75 to 3.50D		-6.00	

HIGH WRAP LENSES:

VARILUX® STYLISTIC®			LENS ID													
			S													
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	C	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75	0.75 to 3.50 D	-6.00
	Essilor Polarized	1.50	✓	✓	✓							☼	-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.25	0.75 to 3.50 D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓	✓							☼	-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.25	0.75 to 3.50 D	-6.00
	Transitions® GEN S™	1.50	✓	✓	✓	✓	✓	✓	✓	✓			-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50 D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50 D	-6.00
	Transitions® XTRActive® Polarized™	1.50	✓										-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50 D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50 D	-6.00
	Transitions® GEN S™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50 D	-6.00
	Transitions® XTRActive® new generation	1.53	✓	✓									-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50 D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50 D	-6.00
	Essilor Polarized	1.59	✓	✓	✓					✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50 D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓	✓							☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50 D	-6.00
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓			-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00	0.75 to 3.50 D	-6.00
	Transitions® XTRActive® new generation	1.59	✓	✓	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50 D	-6.00
	Transitions® XTRActive® Polarized™	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50 D	-6.00
THIN & LITE®	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50 D	-6.00
	Essilor Polarized	1.67	✓	✓	✓							☼	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50 D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓	✓							☼	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50 D	-6.00
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓			-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50 D	-6.00
	Transitions® XTRActive® new generation	1.67	✓	✓	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50 D	-6.00
	Transitions® XTRActive® Polarized™	1.67	✓										-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50 D	-6.00

MFH=14mm

NOTES: *ADD Power in steps of 0.25D.

★ Available in [Essential Blue Series®](#) lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby, C: Copper, ESC: Essilor Sun Colors

✓ Currently Available

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue).

To experience the full benefit of Varilux® Stylistic® lenses, it's recommended (not required) to provide the frame curve, wrap, panto and vertex when ordering.



For the complete list of Transitions® color availability please refer to the *Transitions* Availability Chart.
For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

©2026 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries. *Transitions*, *Transitions Signature*, *XTRActive*, and the *Transitions* logo are registered trademarks and *Light Intelligent Lenses* is a trademark of Transitions Optical, Inc., used under license by Transitions Optical Ltd. *GEN S* is a trademark of Transitions Optical Limited. **Photochromic performance is influenced by temperature, UV exposure, and lens material.** These products may be protected by one or more patents. 485203_PRO_VAR_SHK/HB 06/26

[essilor.com](https://www.essilor.com)