

VARILUX®

LENSES AVAILABILITY

for Eyecare Professionals

Effective JULY 2021

VARILUX® X DESIGN		VARILUX® X FIT		VARILUX® X 4D																											
LENS ID		LENS ID		LENS ID																											
MATERIAL		PRODUCT		INDEX		G		B		GG		COLORS		G1		B1		G-G		C		R		ESC		POWER RANGE		ADD POWER*		MAX CYLINDER	
PLASTIC	★ Clear	1.50																													
	Essilor Polarized	1.50	✓	✓																											
	Xperio UV™ Polarized	1.50	✓	✓																											
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓																										
	Transitions® XTRActive®	1.50	✓	✓																											
TREXA® (TRIVEX®)	★ Clear	1.53																													
	Transitions® Signature® GEN 8™	1.53	✓	✓																											
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59																													
	Essilor Polarized	1.59	✓	✓																											
	Xperio UV™ Polarized	1.59	✓	✓																											
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓																									
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓																										
THIN & LITE®	★ Clear	1.67																													
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓																									
	Transitions® XTRActive®	1.67	✓SM	✓	✓																										
	★ Clear**	1.74																													
	Transitions® Signature® GEN 8™**	1.74	✓	✓																											

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

★ Available in **Essential Blue Series™** lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.



VARILUX® PHYSIO® W3+ TEMPORAL NASAL

VARILUX® PHYSIO® W3+ FIT TEMPORAL NASAL

VARILUX® PHYSIO® W3+ EYECODE™ TEMPORAL NASAL

MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓			✓						-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.67	✓	✓			✓						-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓			✓						-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.67	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74											-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.74	✓										-18.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00

MFH=14mm

VARILUX® COMFORT MAX LENS ID

VARILUX® COMFORT MAX FIT LENS ID

MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓			✓						-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Essilor Polarized	1.67	✓	✓			✓						-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓			✓						-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

★ Available in **Essential Blue Series™** lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

VARILUX PHYSIO DRx™

LENS ID



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓								-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Vantage®	1.50	✓										-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Vantage®	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.74	✓										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

***For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)

VARILUX PHYSIO SHORT DRx™

LENS ID



RX AVAILABILITY FOR VARILUX PHYSIO SHORT DRx MATCHES VARILUX PHYSIO DRx, EXCEPT ADD POWER RANGE IS 0.75 TO 3.00D ACROSS ALL MATERIALS.

MFH=14mm

VARILUX COMFORT DRX™



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	0.75 to 3.50D	-4.00
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓								-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00
	Transitions® Vantage®	1.50	✓										-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓	✓							-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Transitions® Vantage®	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓		✓							-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	★ Clear*	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00
Transitions® Signature® GEN 8™	1.74	✓			✓							-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00	

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

*Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

Please see Xperio UV™ availability chart for full range of Xperio UV colors including Essilor® Sunwear Colors™.

Trivex is a registered trademark of PPG Industries Ohio, Inc.

★ Available in Essential Blue Series™ lenses.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)

VARILUX COMFORT SHORT DRX™



RX AVAILABILITY FOR VARILUX COMFORT SHORT DRX MATCHES VARILUX COMFORT DRX, EXCEPT ADD POWER RANGE IS 0.75 TO 3.00D ACROSS ALL MATERIALS.

MFH=14mm

